

OPERATING MANUAL

Gilson Tapping Sieve Shaker

SS-8R, SS-12R



SS-8R

Rev: 6/2022

Gilson Tapping Sieve Shaker SS-8R & SS-12R Operating Manual

CONTENTS

Safe	ety Instructions	2
1.0	Introduction	3
2.0	Unpacking & Setup	3
3.0	Operating Instructions	3
	3.1 Sieve Stack & EZ-Clamp Assembly	3
	3.2 Timer Setup & Operation	3
4.0	Maintenance	4
	4.1 Routine Annual Maintenance	4
	4.2 Disassembly	4
5.0	Troubleshooting	5
6.0	Technical Support	5
7.0	Parts Diagram	6
8.0	Replacement Parts	7
9.0	Accessories	8

SAFETY INSTRUCTIONS

Please read these instructions thoroughly to familiarize yourself with the operation of the SS-8R or SS-12R before attempting to run it.

The buyer is responsible for ensuring that users are properly trained, that they are aware of all of the information and instructions in this document, and that they are aware of the potential risks of operating the apparatus. The manufacturer will not be responsible for any damage to people and/or property caused by noncompliance with any instructions in this manual.

NOTE: These instructions are intended only as a guide for general operation of this device and should not be used in place of test protocol. Refer to current ATSM standards for complete and detailed test procedures.

Always ensure the motor and other electrical components are properly configured for your intended use and available power source. The standard Gilson Tapping Sieve Shaker comes with a 1/3hp motor wired for 115V,60Hz. Sieve Shakers can also be ordered with special wirings: 230V,60Hz; 110V,60Hz; and 220V,50Hz. Motors are NOT explosion-proof.

Always use a properly-wired, three-pronged plug, or otherwise ground the machine. Make sure the cord is located where no one will trip or get tangled in it.

Always check electrical wiring for loose connections and for pinched or frayed wiring.

Always disconnect and lock out power supply when the machine is not in use, especially before performing maintenance and repairs.

- ▲ WARNING: This machine operates on an electric current. Improper operation could result in electric shock, electrocution, or an explosion! Motors are NOT explosion-proof!
- ▲ WARNING: Do not wear loose clothing that might be caught in the machine and keep all body parts away from moving parts of the machine. ALWAYS wear safety glasses, hearing protection, and other personal protective equipment while operating, maintaining, or repairing this machine.
- **WARNING:** DO NOT operate the machine without having all guards and covers in place.
- ▲ WARNING: DO NOT perform tasks on the machine other than those for which it was designed. Only use the machine in the manner for which it was intended, as described in this instruction manual.





1.0 INTRODUCTION

This manual applies to units with serial number SS-1812 or higher. For older units, request manual for your serial number; or convert to current design by ordering parts numbered 4P through 8P from this manual.

2.0 UNPACKING & SETUP

2.1 Inspection of Shipment

- 1. The SS-8R and SS-12R weigh approximately 135lb and 162lb, respectively. Use appropriate equipment and manpower to uncrate the Sieve Shaker. Wear safety glasses and work gloves.
- 2. In most instances, the simplest method of uncrating is to cut the carton away from the machine. Leave the carton intact as much as possible so that it can be used to return the machine if necessary.
- Lift machine from skids and carton, using appropriate equipment to lift the machine onto a substantial, solid work surface.
- 4. Lay the Sieve Shaker on its side to install the leveling legs. Turn the legs into the threaded holes in the bottom corners of the outer case. Turn the legs all the way in at this point. Set the Sieve Shaker upright. Level the machine by adjusting the legs until the machine is level and rigid.
 - Further leveling will likely be required once the Sieve Shaker is fully assembled and ready for operation. Relevel when you move the Sieve Shaker, when the material in the test sieves does not remain evenly distributed, or when the machine becomes unstable in operation.
- 5. Mount threaded clamp rods (2R) by placing jam nut ends into the tapped holes in the platform until threads are flush with the platform bottom. Tighten jam nut to platform to hold clamp rods.
- 6. Remove the round spacers from the knobs by slightly pushing on push buttons and removing the spacers. Holding both push buttons, install lid assembly by placing the knobs upon the clamp rods. Once the lid assembly is in desired location, release the push buttons and the assembly will remain in that position.
- 7. Place the acorn nuts (1) on the top of the clamp rods (2R) and tighten until secure.
- 8. Place the Sieve Shaker on a substantial, solid work surface sufficient to provide support for this size machine.
- 9. Re-level prior to operating.

3.0 OPERATING INSTRUCTIONS

Please read understand all safety and operating instructions for the Gilson Tapping Sieve Shaker before putting it into service.

The Gilson Tapping Sieve Shaker efficiently separates most free-flowing materials with particle sizes from No.4 to No.635 (4.75mm to 20µm). Performance on extended size ranges can be determined by experimentation. The test specimen should be large enough to be representative, without overloading any individual test sieve. Overloading may result in incomplete separation or damage to the test sieve cloth. Maximum loading for individual test sieves No.4 and smaller should be about 200g for 8in test sieves and about 450g for 12in test sieves. Test sieves larger than No.4 should be limited to about one particle of material for each available opening.

Gilson's innovation EZ-Clamp System quickly secures sieve stacks of varying heights. This system with integral test sieve cover greatly enhances efficiency when testing multiple samples or using more than one sieve stack.

NOTE: Before conducting a test, confirm that the Sieve Shaker has been setup and leveled properly in accordance with step 4 of Section 2.0. Readjust if necessary.

3.1 Sieve Stack & EZ-Clamp Assembly

- Push the buttons on the knobs. The EZ-Clamp assembly will slide freely up or down the clamping rods. When the buttons are released, the assembly will remain in that position.
- Position the assembly slightly higher than the height of the sieve stack.
- Place the sieve stack on the platform. Push the buttons on the knobs and slide the EZ-Clamp assembly down firmly against the top test sieve.
- Tighten both knobs to securely clamp the sieve stack in place.
- When the test is complete, unscrew the knobs enough to release tension on the stack, then push the buttons and raise the assembly to the desired height to remove the test sieves.

3.2 Timer Setup & Operation

This unit is equipped with an easy-to-operate Gilson interval countdown timer. The timer has a large, 0.6in LED display and will operate in four different modes. It is powered by line voltage and will work on power supplies 100–265VAC, 50/60Hz, with up to 20 amps Inductive or Resistive current.

NOTE: The main device controlled by the timer may be restricted to operating on a more limited electrical supply range. Check the device carefully to ensure compatibility with your electrical supply.

Current timer mode is indicated by the four red LEDs on the timer face:

- A = MMSS (99min:59sec x 1sec)
- $B = HHMM(99hr:59min \times 1min)$
- C = SSSS (9999sec x 1sec)
- D = MMMM (9999min x 1min)

(H is for hours, M for minutes, and S for seconds.) To adjust the timer mode, press and hold both the up and down arrow keys at the same time until the display shows the mode. Once the mode letters are displayed, press the up or down arrow to change modes. Press "START"/"STOP" to accept new mode.

To set the run time, press either the up or down arrow. The first digit on the right hand side will flash in half-second intervals. Press either arrow key to adjust to the desired value. To enter the displayed digit and move to the next, press "START"/"STOP". Once the last digit on the left is entered, the timer is ready to start.

Press "START"/"STOP" to initiate the current run program. Once running, pressing "START"/"STOP" again will pause the timer with the current amount of time remaining on screen. When allowed to time-out, the timer beeps and displays "DONE". Press any key to continue. Setting and Mode values are saved automatically and restored on power-up.

Perform your test. Time required to complete a test will vary depending upon the physical characteristics of the test material. Most separations will be complete in ten minutes or less. Refer to your test specifications, and be consistent.

Disconnect and lock out power to the Sieve Shaker when it is not in use.

NOTE: The SS-8R and SS-12R Sieve Shakers are equipped with a tapping mechanism which makes a loud noise. If the noise is unacceptable, consider using the SSA-805R Sound Enclosure.

SS-8R and SS-12R are counterbalanced to permit free-standing operation with most common test sieve loadings. With very tall or very short sieve stacks, the units may be unstable and move around during operation. Normally this condition occurs only with sieve stacks exceeding about 20in (508mm) in height. This is equivalent to using more than eight test sieves and pan of 2.125in stacking height. If you need to use tall sieve stacks, and your Sieve Shaker moves around, try re-leveling. If stability is a problem with short sieve stacks, add extra nonfunctional test sieves below an extended rim pan.

NOTE: Bolting the machine down is not a substitute for leveling. Bolting the shaker down will cause forces from the unbalanced sieve stack to damage the drive mechanism. This damage is not covered by warranty.

When using 8in test sieves on the SS-12R, always use the platform adapter to compensate for the difference in test sieve weights. This adapter can be used with either a regular or an extended rim pan. When the adapter is positioned with larger ID side up, you can use either pan. Invert the adapter for a tighter fit if using only an extended rim pan. Secure the ring to the test sieve platform with the two cap screws provided.

4.0 MAINTENANCE

Before performing maintenance or repairs on the Sieve Shaker, **ALWAYS** read and understand the safety, operating, and maintenance instructions.

Please provide the serial number and model number of the unit when ordering replacement parts.

4.1 Routine Annual Maintenance

- Apply of a few drops of oil to the motor end bearings.
- Inspect the drive belt for wear, tension, and alignment.

A worn, loose, tight, or misaligned drive belt can affect operation of the Sieve Shaker. The belt should be snug: Neither too tight nor too loose. A snug fit assures longer life, less bearing wear, and quieter operation than a belt which is too tight. A loose drive belt may cause the unit to run too slowly or in spurts. The drive belt should deflect 1/64 of the value of the span of the pulleys. The pulleys should be aligned to avoid excessive edge wear.

Never force or pry the belt over the pulley flanges. Use steps 1–3 in Section 4.2.

• Every two years or whenever disassembled, grease the thrust bearing (13) and the face of the cam (31). Perform steps 1–7 in Section 4.2 to access these parts.

4.2 Disassembly

- 1. Disconnect and lock out the power supply.
- Locate the four cover mounting screws, and remove them. Remove right (timer side) section of the cover.
- Inside the left cover section, locate two more mounting screws on the cover flanges. Remove these screws and the left section of the cover.
- 4. Loosen the four motor mounting bolts. Motor will slide toward the left side of the machine, loosening tension on the belt so that you can remove it.
- 5. Remove the four mounting cap screws (16) from the mounting plate (17), and lift the entire unit out of its case.

- 6. Unhook the two hammer springs (23).
- 7. Remove the connecting link nuts (41) and links (40).
- 8. Remove the hammer post cap screws (25), actuator link screws (36), bumper block cap screws (20), and lock ring (19). Now you can lift the sieve platform with main shaft out of its center housing (15).
- If you remove the cam shaft hangers (27) for any reason, reassemble exactly as removed. The roller clutch block (33) and one cam shaft hanger (27) contain overrunning clutches which must be installed to allow correct rotation of the cam shaft (29), or serious damage will result.
- 10. Main bearings have permanent lubrication and are sealed inside the main housing (15). This main housing with bearings and hub should not be disassembled in the field; replace them as a unit if necessary.

5.0 TROUBLESHOOTING

• Unit Fails to Operate:

Check motor, electrical connections, and timer. Replace or reconnect as necessary. Check drive belt tension; replace belt if worn.

• Unit Runs, but Fails to Give Impact Tapping:

Remove covers, drive belt, and mounting cap screws (16). Replace or reconnect hammer springs (23) as required.

Unit Operates but is Excessively Noisy:

Disassemble through step 4 in Section 4.2. Check resilient faces on bumper block (20) and replace block if necessary.

• Unit is Unstable; Shakes or Walks Too Much:

- Readjust level using the leveling legs, and lock with lock nuts.
- b. For SS-12R, be sure that weighted platform adapter is used when using 8in test sieves.
- c. Loosen and re-clamp the sieve stack.
- d. Check sieve stack height. If over 20in, see Section 3.0.

6.0 TECHNICAL SUPPORT


Contact Gilson Technical Support for assistance with applications, operation, maintenance, or repair.

• Telephone: 800-444-1508

Email: techsupport@gilsonco.com

• Web: globalgilson.com

7.0 PARTS DIAGRAM



8.0 REPLACEMENT PARTS

Item No.	Description	No. Req.	Part Number	Item No.		No. Req.	Part Number
EZ-C	CLAMP SIEVE CLAMPING PARTS	S		HAMI	MER ASSEMBLY PARTS		
1	Clamp Rod Acorn Nut	2	RPSS-8RKEY1	21	Upper Hammer	1	RPSS-12RKEY21
2	EZ-Clamp Rod	2	RPSS-12RKEY2R	22	Lower Hammer	1	RPSS-12RKEY22
3	Clamp Rod Nut	2	RPSS-12RKEY3	23	Hammer Spring	2	RPSS-12RKEY23
4	EZ-Clamp Knob	2	RPSS-8RKEY4R	24	Hammer Post	2	RPSS-12RKEY24
6	EZ-Clamp Compression Test Sieve	1	RPSS-8RKEY6R	25	Hammer Post Cap Screw and Washer	2	RPSS-12RKEY25
	Cover for SS-8R	·	00 0	26	Hammer Post Nonmetallic Bushing	4	RPSS-12RKEY26
6	EZ-Clamp Compression Test Sieve Cover for SS-12R	1	RPSS-12RKEY6R		SHAFT ASSEMBLY PARTS Cam Shaft Hanger	0	DDCC ODVEVOZ
7	EZ-Clamp Washers and Retaining	2	RPSS-8RKEY7R	27	Hanger Cap Screw and Lock Washer	2 4	RPSS-8RKEY27 RPSS-12RKEY28
0	Ring	0	DD00 0D1/E)/0D	29	Cam Shaft	1	RPSS-8RKEY29
8	EZ-Clamp Push Button and Spring	2	RPSS-8RKEY8R	30			
OLD :	STYLE CLAMPING PARTS, SN 5	499 AN	ID LOWER	31	Cam Shaft Hanger Bushing	1	RPSS-12RKEY30 RPSS-12RKEY31
N/A	Clamp Center Knob with Screw and	1	RPSS-12RKEY5P	32	Cam Paller Clutch	1 2	
IN/A	Washers	ı	NF33-12NKE13F	33	Roller Clutch	1	RPSS-12RKEY32
N/A	Clamp End Knob with Washers and	2	RPSS-12RKEY7P	34	Roller Clutch Block Cam Shaft Spacer	1	RPSS-12RKEY33 RPSS-12RKEY34
NI/A	Loose Cap Clamp Bar Tube For SS-8R	1	RPSS-8RKEY4P	34A	Cam Shaft Nylon Washer	4	RPSS-12RKEY34A
N/A N/A	Clamp Bar Tube for SS-12R	1	RPSS-12RKEY4P	35 35	•		
	•			35	Actuator Link with Rod Ends (2) and Nuts (2)	1	RPSS-12RKEY35
N/A	Clamp Bar End Block	2	RPSS-12RKEY8P	36	Actuator Link Screw	2	RPSS-12RKEY36
N/A	Clamp Rod	2	RPSS-12RKEY2	37	Stabilizer Rocker Bracket	1	RPSS-8RKEY37
N/A	Compression Test Sieve Cover for SS-8R	1	RPSS-8RKEY6P	38	Stabilizer Rocker Bracket Nylon Bushin	g 2	RPSS-12RKEY38
N/A	Compression Test Sieve Cover for SS-	1	RPSS-12RKEY6P	39	Stabilizer Rocker Bracket Bolt and Nut	2	RPSS-12RKEY39
	12R			40	Connecting Link	2	RPSS-12RKEY40
EZ-C	LAMP UPGRADE KITS			41	Connecting Link Nut	4	RPSS-12RKEY41
	EZ-Clamp Upgrade for SS-8R	_	SSA-807	42	Connecting Link Washer	4	RPSS-8RKEY42
	EZ-Clamp Upgrade for SS-12R	_	SSA-809	43	Connecting Link Rubber Bushing	4	RPSS-12RKEY43
MAIN	SHAFT PARTS			44	Stabilizer Arm	1	RPSS-12RKEY44
9	Sieve Platform with Main Shaft for	1	RPSS-8RKEY9	OUTE	R CASE PARTS		
5	SS-8R	'	TH OO OF INCLIS	48	Leveling Leg, Lock Nut	4	RPSS-12ROCP6
9	Sieve Platform with Main Shaft for	1	RPSS-12RKEY9	51	Outer Case	1	RPSS-12R-OCP-1
10	SS-12R Drive Pulley, Counterweight, Key,	1	RPSS-8RKEY10	52	Left and Right Case Cover for SS-8R	1	RPSS-8ROCP2
10	Setscrew for SS-8R	'	111 33-011KL110	52	Left and Right Case Cover for SS-12R	1	RPSS-12ROCP3
10	Drive pulley, Counterweight, Key, Setscrew for SS-12R	1	RPSS-12RKEY10	53	Cover Mount Screw	4	RPSS-CMS
11	Neoprene Sponge Washer	1	RPSS-8RKEY11	ELEC	TRIC & DRIVING PARTS		
12	Anti-friction Washer	1	RPSS-8RKEY12	46	Motor 115V,60Hz	1	RPSS-8REDP1
13	Thrust Bearing	1	RPSS-12RKEY13	46	Motor 230V,50Hz	1	RPSS-MOTOR
14	Nonmetallic Bushings	2	RPSS-12RKEY14				230/50
15	Main Housing, Bearing, and Hub	1	RPSS-12RKEY15	47	Motor Cord, Grommet, and Plug (Specify Serial Number)	1	RPSS-12R-M CORD
16	Mounting Cap Screw and Lock Washer	8	RPSS-12RKEY16	N/A	Motor Mounting Bolt, Flat Washer, Lock Washer	1	RPSS-MMBW
17	Mounting Plate	1	RPSS-8RKEY17	49	Motor Pulley 60Hz	1	RPSS-PULLEY/60HZ
18	Lower Counterweight, Key, Set- screws for SS-8R	1	RPSS-8RKEY18	49	Motor Pulley 50Hz	1	RPSS-PULLEY/50HZ
18	Lower Counterweight, Key, Set- screws for SS-12R	1	RPSS-12RKEY18	45 50	Drive Belt Timer (Serial Numbers 3315 and Newer	1	RPSS-12REDP7 RP-Timer-ART
19	Lock Ring	1	RPSS-8RKEY19	30	Timor (Octial Nationals 3313 and Newel	, '	III - IIIIIGI-ANI
20	Bumper Block, Cap Screws (2), Lock Washers (2)	1	RPSS-12RKEY20				

9.0 ACCESSORIES



SSA-807



SSA-809 shown on SS-12R with Test Sieves



ACCESSORIES						
Description	Model					
EZ-Clamp Upgrade Kit replaces the original clamping assemblies on older Gilson Tapping Sieve Shakers. Knurled knobs with push-button release slide freely up and down the clamp rods for smooth, easy clamping. Once in position, a quick twist tightly secures the sieve stack. When the test is complete, push the EZ-Clamp button and raise just enough to remove the stack. Upon release of the button, the clamps stay in place, ready for the next test. EZ-Clamp kits include free-sliding push-button knobs, an integral sieve cover, and new clamp rods. EZ-Clamp Upgrade Kit for SS-8R EZ-Clamp Upgrade Kit for SS-12R	SSA-807 SSA-809					
Gilson Sound Enclosure controls noise and dust associated with SS-8R and SS-12R Sieve Shakers and other lab equipment. Sturdy painted steel case with full-width hinged doors is lined with 1 in (25.4mm) of sound-attenuating foam. Product Dimensions: 31x19x46in (787x482x1,168mm), WxDxH.	SSA-805R					
Clean-N-Stor accessories are handy, time-saving devices for emptying, cleaning, and weighing functions associated with sieving operations. Inverting an 8in or 200mm test sieve on the stainless steel funnel allows quick emptying and cleaning of contents into a receiving scoop or pan. A sieve stack can also be stored on top of the funnel. A scoop and soft-bristle cleaning brush are included with all models. The SSA-801 attaches to the top of the SS-8R case. SSA-802 is a stand-alone model that can be positioned directly over an electronic balance, so sieve fractions can be weighed as the test sieve is being cleaned. OBA-15R is an adjustable-height Clean-N-Stor version designed to fit over taller balances.						
Clean-N-Stor Attachment for SS-8R Stand-Alone Clean-N-Stor Adjustable-Height Clean-N-Stor	SSA-801 SSA-802 OBA-15R					
Platform Adapters permit smaller diameter test sieves to be used in SS-8R and SS-12R Gilson Tapping Sieve Shakers. Each adapter is designed to compensate for weight differences as well as frame diameter. SSA-810 and SSA-811 fit either sieve shaker, and must be used in conjunction with the included adapter when using with the SS-12R. The SSA-812 for 10in (254mm) test sieves is for use with the SS-12R only.						
Platform Adapter for 3in (76.2mm) Test Sieves Platform Adapter for 6in (152.4mm) Test Sieves Platform Adapter for 10in (254mm) Test Sieves	SSA-810 SSA-811 SSA-812					



SSA-801 shown with Test Sieves on SS-8R



SSA-802 shown with Test Sieve



SSA-811 shown with Test Sieves