ROTAVAPOR

RE-111 RE-121 EL-131

95068 A4 3000 185

Inhaltsverzeichnis		Contents	Index		Seite Page
Technische Spezifikationen		<i>Technical specifications</i>	Spécifications techniques		3
Zusammenbau		<i>Assembling</i>	Assemblage		4
Hinweise und Tips		<i>Notes and tips</i>	Conseils et recommandations		6
Glasaufbauten	RE-A	<i>Glass assemblies</i>	<i>RE-A</i>	Montages en verre	RE-A
	RE-B		<i>RE-B</i>		RE-B
	RE-C		<i>RE-C</i>		RE-C
	EL-S		<i>EL-S</i>		EL-S
	EL-E		<i>EL-E</i>		EL-E
	EL-C		<i>EL-C</i>		EL-C
Aggregate	RE-111	<i>Drive units</i>	<i>RE-111</i>	Blocs d'entraînement	RE-111
Elektronik-Schema	RE-111	<i>Wiring diagram</i>	<i>RE-111</i>	Schéma de câblage	RE-111
	RE-121		<i>RE-121</i>		RE-121
	EL-131		<i>EL-131</i>		EL-131
Elektronik-Schema	RE-121/EL-131	<i>Wiring diagram</i>	<i>RE-121/EL-131</i>	Schéma de câblage	RE-121/EL-131
Wasserbad	B-461	<i>Waterbath</i>	<i>B-461</i>	Bain marie	B-461
Ölbad	B-471	<i>Oilbath</i>	<i>B-471</i>	Bain d'huile	B-471
Servo-Schnellheber	B-011	<i>Servo jack</i>	<i>B-011</i>	Elévateur servo	B-011
Hebestativ		<i>Lifting stand</i>		Statif élévateur	30
Temperatursonden		<i>Temperature probes</i>		Sondes de température	32
Zubehör		<i>Accessories</i>		Accessoires	34

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


Technical specifications	R-111	RE-121	EL-131	Waterbath 461	Oilbath 471
Drive unit		Sparkless induction motor, constant torque, electronically regulated		-	-
Rotation		stepless speed control up to 220 rpm		-	-
Vacuum tightness		<1 mbar		-	-
Digital LED display	-	Speed of rotation		-	-
	-	bath temperature		-	-
	-	vapour temperature (Probe is optional)		-	-
Voltages	220 V/50 Hz	240 V/50 Hz	117 V/60 Hz	220 V, 117 V	220 V, 117 V
Power consumption		85 W		1200 W	1200 W
Temperature range				20–100°C	30–180°C
Dry run dual safety switch				standard	standard
Level control				standard	-
Bath diameter				265 mm	265 mm
Bath volume				9 l	9 l

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Assembling

- Loosen lower clamp and place drive unit onto the servo jack or lifting stand.
Tighten clamp.
- Inclination angle of rotation axis should be 25°. Check the marks. If necessary loose upper clamp and tilt until the angle marks are congruent.
- Lift the Jacquard-handle as shown for additional turning.
- Put bath onto the base plate.
Insert the fifth plastic foot (18689) if the bath has to be displaced more than 5 cm to the right.
- Mount glass assembly according to sketches on pages 7 to 10.
- Adjust variable stopper such that no rotating part touches the bath bord.
- Connect the feed water hose to the small and the drain hose to the large connection of the bath.

For small flasks put level adjustment in a higher position:
Turn overflow tube to loose it.

7. Cable


The cable may be fixed at the pins located at the rear of the drive unit.
For access, push the plastic cable plate backwards from the drive unit.
The spiro-hose (00728) is wrapped around the probe cables to secure them additionally. See illustration.

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
18915 O-Ring GL 14



18916 Olive GL 14, gebogen

Assembly B**Glass assembly B compl. for RE***Gasket KD-22**Cock, STJ 19/38**PTFE tube, 210 mm**Distribution piece glass**Condenser**Vacuum connection**Screw cap with insert spring**Introduction stopcock compl.**Distribution piece with gasket and screw cap**KS-Clip**Vapour duct**Scaffolding for Assembly B compl., consisting of:**Rod, 300 x 12**Condenser clamp with rod**Cross piece**Saveclip 29****Hose adapter GL 14 compl.****Screw cap GL 14**O-ring GL 14**Hose stem GL 14, bent***Assemblage B****Assemblage B compl. pour RE***Joint KD-22**Robinet, CN 19/38**Tuyau PTFE, 210 mm**Tête de distribution verre**Réfrigérant**Raccord de vide**Couplage avec ressort**Robinet d'introduction compl.**Tête de distribution avec joint et couplage**KS-Pince**Conduit de vapeur**Support pour assemblage B compl., composé de:**Tige, 300 x 12**Bride pour réfrigérant et tige**Pièce en croix**Agrafe de sécurité 29****Olive GL 14 compl.****Ecrou chapeau GL 14**Joint torique GL 14**Olive GL 14, coudée*

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



Cold Trap 2-parts



Glass assembly C compl. for EL
O-ring, 88,3 × 5,3
<i>Cold finger</i>
<i>Cock, STJ 19/38</i>
<i>STJ 19 clip</i>
<i>Distribution head, glass only</i>
<i>Distribution head compl. with drain valve</i>
<i>Screw cap, di 51,5 mm</i>
<i>Cold trap mantle</i>
<i>Scaffolding for Cold trap compl.</i>
<i>Cold trap 2-parts compl. (without distribution head 11130)</i>

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




Wiring diagram RE-121/EL-131

The printed circuit boards are delivered as a standard exchange item. Full price applies until we receive the faulty print for credit.

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Digital display

300 or a higher number on display in position 'Bath' or 'Temp.' means:


- no temperature probe connected or
- probe line interrupted

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






Caution!
Only authorized service technicians
are entitled to open the 18644
assembly containing a loaded
spring.

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


Evaporating flasks (pear-shaped)
For EL systems only with adapter
flange-STJ

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
Adapter flange-STJ
for EL systems
1 adapter STJ 29/32 is supplied
with the ROTAVAPOR

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



Flanged drying flask
for EL systems

Inhalt/Capacity/Capacité	Code
500 ml	00459
1000 ml	00460
2000 ml	00461

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
Destillierspinnen

**Specila flasks for drying
powdery substances**

Inhalt/Capacity/Capacité	NS/STJ/CN 29/32	NS/STJ/CN 24/
500 ml	00452	115
1000 ml	00453	004
2000 ml	00454	115

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


Cooling coil insert

This coil can be inserted in the cold-trap condenser in place of the cold finger and can be fixed by means of a quick fastener.

With this accessory, the ROTAVAPOR can now be operated with a cold-trap condenser or with a normal condenser, as desired.

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Flanged flasks

for flanged vapour duct

The wide flask neck makes it easy to remove substances from the evaporating flask.

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

Stand assembly adapter

This adapter enables the drive unit to be mounted on a laboratory stand having a rod of diameter 12 to 15 mm.

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
ROTAVAPOR®

BÜCHI

Laboratory Techniques Ltd.
CH - 9230 Flawil/Switzerland

Tel. 071 84 8181, Telex 881230 buchi ch, Telefax 071 83 57 11


Introduction

It is now hardly possible to consider rapid and careful evaporation of solvents in the laboratory without ROTAVAPOR (vacuum rotary evaporator).

The first ROTAVAPOR for professional users was originally produced by BÜCHI in the fifties.

Today's ROTAVAPOR which combines traditional craftsmanship with up-to-date technology, satisfies **your** requirements too. In addition to the conventional evaporation of solvents, the number of possible applications of the ROTAVAPOR has greatly increased.

In particular, the first fully automatic ROTAVAPOR RE-140/EL-141 opens up new horizons.

Advantages of the ROTAVAPOR

When BÜCHI designs a new system in collaboration with the users, very clear goals are first set. These were so ambitious that new standards were applied to many points:

Problem-solving

The ROTAVAPOR, from the robust standard model RE-111 up to the fully automatic RE-140/EL-141, takes full account of **your** personal needs too. In particular, the comprehensive range of accessories make it possible to carry out a large number of applications.

Finest quality

The ROTAVAPOR offers you top performance and reliability at reasonable cost. The BÜCHI philosophy is to give the best and not just the necessary quality. Quality also includes advice for special needs and a first-class customer service.

Safest ROTAVAPOR

All baths are designed as safety water baths. An additional thermal cut out prevents the worst if you have switched the bath on without liquid. A standard feature at no extra cost.

The glass components are provided with threaded safety connectors.

The fully automatic ROTAVAPOR RE-140/EL-141 has a remote control for those who want to remain "at a distance".

All systems are supplied with the BÜCHI STJ combi clip. It is permanently fitted and retains the flask even at normal pressure.

Balanced

A special exclusive feature is the servo quick-action jack. In the B-011 version, it can be raised and lowered almost without effort; Swiss precision work, pioneering in a new direction!

The quick-action jack B-143 belonging to the ROTAVAPOR RE-140/EL-141 is of course motor-driven. Automatic raising or lowering by pressing a key or according to the distillation program.

Application of the ROTAVAPOR

The popularity of the ROTAVAPOR is also based on its versatility. The ROTAVAPOR adapts to changing needs or can be suitably modified. The two different drive units can be provided with a total of six tried-and-tested glass assemblies. Their application can be optimised by the extensive range of accessories.

The following areas of application of the ROTAVAPOR can also be combined, e.g. extraction/concentration with simultaneous solvent regeneration.

- Distillation
- Concentration
- Reaction
- Powder drying
- Coating
- Crystallisation
- Solvent purification
- Extraction
- Degassing

General specifications

Drive unit

Powerful sparkless induction motor, 0–220 rpm, electronically regulated. The torque is proportional to the load. Power input 85W.

Seals

Durable PTFE seal vulcanised on elastic rubber compound. Vacuum-tight to <1 mbar.

Baths

Ø 265 mm, capacity 9 litres, overtemperature cut out with double circuit system, stainless steel body.
Oil bath 30–180° C, 1200W
Water bath 20–100° C, 1200W with top-up, overflow and drain function

Glassware

First-class, chemically inert borosilicate glass.

Flasks

Both evaporating and receiving flasks can be obtained with capacities from 50 ml to 3 litres. The units are supplied as standard with one 1 litre receiving flask SPJ 35/20 and one evaporating flask STJ 29/32. Naturally, you will also find all other standard ground joints in our range (see page 13).

Connectors

- On the glass components:
- threaded safety connectors Ø 9 mm
 - On the water bath:
 - inlet 8 mm
 - outlet 14 mm

Cooling water consumption

Approx. 1.5 litres per minute

Raising distance

The raising/lowering distance of the quick-action jack is 150 mm

Shipping weight




20–25 kg, according to model

For your use: The optimum assembly

Models RE-111, 121 and 140

The RE models have a fixed seal, the KD-22, and a rotating vapour duct. The glass assemblies A, B, C which are possible with this system are suitable for concentration/distillation. Flasks with an internal neck diameter of 50 mm, and a corresponding vapour duct, are obtainable as accessories.

A Assembly

This assembly is characterised by the diagonal condenser. It has widespread use for the **simple distillation** of organic solvents.

The vapours pass from the evaporating flask through the vapour duct into the condenser. The vapours condense on the tap-water-cooled coils. The condensate runs into the receiving flask.

The stopcock on the condenser is used for continuously filling the evaporating flask and for aerating the unit.

B Assembly

The vertical condenser of the B assembly is cooled with tap-water. This configuration requires 10 cm less bench space than the A assembly. Also, the vapour duct is shorter. The B assembly is therefore particularly suitable for **stripping solvents** with high heat of evaporation, for example water and methanol.

The vapours pass from the evaporating flask through the short vapour duct. The vapours already condense in the distribution head or on the condenser coils. The liquid is collected in the round-bottomed flask.

The stopcock on the distribution head is used for feeding solvents and for aerating the system.

C Assembly

The two-piece cold-trap condenser can be filled with dry-ice, cooling brine or liquid air. This enables **solvent freezing**. The applications of the cold-trap condenser are: Distillation at low temperatures, freeze-drying.

The cold finger can be removed without tools for recovery of the distillate and cleaning.




There is no provision for continuous feeding of the evaporating flask.

The vapours pass through the short vapour duct into the cold-trap condenser. The distillate can freeze out on the cold finger or run as liquid into the receiving flask.

Models EL-131 and 141

The EL glass assemblies S, E, C permit the following applications in addition to distillation and concentration: reactions under reflux, degassing of samples, Soxhlet extractions. In the EL models, the vapour duct is fused to the distribution head and stationary. The KD-26 seal

rotates. Evaporating flasks with a flange of internal 35 mm are obtainable as accessories and can be fixed directly to the unit.

S Assembly

The vertical condenser is cooled with tap-water. A PTFE valve is fitted to the distribution head. The valve is opened for distillation. The vapours pass through the distribution head, condense in the condenser and run freely into the receiving flask.

If the receiver valve is closed, the condensate runs back into the evaporating flask; this is called **distillation under total reflux**.

With total reflux, solutions can be degassed and a very large number of reactions can be carried out under vacuum or at normal pressure. This offers the possibility of evaporating off solvent at any time by opening the receiver valve. Using the glass stopcock, reactant or solvent can be fed into the evaporating flask.

A Soxhlet extraction unit can be supplied as an accessory. This is inserted between the distribution head and the condenser. Extraction takes place when the receiver valve is closed and the extract is concentrated when the valve is opened.

E Assembly

An **expansion vessel** is fitted in front of the descending condenser. The vapours pass from the evaporating flask into the distribution head and rise into the expansion vessel, which acts as a buffer vessel for any foam reaching this point. By widening the cross-section, the vapours are expanded and their speed reduced. Aerosols and drops carried over are separated off and run back into the evaporating flask.

The vapours must pass through the whole length of the condenser so as to reach the receiving flask as condensate. The distribution head has a stopcock for continuous feed and for aerating the unit.

The main applications of the E assembly are: **Distillation of solvents with low heat of evaporation, high-foaming solvents and solvent recovery.**

C Assembly

In this assembly, the **cold-trap condenser** replaces the conventional coil condenser. The inner part of the cold-trap condenser, namely the **cold finger**, receives the coolant, for example dry-ice, liquid air or cooling brine.

The applications of this C assembly are: **Distillation at low temperatures, freeze-drying, condensation of low-boiling solvents.**

The distribution head has a feed/aerating stopcock next to the PTFE receiver valve.


The vapours pass from the evaporating flask through the distribution head into the cold-trap condenser. Products which freeze out deposit on the cold finger. Liquid condensate run as into the receiving flask or, if the valve is closed, back into the evaporating flask.

ROTAVAPOR RE-111

The most popular ROTAVAPOR is the Model RE-111. It forms part of the standard equipment in modern chemical laboratories throughout the world. The ROTAVAPOR RE-111 is the result of decades of BÜCHI experience in the construction of such systems; in other words: it was shaped through distillations.


The mechanical construction of the drive unit is made to suit the RE glassware A, B, C. The ROTAVAPOR RE-111 can be used either with the water bath B-461 or with the oil bath B-471. Both are suitable for heating evaporating flasks from 50 ml to 3 litres.

The advantages of the ROTAVAPOR RE-111 are enhanced by the servo quick-action jack B-011; safe, problem-free distillation in regular daily use is guaranteed.

Ordering Information

The units are supplied with one 1 litre evaporating flask STJ 29/32 and one receiving flask SPJ 35/20.


Accessories

The accessories suitable for this system are listed on pages 12-15.

ROTAVAPOR RE-121

The large green LED show you, as required, the rpm, the bath temperature and, with an additional sensor, the vapour temperature. A precision of 1% is achieved by modern electronics in a compact structure. The ROTAVAPOR RE-121 offers quality at reasonable cost. The drive unit is suitable for the RE glass assemblies A, B and C. Either the water bath B-461 or the oil bath B-471 can be used, according to the application. For

this unit, we recommend especially the servo jack B-011. ROTAVAPOR RE-121, digitised distillation!

The units are supplied with one 1 litre evaporating flask STJ 29/32 and one receiving flask SPJ 35/20. The sensor for measuring the vapour temperature must be ordered separately: see Accessories.

Ordering Information

ROTAVAPOR RE-121	Glass Assembly	220-240 V 50-60 Hz	100-117 V 50-60 Hz	Dimensions mm
With servo jack B-011 and water bath B-461	A	19155	19157	Width 75
With servo jack B-011	A	19152	19152	Height 65
With lifting stand	A	19143	19151	Depth 35
With servo jack B-011 and water bath B-461	B	19164	19166	Width 62
With servo jack B-011	B	19161	19163	Height 100
With lifting stand	B	19158	19160	Depth 35
With servo jack B-011 and water bath B-461	C	19173	19175	Width 58
With servo jack B-011	C	19170	19172	Height 80
With lifting stand	C	19167	19169	Depth 35

Accessories

The accessories suitable for this system are listed on pages 12-15.


ROTAVAPOR EL-131

We recommend the ROTA VAPOR EL-131 for applications which are not restricted only to distillation and concentration. The glass assemblies S, E and C can be used with the ROTA VAPOR EL-131.

The rpm, the bath temperature and, with an optional sensor, the vapour temperature are reversibly displayed with clearly visible, green LED elements.

The most modern electronics — only components of «process quality» are fitted — offers you a precision of 1%.

The ROTA VAPOR EL-131, with its variety of possible applications, is the first choice for synthesis laboratories. Oil bath or water bath can be added to the ROTA VAPOR EL-131. For routine use, the unique servo quick-action jack B-011 is of course a must!

**Reaction and distillation:
ROTA VAPOR EL-131
with digital display at an
affordable cost!**

The units are supplied with one 1 litre evaporating flask STJ 29/32 and one receiving flask SPJ 35/20.
The sensor for measuring the vapour temperature must be ordered separately: see Accessories.

Ordering Information

ROTA VAPOR EL-131

With servo jack B-011 and water bath (B-040)

With servo jack B-011

With lifting stand

With servo jack B-011 and oil bath (B-040)

With servo jack B-011

With lifting stand

With servo jack B-011 and water bath (B-040)

With servo jack B-011

With lifting stand

Accessories


The accessories suitable for this system are listed on pages 12–15.

ROTAVAPOR RE-140/EL-141

These two models differ in the glass components which they use. The RE-140 can be combined with the RE assemblies A, B or C and the EL-141 can be combined with the EL glassware S, E or C.

The ROTA VAPOR RE-140/EL-141 is the first rotary evaporator to permit fully automatic distillation. The use of the microprocessor makes it possible to carry out distillations according to fixed programs. Particular

attention has been paid to make it easy to use: rotation speed and distillation pressure are selected via rotary potentiometers without anything having to be keyed in. In this instance, it is more practical than keying-in. Other functions, for example motor-driven raising/lowering of the quick-action jack, are effected by keying-in. The keys do not have any confusing dual functions and, when keying-in, the confirmatory click is felt with the fingertip.

Main features of ROTA VAPOR RE-140/EL-141

- Fully automatic distillation
- Integrated vacuum controller
- User-friendly electronics
- Motor-driven quick-action jack
- Six possible glass assemblies

Automatic distillation

An automatic distillation is best described with the aid of a practical example: a substance dissolved in water loses its biological activity at 60°C. The substance must be concentrated to complete dryness by evaporation, quickly and under mild conditions.

a) Bath temperature

In the example, you set this to 55°C which is below the temperature which is harmful to the substance.

On pressing **Desired bath temperature [°C]**, the previously used temperature appears in the display. If this value is other than 55, you press **Set bath temperature [°C]** and key in 55.

The diode **Heating** lights up when current is flowing through the heating system. The present temperature is displayed by means of the key **Actual bath temperature [°C]**.

If for any reason the bath temperature rises 5°C above the desired temperature, the system automatically switches off and a fault is displayed.

b) Distillation vacuum

A difference between bath temperature and boiling point of 20°C ensures a sufficiently large evaporative capacity. Accordingly, with a bath temperature of 55°C minus 20°C, the boiling point must be 35°C. This is obtained at a pressure of 30 mbar.

With the rotation switched on, you actuate the key **Controlled vacuum** and the pressure is reduced to the preset value by means of the potentiometer **Vacuum**. To display the present vacuum, the key **Actual vacuum [mbar]** is pressed.

The key **Aerate** aerates the unit while it is being pressed. With this key, you can at any time immediately suppress foaming which occurs during evacuation.

The desired pressure of 30 mbar is now kept constant automatically.

The Distillation


The object is to evaporate the sample to complete dryness.

By pressing the key the flask is lowered into the bath by means of the motor. Instead of the bath temperature, we can now display the vapour temperature with the key **Vapour temperature [°C]**. This will adjust itself to the desired 35°C. As soon as this temperature has been reached, we now select the mode of operation. Constant vapour temperature with the key **Controlled vapour temperature**. The ROTAVAPOR now operates not under constant vacuum, but changes the vacuum, as required, in order to keep the vapour temperature constant.

c) End of Distillation

As soon as all the solvent has evaporated off, the ROTAVAPOR must stop operating. As required, the sample has been evaporated rapidly, under mild conditions and to dryness.

In every distillation, the vapour temperature drops when no more solvent evaporates. The microprocessor therefore monitors the vapour temperature continuously. If this drops over a certain time, the ROTAVAPOR is switched off. Switching-off procedure consists of: stop rotation, aerate, raise evaporating flask out of the bath, short acoustic signal. You engage this function by pressing the key **Automatic stop**.


Necessary Accessories: These items are not included as standard item and must be ordered separately.

1 Water-jet pump unit, Code 23349, or when using a mechanical vacuum pump: 1 Vacuum valve with 5 mm orifice, Code 27259.

Ordering Information

ROTAEVAPOR RE-140	Dimensions mm (in.)
With motor-driven jacket B-142 and water bath B-143	Width 600 (23.6) Height 600 (23.6) Depth 500 (19.7)
With motor-driven jacket B-142 and oil bath B-144	Width 600 (23.6) Height 600 (23.6) Depth 500 (19.7)
With motor-driven jacket B-142 and oil bath B-144	Width 600 (23.6) Height 600 (23.6) Depth 500 (19.7)
With motor-driven jacket B-142 and water bath B-143	Width 600 (23.6) Height 600 (23.6) Depth 500 (19.7)
With motor-driven jacket B-142 and oil bath B-144	Width 600 (23.6) Height 600 (23.6) Depth 500 (19.7)
ROTAEVAPOR RE-140	Dimensions mm (in.)
With motor-driven jacket B-142 and water bath B-143	Width 600 (23.6) Height 600 (23.6) Depth 500 (19.7)
With motor-driven jacket B-142 and oil bath B-144	Width 600 (23.6) Height 600 (23.6) Depth 500 (19.7)
With motor-driven jacket B-142 and oil bath B-144	Width 600 (23.6) Height 600 (23.6) Depth 500 (19.7)
With motor-driven jacket B-142 and water bath B-143	Width 600 (23.6) Height 600 (23.6) Depth 500 (19.7)
With motor-driven jacket B-142 and oil bath B-144	Width 600 (23.6) Height 600 (23.6) Depth 500 (19.7)

Accessories for ROTAVAPOR 111, 121/131



Water bath B-461

The water bath B-461 has a diameter of 265 mm and can take flasks up to 3000 ml. The 1200 Watt heating system keeps the temperature constant, even with large heat extraction. The temperature range is 20 to 100° C. Temperature accuracy $\pm 1^\circ \text{C}$. The stainless steel body has an incurved edge so that the bath liquid does not splash over the edge. The water bath B-461 has a continuous feed as well as an overflow and drain device. The latter is constructed so that the level cannot fall below the heating coil. As an active safety feature, the B-461 uses a double circuit system, with two sensors and two separate switching circuits.


If distilled or deionised water is used in the water bath, borax must be added.

Ordering information

Oil bath B-471

The oil bath can be used in the range from 30 to 180° C. Temperature accuracy $\pm 1^\circ \text{C}$. Thermal oil, water-soluble polyalkylene glycol or PEG can be used as the bath liquid. The oil bath is of course equipped with a safety circuit. This consists of a second probe and a second switching circuit. The diameter of the bath is 265 mm, the body is made of stainless steel and the heating power is 1200 Watt. Flasks of 30 to 3000 ml can be used with the oil bath B-471.

Servo Jack B-011

The servo jack B-011 is the ideal complement for the 111, 121 and 131 units.


The new servo jack is exceptionally easy to use and extremely stable. This masterpiece of Swiss precision is designed for hard daily use.

The glass assemblies are raised and lowered with precision, without effort and without jolts. This unique solution is afforded by the patented energy accumulator. The vertical adjustment is 150 mm. The water or oil bath, added as an autonomous module, makes the ROTAVAPOR compact.

Ordering information

Ordering code


Lifting stand

This robust lifting stand is suitable for raising/lowering the 111, 121/131 units. The vertical adjustment distance is 150 mm. The BÜCHI heating baths naturally fit into the stand. The lifting stand has no servo assistance for power boosting.

Ordering code


Accessories

Vapour duct

Assembly	mm	Standard joints	with Cambi-clips			
	29,2/32	24/40	34,5/35	29,2/42	24/29	19,8/26
A	280	23733				
A	280	17484				
A	290		39733			
A	290		17486			
A	286			10981		
A	290				27079	
A	290				17488*	
A	280					27077
A	280					17490*
A	276					10984
B&C	165	23745				
B&C	165	17485*				
B&C	176		23738			
B&C	175		17487*			
S&C	165			10986		
B&C	176				27078	
B&C	175				17489*	
B&C	180					27076
B&C	180					17491*
B&C	186					10989

Evaporating flasks (pear-shaped)

For EL systems only with adapter flange-STJ

STJ	29/32	24/40	34,5/35	29/42	24/29	19/26
50 ml	00431	08750	08740	08736	00472	08743
100 ml	00432	08751	08741	08737	08479	08744
250 ml	00433	08754	08755	08738	08753	08745
500 ml	00434	08758	08759	08739	08757	08746
1000 ml	00435	00440	08763	08762	08761	08747
2000 ml	00436	08765	08766	08769	08764	08748
3000 ml	00437	08767	08742	08770	08735	08749

Adapter flange STJ for EL systems


1 adapter STJ 29,2/32 is supplied with the ROTAVAPOR

STJ	29,2/32	24/40	34,5/35	29,2/42	24/29	19/26
	23733	23747	00921	27101	27100	00916


BÜCHI seal

	Code
KD-22 for RE systems	00636
KD-26 for EL systems	10179


Both seals can be used without grease.


Accessories

Receiving flask
round with SPJ 35/20

Flanged evaporating flasks
for EL systems


Special flasks for drying
powdery substances


Capacity	Code No.	Code No.
500 ml	064404	11579
1000 ml	064405	064420
2000 ml	064406	11580

Capacity	Code No.	Code No.
500 ml	064407	11581
1000 ml	064408	064421
2000 ml	064409	11582


Capacity	Code No.	Code No.
500 ml	064404	11579
1000 ml	064405	064420
2000 ml	064406	11580

Spider evaporators
for RE models

Spider evaporators, for RE models

Cooling coil insert

This coil can be inserted in the cold-trap condenser in place of the cold-finger.


With this accessory, the ROTAVAPOR can now be operated with a cold-trap condenser or with a normal condenser, as desired.

Capacity	Code No.	Code No.
500 ml	064404	11579
1000 ml	064405	064420
2000 ml	064406	11580

Capacity	Code No.	Code No.
500 ml	064407	11581
1000 ml	064408	064421
2000 ml	064409	11582


Capacity	Code No.	Code No.
500 ml	064404	11579
1000 ml	064405	064420
2000 ml	064406	11580

Accessories


Flanged flask for flanged vapour duct, for RE models
The wide flask neck makes it easy to remove substances from the evaporating flask.

Capacity	Code
500 ml	00444
1000 ml	00442
2000 ml	00443
3000 ml	00444


Flanged vapour duct
for RE models

«Soxhlet» extraction unit
for EL models


Code	Item
00509	Top assembly
00510	Bottom assembly and solid/liquid compartment
00520	Unfilled compact compartment
01541	Cleaning seal

Code	Item
11744	Solid/liquid complete
11745	Capacity 200 ml
08569	Capacity 500 ml
00569	Extracting cartridges
00989	Capacity 200 ml
00989	Capacity 500 ml


Vapour temperature measuring sensors

Recommended accessories for the RE-121/EL-131 (for the RE-140/EL-141 included when supplied)

0-1000 mbar vacuummeter
for the RE-111/121 and EL-131

Stand assembly adapter
This adapter enables the drive unit to be mounted on a laboratory stand having a rod of diameter 12 to 15 mm.
Suitable for units 111/121/131

Code	Item
13793	Stand assembly adapter


Code	Item
13793	Stand assembly adapter

Code	Item
13793	Stand assembly adapter

ROTAVAPOR R-152

Introduction


The ROTAVAPOR R-152 has been designed to process batches between 1 and 15 litres of solvent or for continuous feed. The R-152 is predominantly used for the gentle evaporation of solvents at normal pressure or under vacuum. The ROTAVAPOR R-152 is suitable for many applications in addition to the distillation and the stripping of solvents. The design concept of the glass configurations combined with the accessories provided, make it possible to carry out the following operations: powder and granule drying, reflux reactions, recrystallization, solid/liquid extractions, automatic concentration by evaporation, and many other operations.


Technical data

Rotary drive	Motor	Enclosed induction motor with cooling, 100 Watt, protection class IP 54			
	Speed regulation	Electronic, infinitely variable from 15 to 130 rpm, high torque over the entire range			
Bath	Diameter Depth Capacity Bath pan	412 mm 220 mm 29 l Made of stainless steel, water bath with continuous replenishing with needle valve for constant level			
	Heater regulation	Oil bath with protective surround eliminating risk of accidental contact with hot bath wall			
Power rating in Watt	Water bath Temperatur range °C	3750 20–100	Oil bath 20–180		
Over temperature		Cut-out switch which responds whenever the required temperature is exceeded by 15 °C for any reason.			
Raising/Lowering		A warning light is illuminated			
		Motorized, with automatic lowering in the event of a power failure			
Flasks	Receiving flask Evaporation flask	10 l with built-in drain valve for emptying and volume scale 10 l as standard, 6 and 20 l as option			
Display/ Operating elements	Bath temperature Vapour temperature Speed of rotation	with large green LED's, switchable by			
Aeration		Built-in aeration valve, aerates when the aeration key is pressed or in event of a power failure. Also suitable for use with inert gas			
Vacuum	Thightnes	Better than 1 mbar			
Pump	Recommended pump rate	3–5 m ³ /hour			
Cooling water consumption	Approx. 200 l/hour				
Weight	Approx. 70 kg, depending on glass configuration				
Dimensions	Width R configuration	90 cm	Depth 49 cm		
All other configurations	105 cm	58 cm	Height 172 cm		


The six different glass configurations


Rotavapor R-152/R
Reflux


Rotavapor R-152/RW
Reflux with interchangeable receiver


Rotavapor R-152/D
Descending

Rotavapor R-152/DW
Descending with interchangeable receiver

Rotavapor R-152/D2
Descending with two condensers


Rotavapor R-152/D2W
Descending with two condensers and with interchangeable receiver

Ordering Information

Glass configuration/Bath-type	Voltage (Volt) Frequency (Hertz)	3 x 380/220 50-60	1 x 220 50-60	3 x 220 50-60	1 x 240 50-60
Reflux					
Water bath	152/R	27192	27194	27196	27198
Water bath, change overreceiver	152/RW	27210	27211	27213	27212
Oil bath	152/R	27193	27195	27197	27197
Oil bath, change overreceiver	152/RW	27230	27231	27233	27232
Descending					
Water bath	152/D	27214	27215	27217	27216
Water bath, change overreceiver	152/DW	27222	27223	27225	27224
Oil bath	152/D	27234	27235	27237	27236
Oil bath, change overreceiver	152/DW	27242	27243	27245	27244
Descending with two condensers					
Water bath	152/D2	27218	27219	27221	27220
Water bath, change overreceiver	152/D2W	27226	27227	27229	27228
Oil bath	152/D2	27238	27239	27241	27240
Oil bath, change overreceiver	152/D2W	27246	27247	27249	27248


ROTAVAPOR M

ROTAVAPOR M Standard

With a vertical condenser, this ROTAVAPOR is suitable for working under reflux. If the condenser is swivelled downwards, the condensate runs into the receiving flask.

The many accessories make this system particularly suitable for educational purposes. Ask for the guide for carrying out the experiments!

Special accessories (only for ROTAVAPOR M Standard!) for carrying out the following operations:

The ROTAVAPOR, Standard and Compact are designed for working with small volumes. Evaporating flasks with capacities of 30, 50, 100 and 200 ml can be firmly screwed to the drive unit.

The volumes of the receiving flasks are 50, 100 and 250 ml.

The ROTAVAPOR M is supplied with a water bath or oil bath, as desired.

General specifications of ROTAVAPOR M


Drive unit: Sparkless induction motor
Power input 24 W
Rpm fixed at 150

Baths: Ø 140 mm external and Ø 110 mm internal
Bath capacity 1 litre, power input 400 W
Water bath 20–110°C
Oil bath 60–200°C


Weight: 5 kg net, shipping weight 7 kg.

ROTAVAPOR M Compact

With its glass assembly, the ROTAVAPOR M Compact is suitable for applications restricted exclusively to the distillation of solvents. The ROTAVAPOR M Compact has a structure requiring a particularly small amount of space.

ROTAVAPOR 176 Ex/186 Ex

ROTAVAPOR 176 Ex

Rotary evaporator with 20 litre flask for use in rooms with explosion hazard. Complies with SEV and VDE codes. Steam-heated water bath, infinitely variable speed control with mechanical gear unit.


Alternation receiving device with two receiving flasks of 10 litres each.

ROTAVAPOR 186 Ex

Rotary evaporator with 50 litre evaporating flask, explosion-proof in compliance with SEV and VDE codes. Steam-heated water bath adjustable from 30 to 100°C. Infinitely variable mechanical speed control from 5 to 95 rpm. 2 alternative setups as with ROTAVAPOR 176.

Alternation receiving device with two receiving flasks of 20 litres each.

B-168 Vacuum/Distillation Controller

The B-168 Vacuum/Distillation Controller keeps a vacuum constant and can also be used as an apparatus for automatic distillation control of a ROTAVAPOR by the addition of the Dual-Temp probe.

If the Dual-Temp probe for distillation control is not connected, the B-168 performs the functions of a vacuum controller: a preselected pressure in the range of atmospheric pressure to 1 mbar is kept constant by switching a solenoid valve in the vacuum line or by controlling a water jet pump equipped with valves. The band width of the switching points (hysteresis) is matched automatically.

All B-168 units are equipped as standard with an aeration valve.

Additional applications for the B-168 include a wide range of vacuum systems such us drying ovens, desiccators, filtrations and many other equipment requiring a defined, precisely controlled vacuum.

Ordering Information

B-168 Vacuum/Distillation Controller

Voltage	Part Number
230 V	27053
100-115 V	27052
240 V	26564

B-167 water-jet pump unit

Part Number	27054

B-164 24 V water-jet pump unit

Part Number	28348

Dual-Temp probe, complete

Condenser Type	Part Number
For diagonal condenser with SG 119/38	24919
For vertical condenser with SG 929/32 and inside diameter 0.15-22 mm until 200 mm	24910

Manufacturing program

Evaporation of Solvents

- Rotary Evaporators ROTAVAPOR®
- Vacuum Controller
- PTFE pump and controller
- Recirculation Chiller
- Circulating Aspirator
- Heating Baths
- Spray Dryer
- Tube Ovens

Preparative Liquid-Chromatography

- Medium Pressure Liquid Chromatography System
- Fraction Collector

Determination of Nitrogen and Fat

- Kjeldahl Nitrogen and Protein Analysis
- Total Fat Determination System
- Extraction and Gas Scrubbing System

Various Instruments

- Melting Point, Boiling Point and Dropping Point
- Water Distillation Units
- Stereo Models

Laboratoriums-Technik

Buchi Laboratoriums-Technik AG
P.O. Box
CH-9230 Flawil/Switzerland

Telephone 071 848181
Telex 881230 buchi ch
Telefax 071 835711