

Quick Reference Guide

for

Troubleshooting 7673 Autosampler Errors

A16008

This document is believed to be accurate and up-to-date. However, Agilent Technologies, Inc. cannot assume responsibility for the use of this material. The information contained herein is intended for use by informed individuals who can and must determine its fitness for their purpose.

A16008.doc http://www.chem.agilent.com Page 1 of \$

The Quick Reference Guide will help to isolate problems associated with the 7673 Autosampler System to help determine the site of the problem, i.e. in the injector (tower), controller, or 100 vial tray (if present).

Common problems stem from using vials that are not designed to work with the autosampler, syringes that have become "sticky" or "frozen" so that the plunger movement is restricted, or a build-up of film on the surface of the sample turret in the tower. It is recommended that these issues be fully evaluated before moving to in-depth troubleshooting.

If further diagnostics are desired, please see associated documents:

Injector Diagnostics-18593A/B or G1513A

Tray Diagnostics with 18594A/B Controllers

Tray Diagnostics with the G1512A Controller

This document is believed to be accurate and up-to-date. However, Agilent Technologies, Inc. cannot assume responsibility for the use of this material. The information contained herein is intended for use by informed individuals who can and must determine its fitness for their purpose.

LED Error Codes

INJECTOR ERRORS

Injector LedDescriptionRed Steady OnSafety Error

Red 2 Flashes Injector Syringe Error
Red 3 Flashes Injector Turret Error
Red 4 Flashes Injector Plunger Error
Red 5 Flashes Incomplete Injection Error
Red 6 Flashes Injector Memory Error

Injector Front Panel LEDs Controller LED Description

Red, yellow, green flashing in rotation Green On Injector Alignment Error

Yellow "ALIGN" LED on

The yellow "ALIGN" LED is located inside the tower door where the control switches are located.

TRAY ERRORS

Controller LED Description

Red Steady On Tray Not Connected Error
Red 2 Single Flashes* Vial Lost/Left in Gripper Error
Red 3 Single Flashes* Theta Axis Error

Red 4 Single Flashes*
Red 5 Single Flashes*
Red 5 Single Flashes*
Red 6 Single Flashes*
Red 7 Single Flashes*
Red 7 Single Flashes*
Red 7 Single Flashes*
Red 7 Single Flashes*
Red 8 Single Flashes*
Red 9 Single Flashes*
Red 1 Single Flashes*
Red 1 Single Flashes*
Red 2 Single Flashes*
Red 3 Single Flashes*
Red 3 Single Flashes*
Red 3 Single Flashes*
Red 4 Single Flashes*
Red 5 Single Flashes*
Red 5 Single Flashes*
Red 5 Single Flashes*
Red 6 Single Flashes*
Red 7 Single Flashes*

CONTROLLER ERRORS

Controller LED Description

Red, 2 Double Flashes

Red, 3 Double Flashes

Red, 4 Double Flashes

Red, 5 Double Flashes

Controller RAM Error

Controller ROM Error

Configuration Error

MISCELLANEOUS ERRORS

Injector LED Controller LED Description

Red, Yellow, Green On Green On Injector Processor Error
Red, Yellow, Green On Red 18-Volt Supply Error
None None Power Supply Error

^{*}The Controller/GPIB board also uses single red flashes to communicate error conditions found during its power-up selftest. See "Controller/GPIB Board Errors."

Troubleshooting—Automatic Liquid Sampler with G1512A Controller

LED Error Codes

CTOR ERRORS Injector LED	Controller LED	Description
Red On Steady	Green On Steady	Injector Door Open/
		Not Mounted Error
Red, 2 Single Flashes	Green On Steady	Syringe Error
Red, 3 Single Flashes	Green On Steady	Turret Error
Red, 4 Single Flashes	Green On Steady	Plunger Error
Red, 5 Single Flashes	Green On Steady	Incomplete Injection Error
Red, 6 Single Flashes	Green On Steady	Injector Memory Error
None	None	Power Error

TRAY ERRORS Injector LED	Controller LED	Description
Green On Steady	Yellow, 3 Single Flashes/ Green On Steady/3 Beeps	Bottle In Gripper Error
Green On Steady	Yellow, 4 Single Flashes/ Green On Steady/4 Beeps	Hard Tray/Injector Error

SCELLANEOUS LED PATTERNS Injector LED	Controller LED
Red, Yellow, Green On Steady	None
Red, Yellow, Green On Steady	Green On Steady
Red, Yellow, Green On Steady	Yellow, Green On Steady
None	Yellow, 1 Single Flash/1 Beep
None	Yellow, 2 Single Flashes/2 Beeps
None	Yellow, 3 Single Flashes/3 Beeps
None	Green, 1 Single Flash
None	Yellow, Green, 1 Alternating Flash/Continuous Beep

Troubleshooting—Bar Code Reader with 18594A/B Controller

Error Condition Beeps

If an error condition exists, there will be a series of beeps at power-on to draw your attention, a pause of approximately half a second, and the series of beeps will then be repeated once.

- 1. Count the beeps.
- 2. Decode the error message (see below).

Beeps	Indication	Action	
10 beeps	ROM error	Replace BCR board	
11 beeps	RAM error	Replace BCR board	
13 beeps	Reader cable is disconnected or damaged	Connect cable or replace	

Configuration Beeps

Following the error tests there will be a series of beeps to indicate the configuration of the switch settings. This series will not be repeated.

- 1. Count the beeps.
- 2. Decode the configuration message (see below).

Beeps	Indication	Beeps	Indication
1 beep	board is set in TTL mode	10 beeps	read mixed, zero suppressed,
2 beeps	board is set in command-only mode		checksum off
3 beeps	read Code 39 with checksum off	11 beeps	read mixed, zero suppressed,
	—Normal mode		checksum on
4 beeps	read Code 39 with checksum on	12 beeps	read mixed, checksum off
5 beeps	read Code 2 of 5 with checksum off	13 beeps	read mixed, checksum off
6 beeps	read Code 2 of 5 with checksum on	14 beeps	read mixed, checksum off
7 beeps	read UPC code	15 beeps	read mixed, checksum off
8 beeps	read mixed, checksum off	16 beeps	read mixed, checksum off
9 beeps	read mixed, checksum on		