IKA

designed for scientists

IKA ETS-D5

Betriebsanleitung	DE	4
Ursprungssprache		
Operating instructions	EN	11
Mode d'emploi	FR	18
Instrucciones de manejo	ES	25
Instruções de serviço	PT	32
使用说明	ZH	39
Veiligheidsaanwijzingen	NL	45
Avvertenze per la sicurezza	IT	46
Säkerhetsanvisningar	SV	48
Sikkerhedshenvisninger	DA	49
Sikkerhetsveiledning	NO	51
Turvallisuusohjeet	FI	52
Wskazówki bezpieczeństwa	PL	54
Bezpečnostní pokyny a informace	CS	55
Biztonsági óvintézkedések	HU	57
Varnostní napotki	SL	58
Bbezpečnostné pokyny	SK	60
Ohutusjuhised	ET	61
Drošības norādījumi	LV	63
Saugos nurodymai	LT	64
Указания за безопасност	BG	66


Fig. 1

Fig. 2


Fig. 3

Fig. 4

Fig. 5

DE	
A	Messfühler
B	Abdeckkappe
C	LCD-Display
D	Taste „▲“
E	Taste „▼“

EN	
A	Sensor
B	Cap
C	LCD display
D	Button "▲"
E	Button "▼"

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Safety instructions

For your protection

- Read the operating instructions in its entirety before using the device and follow the safety instructions. If this accessory is used with another device, observe also its operating instructions.
- Keep the operating instructions in a place where it can be accessed by everyone.
- Ensure that only trained staff work with the device.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- To avoid body injury and property damage, observe the relevant safety and accident prevention measures when processing hazardous materials.
- Wear your personal protective equipment in accordance with the hazard category of the media to be processed.
- **Caution!** The supply voltage (Limited Power Supply) for the **ETS-D5** must be provided using an IEC 61558 compliant transformer with double or reinforced insulation. Alternatively, laboratory equipment that meets the requirements of IEC 61010 or NEC Class 2 may be used.
- **Caution!** In unsupervised and safe operation, this equipment must only be used with media that have a flash point higher than the temperature safety limit set for the magnet agitator and/or heating system used. The safety temperature limit must always be set to at least 25 °C lower than the flash point of the media used. (acc. to EN 61010-2-010)
- **Burn hazard!** Take care when touching the temperature sensor! The temperature sensor may heat up to 450 °C when inserted in the media!
- Do not touch the temperature sensor while measurements are being taken. This will prevent incorrect results.
- Ensure that the external temperature sensor is inserted into the medium to a depth of at least 20 mm when connected.
- Please ensure that the spiral cable does not come into contact with the heating plate.

Contents

Page

EU Declaration of conformity	11
Safety instructions	11
Intended use	12
Unpacking	12
Factory settings	12
Commissioning	14
Operating modes	14
Setting the HI-TEMP	15
Adjusting the set-point temperature	15
Replacing extension cable and sensor	15
Error codes	16
Accessories	16
Technical Data	17
Warranty	17

EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61326-1, EN 60529 and EN ISO 12100.

A copy of the complete EU Declaration of Conformity can be requested at sales@ika.com.

- The stainless steel temperature sensor must not be used with aggressive media such as acids, caustic solutions or distilled water, due to the risk of corrosion. The H 66 glass sensor should be used in such cases.
- Only use glass encapsulated temperature sensors for electrolysis procedures.
- Always use the extension cable H 70 when the media being processed produces vapour. This ensures that the control unit does not come into contact with the vapour.
- Beware of hazards due to:
 - flammable materials
 - combustible media with a low boiling temperature
- Do not use the device in explosive atmospheres, it is not EX-protected.
- Danger!** The user must check that all conditions have been met for the safe and correct operation of the equipment before operating the ETS-D5 with a heater that does have a DIN 12878 compliant connector, and therefore does not provide cut-out protection. [The cut-out protection is a mechanism that ensures that the electrical heating circuit is switched off should the contact thermometer (ETS-D5) be broken or if the connector is removed.]

For protection of the equipment

- Protect the device and accessories from bumps and impacts.

Use of "H+P" magnetic stirring systems (no break fuse protection included)

- When the supply current for the ETS-D5 is interrupted, the electronics for the heatable magnetic stirrer system will take over the temperature control. This means that the magnetic stirrer system will heat up to the temperature set (e.g. 300 °C). Readjust the set-point temperature for the magnetic stirrer system or switch the heating off.

Use of "Corning" magnetic stirrer PC-400, PC-420, PC-600, PC-620

- When using the **ETS-D5** in conjunction with a "Corning" Magnetic Stirrer it is absolutely essential that the rotary "heat" dial on the magnetic stirrer is set to the "O" (off) position. The ETS-D5 will only control the set temperature when it is in this position.
- Warning!** The heating plate temperature limiter is not active in this operating mode.
- If the rotary "heat" dial on the magnetic stirrer is not in the "O" position then the magnetic stirrer will heat up to the temperature set, regardless of the settings on the ETS-D5.


Warning! ETS-D5 is inactive.

The safety advices given above represent the current status of the known risks. Nonetheless, the user must check equipment for correct and safe operation before operating any system comprising the ETS-D5 and a magnetic stirring/heating system.

Intended use

Use

- The ETS-D5 is an instrument for the precise control of temperature. The sensor on this instrument must be inserted at least 20 mm into the medium to be tempered. The instrument can be connected to any magnetic stirrer or heating plate that has a contact thermometer connection, provided this conforms with the requirements listed in the technical data (see section Technical Data).

Area of use

- Indoor environments similar to that a laboratory of research, teaching, trade or industry area.
- The safety of the user cannot be guaranteed:
 - if the device is operated with accessories that are not supplied or recommended by the manufacturer,
 - if the device is operated improperly or contrary to the manufacturer's specifications,
 - if the device or the printed circuit board are modified by third parties.

Unpacking

Unpacking

- Unpack the device carefully. Any damage should immediately be reported to the carrier (mail, rail or freight forwarding company).

Scope of delivery

- IKA ETS-D5
- Operating instructions

Factory settings

The ETS-D5 is supplied with the following factory settings

- Operating mode: A
- Set-point temperature: - 10 °C
- HI-TEMP: 450 °C

Restoring parameter to factory settings (RESET)

The following steps must be followed to restore factory settings:

- Switch off the equipment
- Press and hold down the membrane keys "▲" and "▼" and switch on the equipment
- The factory settings will now be restored (see above)

Commissioning

Plug the instrument into the contact thermometer jack on the magnetic stirrer (e.g. IKA RCT basic, IKA RET basic, RH basic, RH digital KT/C, Heidolph MR 3001K, Heidolph MR 3002, VWR VMS-D, VWR VMS-A, H+P Variomag, etc.). Please observe all safety instructions in this respect.

Procedure for connecting a Corning Magnetic Stirrer:
Turn the rotary "heat" dial on the magnetic stirrer to the "O" (off) position. Next, connect the ETS-D5 to the magnetic stirrer using an accessory AD-C1 adapter.
Please observe all safety instructions in this respect.
The device will be activated when the magnetic stirrer is switched on.

If the supply voltage at the contact thermometer connector is not between 8 V and 16 V DC then it will be necessary to power the device using the H 52 power supply set (H 50 double connector adapter and H 51 plug-in power supply unit). The power supply and/or double connector adapter will provide the necessary supply voltage for the device. An automatic self-test will be performed when the device is activated. All segments on the LCD display will be visible during this test (Fig. 1).

Operating modes

The ETS-D5 can be run in three different operating modes (A, B, C).

The selected operating mode will be permanently displayed (Fig. 4).

Operating mode "A"

This is the factory-setting mode for this instrument. When the instrument is switched on, the set-point temperature is always set to - 10 °C. The set-point temperature can be adjusted using the "▲" and "▼" keys. The HI-TEMP (maximum set-point temperature) can only be set in this mode.

Operating mode "B"

In this mode the most recent set-point temperature is used, and will become active when the instrument is switched off and on again. The settings can be adjusted using the "▲" and "▼" keys.

Operating mode "C"

In this mode the current settings are used and will become active when the instrument is switched off and on again. The setting cannot be adjusted.

Holding down the "▲" key when switching on the instrument switches to the next operating mode in the sequence A-B-C-A-B-C-A etc.

Setting the HI-TEMP

The "HI-TEMP" (maximum set-point temperature) can only be set when in operating mode "A". After the self-test, HI-TEMP will be visible on the LCD display (Fig. 2). The HI-TEMP can now be set between 0 and 450 °C using the "▲" and "▼" keys. The displayed HI-TEMP value will be set and stored if no key is pressed for approximately 5 seconds. The LCD display will then appear as shown in Fig. 3.


Adjusting the set-point temperature

The set-point temperature can be adjusted in 0.1 K steps by pressing the "▲" and "▼" keys. If the key is held down for longer than 5 seconds then the set-point temperature will change in 1K or 10 K steps.

The LCD display (Fig. 4) shows: the instantaneous actual temperature with physical units; the set-point temperature (underneath actual temperature);  to indicate that the heating for the connected device is switched on; the selected operating mode (top left).

Replacing extension cable and sensor

To replace the sensor or the extension cable, pull the sensor and protective cap downwards until the plastic snap fastener is released. Then push the sensor or extension cable and protective cap over the fitting on the instrument to re-establish the connection.

Error codes

- The fault is shown by an error code on the display as following if the error occurs. (Fig. 5)

- Proceed as follows in such cases:
 - Disconnect the power supply for at least 10 sec.
 - Carry out corrective measures.
 - Restart the device.

Error code	Cause	Effect	Solution
Er 1	Sensor not connected or defective (Sensor rupture)	Device off	<ul style="list-style-type: none"> - Connect temperature sensor - Replace temperature sensor
Er 5	Sensor not in media (no change in temperature of the media after heater active for 3 minutes)	Device off	<ul style="list-style-type: none"> - Insert the temperature sensor in the media to a depth of at least 20 mm
Er 6	Temperature sensor short-circuit	Device off	<ul style="list-style-type: none"> - Switch off the device and replace temperature sensor
Er 7	The temperature of the medium is higher than the maximum target temperature that can be set	Device off	<ul style="list-style-type: none"> - Switch off the device and allow the media to cool down - Increase the maximum target temperature that can be set
Er 8	Calibration error	Device off	<ul style="list-style-type: none"> - Switch off the device - Re-calibration by authorised service personnel

- If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department,
- Send the device for repair, including a short description of the fault.

Accessories

- For accessories see www.ika.com.

Technical Data

Temperature measuring range	°C	- 50 ... 450
Temperature measurement resolution	K	0.1
Accuracy of temperature measurement	K	± 0.2 + sensor tolerance PT 1000 DIN IEC 751 Class A
Set temperature resolution	K	0.1
Type of control		PID
Control deviation (500 ml water in 600 ml glass beaker, magnetic stirring bar 40, at 600 rpm, 50 °C)	K	± 0.5
Immersion depth max.	mm	200
Supply voltage	V DC	8 ... 16 (Limited Power Supply)
Current consumption	mA	15
Permissible ambient temperature	°C	0 ... 60
Permissible relative humidity	%	80
Permissible duration of operation	%	100
Protection class according to DIN EN 60529		IP 54
Cable length	mm	350
Male connector		6-pin DIN 45322
Socket		DIN 12878 Class 2
Dimensions (W x D x H)	mm	82 x 22 x 83 (without sensor)
Weight	kg	0.2

Subject to technical changes!

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

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IKA Vietnam Company Limited
IKA India Private Limited