


Instruction Manual

HI 83224

Multiparameter Bench Photometer



www.hannainst.com

Dear Customer,

Thank you for choosing a Hanna Instruments product. Please read this instruction manual carefully before using this instrument. This manual will provide you with the necessary information for correct use of this instrument, as well as a precise idea of its versatility. If you need additional technical information, do not hesitate to e-mail us at tech@hannainst.com or view our worldwide contact list at www.hannainst.com.

TABLE OF CONTENTS

PRELIMINARY EXAMINATION	3
ABBREVIATIONS	3
GENERAL DESCRIPTION	3
SPECIFICATIONS	4
PRECISION AND ACCURACY	4
PRINCIPLE OF OPERATION	4
FUNCTIONAL DESCRIPTION	6
TIPS FOR AN ACCURATE MEASUREMENT	8
HEALTH & SAFETY	12
METHOD REFERENCE TABLE	13
OPERATIONAL GUIDE	15
SETUP	19
HELP MODE	21
AMMONIA LOW RANGE	22
AMMONIA HIGH RANGE	24
CHLORINE, FREE	26
CHLORINE, TOTAL	28
NITRATE	30
NITROGEN, TOTAL LR	32
NITROGEN, TOTAL HR	37
OXYGEN DEMAND, CHEMICAL LR	42
OXYGEN DEMAND, CHEMICAL MR	45
OXYGEN DEMAND, CHEMICAL HR	48
PHOSPHORUS, REACTIVE	51
PHOSPHORUS, ACID HYDROLYZABLE	53
PHOSPHORUS, TOTAL LR	57
PHOSPHORUS, REACTIVE HR	61
PHOSPHORUS, TOTAL HR	64
ERRORS AND WARNINGS	68
DATA MANAGEMENT	69
STANDARD METHODS	69
ACCESSORIES	70
WARRANTY	71
HANNA LITERATURE	71

All rights are reserved. Reproduction in whole or in part is prohibited without the written consent of the copyright owner, Hanna Instruments Inc., Woonsocket, Rhode Island, 02895, USA.

PRELIMINARY EXAMINATION

Please examine this product carefully. Make sure that the instrument is not damaged. If any damage occurred during shipment, please notify your Local Hanna Office.

Each Meter is supplied complete with:

- Ten Sample Vials
- Cloth for wiping vials (1 pcs.)
- Scissors
- Instruction Manual

Note: Save all packing material until you are sure that the instrument works correctly. Any defective item must be returned in its original packing with the supplied accessories.

ABBREVIATIONS

°C:	degree Celsius
°F:	degree Fahrenheit
COD:	Chemical Oxygen Demand
EPA:	US Environmental Protection Agency
µg/L:	micrograms per liter (ppb)
mg/L:	milligrams per liter (ppm)
g/L:	grams per liter (ppt)
mL:	milliliter
LR:	low range
MR:	medium range
HR:	high range

GENERAL DESCRIPTION

HI 83224 is a multiparameter bench photometer. It measures 15 colorimetry based methods. Furthermore, it can identify the samples via a bar code label placed on the vials.

The reagents are in liquid or powder form and are supplied in bottles, ready-to-use vials or in packets. The amount of reagent is precisely dosed to ensure maximum repeatability.

HI 83224 bench photometer can be connected to a PC via an USB cable. Its software companion is **HI 92000** Windows® Compatible Software that helps the user to manage all the results.


HI 83224 has a very important feature. The samples are identified via bar code labeled vials. The corresponding bar codes for different types of analysis are shown in a table (see page 14). For those methods that don't have predosed reagents the vials supplied with the instrument should be used. The bar code has 2 fields of 2 digits. The first field is for method identification and the other one is for reagent lot identification.

HI 83224 has a powerful interactive user support that assists the user during the analysis process.

Each step in the method is help supported on the LCD. Help functions are accessed easily by pressing the **HELP** key. A tutorial mode is also available.

FUNCTIONAL DESCRIPTION










INSTRUMENT DESCRIPTION



- 1) Vial holder
- 2) Liquid Crystal Display (LCD)
- 3) Keypad
- 4) Bottom lid
- 5) USB connector
- 6) Fuse holder
- 7) ON/OFF power switch
- 8) Power supply connector

KEYPAD DESCRIPTION

The keypad contains 8 direct keys and 3 functional keys with the following functions:

	Press to perform the function displayed above it on the LCD.
	Press to exit the current screen.
	Press to access the method menu.
	Press to move up in a menu or a help screen, to increment a set value, or to access second level functions.
	Press to move down in a menu or a help screen, to decrement a set value, or to access second level functions.
	Press to identify the bar code on the vial.
	Press to recall the log.
	Press to display the help screen.
	Press to access the setup screen.

TIPS FOR AN ACCURATE MEASUREMENT

The instructions listed below should be carefully followed during testing to ensure most accurate results.

COLLECTING AND MEASURING SAMPLES

- For adding the exact amount of sample or liquid reagent to the reagent vials it is strongly recommended to use the available Hanna automatic fixed-volume pipettes or class A laboratory pipettes (symbolized like a generic pipette tip in the following method related chapters):

pipette code	volume
HI 731340	200 μ L
HI 731341	1000 μ L
HI 731342	2000 μ L




For correct use of the automatic pipette, please follow the related Instruction Sheet.


Alternatively, the optional **HI 740142** 1 mL graduated syringe or **HI 740226** 5 mL graduated syringe can be used.

For the correct use of the syringes, see instructions below.

- In order to measure exactly 5 mL of reagent with the 5 mL syringe:
 - push the plunger completely into the syringe and insert the tip into the solution
 - pull the plunger up until the lower edge of the seal is exactly on the 5 mL mark
 - take out the syringe and clean the outside of the syringe tip. Be sure that no drops are hanging on the tip of the syringe, if so eliminate them. Then, keeping the syringe in vertical position above the vial, push the plunger completely down into the syringe. Exactly 5 mL has been added to the vial



2. In order to measure exactly 0.5 mL of reagent with the 1 mL syringe:
- push the plunger completely into the syringe and insert the tip into the solution
 - pull the plunger up until the lower edge of the seal is exactly on the 0.0 mL mark
 - take out the syringe and clean the outside of the syringe tip. Be sure that no drops are hanging on the tip of the syringe, if so eliminate them. Then, keeping the syringe in vertical position above the vial, push the plunger down into the syringe until the lower edge of the seal is exactly on the 0.5 mL mark. Now the exact amount of 0.5 mL has been added to the vial, even if the tip still contains some solution




- Color or suspended matter in large amounts may cause interference, therefore, these should be removed by treatment with active carbon and by prior filtration.

USING LIQUID AND POWDER REAGENTS

- Proper use of the dropper:
 - to get good reproducible results, tap the dropper on the table for several times and wipe the outside of the dropper tip with a cloth
 - always keep the dropper bottle in a vertical position while dosing the reagent




- Proper use of the powder reagent packet:
 - (a) use scissors to open the powder packet
 - (b) push the edges of the packet to form a spout
 - (c) pour out the content of the packet




USING VIALS


- Never insert hot vials in the instrument, the vial holder may become damaged.
- In order to avoid reagent leaking and to obtain most accurate results, it is recommended to close the vial tightly with the supplied cap after addition of reagents or sample.
- Shaking the vial can generate bubbles in the sample, causing higher readings. To obtain accurate measurements, remove such bubbles by swirling or by gently tapping the vial.
- Proper mixing is very important for reproducibility of the measurements. The right way of mixing a vial is specified for each method in the related chapter.
 - (a) invert the vial a couple of times or for a specified time: hold the vial in the vertical position with the cap up. Turn the vial upside-down and wait for all of the solution to flow to the cap end, then return the vial to the upright vertical position and wait for all of the solution to flow to the vial bottom. This is one inversion. The correct speed for this mixing technique is 10 complete inversions in 30 seconds.




This mixing technique is indicated with "invert to mix" and the following icon:



- (b) **shake the vial**: move the vial up and down. The movement may be gentle or vigorous.
This mixing technique is indicated with “shake gently” or “shake vigorously” and one of the following icons:



Pay attention to push the vial completely down in the holder.



- Whenever the vial is placed into the measurement cell, it must be dry on the outside, and completely free of fingerprints, oil or dirt. Wipe it thoroughly with **HI 731318** or a lint-free cloth prior to insertion.
- Do not let the reacted sample stand too long after reagent is added. For best accuracy, respect the timings described in each specific method.
- It is possible to take multiple readings in a row, but it is recommended to take a new zero reading for each sample and to use the same vial for zeroing and measurement where possible (for most precise results follow the measurement procedures carefully).
- All the reaction times reported in this manual are at 25 °C (77 °F). In general, the reaction time should be increased for temperatures lower than 20 °C (68 °F), and decreased for temperatures higher than 25 °C (77 °F).



DIGESTION

- Some analytical methods require digestion of the sample. Use the **HI 839800 test tube heater** for digestion. The optional **HI 740217 safety shield** is strongly recommended. For correct use of the reactor follow the Reactor Instruction Manual. At the end of the digestion period the vials are still hot, allow the vials to cool to room temperature in the optional **HI 740216 test tube cooling rack**.

REAGENT BLANK CORRECTION

- Some methods require a “reagent blank correction”. The blank and the sample are prepared exactly in the same way, with the only difference that for the blank deionized water is used instead of sample.
A blank vial may be used more than once: stability and storing conditions are described for each method in the related chapter.

INTERFERENCES

- In the method related measurement sections we have reported the most common interferences that may be present in an average wastewater matrix. It may be that for a particular treatment process other compounds do interfere with the method of analysis.

HEALTH & SAFETY



- The chemicals contained in the reagent kits may be hazardous if improperly handled.
- Read the Material Safety Data Sheet (MSDS) before performing tests.
- Safety equipment: Wear suitable eye protection and clothing when required, and follow instructions carefully.
- Reagent spills: If a reagent spill occurs, wipe up immediately and rinse with plenty of water.
If reagent contacts skin, rinse the affected area thoroughly with water. Avoid breathing released vapors.
- Reagent vial disposal: Reagents vials may contain different waste pollutants. After use dispose of the reagent vials according to the local regulations.

METHOD REFERENCE TABLE

Method reagent set code	Method description	Method	Page
1	HI 94764A-25	Ammonia LR	22
2	HI 94764B-25	Ammonia HR	24
3	HI 93701-01 HI 93701-03	Chlorine, Free	26
4	HI 93711-01 HI 93711-03	Chlorine, Total	28
5	HI 94766-50	Nitrate	30
6	HI 94767A-50	Nitrogen, Total LR	32
7	HI 94767B-50	Nitrogen, Total HR	37
8	HI 94754A-25	Oxygen Demand, Chemical (COD) LR	42
9	HI 94754B-25	Oxygen Demand, Chemical (COD) MR	45
10	HI 94754C-25	Oxygen Demand, Chemical (COD) HR	48
11	HI 94758A-50	Phosphorus, Reactive	51
12	HI 94758B-50	Phosphorus, Acid Hydrolyzable	53
13	HI 94758C-50	Phosphorus, Total LR	57
14	HI 94763A-50	Phosphorus, Reactive HR	61
15	HI 94763B-50	Phosphorus, Total HR	64

VIAL IDENTIFICATION

Predosed vials are related to the different methods and can be identified by a bar code printed on the vial label (all methods except Free and Total Chlorine):



Method reagent set code	Method description	Vial bar code
HI 94764A-25	Ammonia LR	01xx
HI 94764B-25	Ammonia HR	02xx
HI 93701-01, HI 93701-03	Chlorine, Free	-
HI 93711-01, HI 93711-03	Chlorine, Total	-
HI 94766-50	Nitrate	05xx
HI 94767A-50	Nitrogen, Total LR	06xx
HI 94767B-50	Nitrogen, Total HR	07xx
HI 94754A-25	Oxygen Demand, Chemical (COD) LR	12xx
HI 94754B-25	Oxygen Demand, Chemical (COD) MR	13xx
HI 94754C-25	Oxygen Demand, Chemical (COD) HR	24xx
HI 94758A-50	Phosphorus, Reactive	30xx
HI 94758B-50	Phosphorus, Acid Hydrolyzable	31xx
HI 94758C-50	Phosphorus, Total LR	32xx
HI 94763A-50	Phosphorus, Reactive HR	33xx
HI 94763B-50	Phosphorus, Total HR	34xx

Note: xx represents the reagent lot code.

OPERATIONAL GUIDE

POWER CONNECTION

Connect the instrument to the AC power supply via the power cord that is supplied with the instrument.

Note: Ensure the main line is surge protected.




Note: Always turn the meter off before unplugging it to ensure no data is lost.

METHOD SELECTION


The instrument can use three operating modes in order to perform measurements: Automatic, SemiAutomatic and Manual. Select the desired operating mode in Setup. The operating mode is related to how the desired method is selected.

Automatic

- Turn the instrument ON via the ON/OFF power switch.
- The meter will perform an autodiagnostic test. During this test, the Hanna Instrument logo will appear on the LCD. After 5 seconds, if the test was successful, the last selected method will appear on the display.
- Take a vial from the set of vials for the method to be measured. Insert the vial into the meter and press the **Zero** key. The instrument will check if the vial corresponds to the current method and will perform the zero sequence. If the meter detects that the vial is for another method, it will automatically change the method and will perform the zero sequence.
- After the desired method was selected, follow the measurement procedure described in the related chapter.




- Before performing a test carefully read all the instructions related to the selected method.
- If the method identification wasn't performed successfully the following error messages will appear:




SemiAutomatic

- Turn the instrument ON via the ON/OFF power switch.
- The meter will perform an autodiagnostic test. During this test, the Hanna Instrument logo will appear on the LCD. After 5 seconds, if the test was successful, the last selected method will appear on the display.

- Take a vial from the set of vials for the method to be measured. Insert the vial into the meter and press the **Zero** key. The instrument will check if the vial corresponds to the current method and will perform the zero sequence. If the meter detects that the vial is for another method, it will display a warning message which allows the user to change the method to the newly detected one or keep the existing method.
- After the desired method was selected, follow the measurement procedure described in the related chapter.




- Before performing a test carefully read all the instructions related to the selected method.
- If the method identification wasn't performed successfully the following error message will appear:




Manual

- Turn the instrument ON via the ON/OFF power switch.
- The meter will perform an autodiagnostic test. During this test, the Hanna Instrument logo will appear on the LCD. After 5 seconds, if the test was successful, the last method used will appear on the display.
- In order to select the desired method press the **METHOD** key and the screen with the list of methods will appear.
- Press the **▲ ▼** keys to highlight the desired method. Press **Select**.



- Take a vial from the set of vials for the method to be measured. Insert the vial into the meter and press the **Zero** key. The instrument will check if the vial corresponds to the selected method and will perform the zero sequence. If the meter detects that the vial is for another method, it will display a warning message. The measurement can continue with the selected method by pressing the **Continue** key or to be stopped by pressing the **Cancel** key.
- After the desired method is selected, follow the measurement procedure described in the related chapter.




- Before performing a test carefully read all the instructions related to the selected method.


DATA MANAGEMENT

You can keep track of your results using the instrument's data log function. Up to 200 individual measurements can be stored in the data log. Storing, viewing and deleting the data is possible using ▼ or ▲ and RCL keys.

Storing data: You can store only a valid measurement. To store a valid measurement press ▼ or ▲ to access the second level functions and then press **Log**. The measurement will be stored with date and time stamps.



Viewing and deleting: You can view and delete the data log by pressing the RCL key. You can delete the last saved measurement. Additionally, you can delete the data records all at once.



CHEMICAL FORM

Chemical form conversion factors are pre-programmed into the instrument and are method specific. In order to view the displayed result in the desired chemical form press ▲ or ▼ to access the second level functions and then press the **Chem Frm** key to toggle between the available chemical forms for the selected method.




METHOD CHECK



Methods associated with one vial can be verified using the Method Check feature.

This feature can be used at anytime in order to avoid mistakes during the measurement procedure.


- Place the vial into the holder and push it completely down.



- Press **METHOD CHECK** to start the identification procedure.
- If the method identification isn't performed successfully, the following warning messages will appear:



- If the method is successfully identified, the name is displayed:



SETUP

In the Setup mode the instrument's parameters can be changed. Some parameters affect the measuring sequence and others are general parameters that change the behaviour or appearance of the instrument.

Press **SETUP** to enter the setup mode.

Press **ESC** or **SETUP** to return to the main screen.

A list of setup parameters will be displayed with currently configured settings. Press **HELP** for additional information.

Press the **▲ ▼** keys to select a parameter and depending to the parameter type, select the new value as follows:



Operating Mode

The instrument has 3 operating modes.

In Automatic mode, no user involvement is required in the method selection.

In SemiAutomatic mode, user involvement is required when a different vial is identified related to the existing method.

In Manual mode, user involvement is required at every level in the method selection.

Press the functional key to select the desired operating mode.





Backlight

Values: 0 to 8.

Press the **Modify** key to access the backlight value.

Use the **◀ ▶** functional keys or the **▲ ▼** keys to increase or decrease the value.

Press the **Accept** key to confirm or **ESC** to return to the setup menu without saving the new value.



Contrast



Values: 0 to 20.

This option is used to set the display's contrast.

Press the **Modify** key to change the display's contrast.

Use the **◀ ▶** functional keys or the **▲ ▼** keys to increase or decrease the value.

Press the **Accept** key to confirm the value or **ESC** to return to the setup menu without saving the new value.




Average

Option: Enable or Disable.

This option is used to enable/disable averaged measuring mode. If enabled, the instrument takes 180 readings and displays the resulting average value. The partial average is displayed during measurement.

Press the functional key to enable or disable this option.





Date / Time

This option is used to set the instrument's date and time.

Press the **Modify** key to change the date/time.



Press the ◀ ▶ functional keys to highlight the value to be modified (year, month, day, hour, minute or second). Use the ▲ ▼ keys to change the value. Press the **Accept** key to confirm or **ESC** to return to the setup without saving the new date or time.



Time format

Option: AM/PM or 24 hour.

Press the functional key to select the desired time format.




Date format

Press the **Modify** key to change the Date Format.

Use the ▲ ▼ keys to select the desired format.

Press the **Accept** key to confirm or **ESC** to return to the setup menu without saving the new format.



Language

Press the corresponding key to change the language.

If the new language cannot be loaded, the previously selected language will be reloaded.





Tutorial

Option: Enable or Disable.

If enabled this option will provide the user short guides related to the current issue, on the screen.

Press the functional key to enable/disable the tutorial mode.




Beeper

Option: Enable or Disable.

When enabled, a short beep is heard every time a key is pressed.

A long beep alert sounds when the pressed key is not active or an error condition is detected.

Press the functional key to enable/disable the beeper.




Instrument ID

Option: 0 to 9999.

This option is used to set the instrument's ID (identification number). The instrument ID is used while exchanging data with a PC.

Press the **Modify** key to access the instrument ID screen. Press the ▲ ▼ keys in order to set the desired value.


Press the **Accept** key to confirm the value or **ESC** to return to the setup menu without saving the new value.



Meter information

Press the **Select** key to view the instrument model, firmware version, language version and instrument serial number.

Press **ESC** to return to the Setup mode.



HELP MODE

HI 83224 offers an interactive contextual help mode that assists the user at any moment.

To access the help screens press **HELP**.




The instrument will display additional information related to the current screen. To read all the available information, scroll the text using the ▲ ▼ keys.

Press the **Support** key to access a screen with Hanna service centers and their contact details.

Press the **Accessories** key to access a list of instrument reagents and accessories.

To exit support or accessories screens press **ESC** and the instrument will return to the previous help screen.

To exit help mode press **HELP** key again and the meter will return to the previously selected screen.



AMMONIA LOW RANGE

SPECIFICATIONS

Range	0.00 to 3.00 mg/L (as $\text{NH}_3\text{-N}$)
Resolution	0.01 mg/L
Accuracy	± 0.10 mg/L or $\pm 5\%$ of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 420 nm
Method	Adaptation of the <i>ASTM Manual of Water and Environmental Technology</i> , D1426-92, Nessler method. The reaction between ammonia and reagents causes a yellow tint in the sample.

REQUIRED REAGENTS

Code	Description	Q.ty/test	Q.ty/set
HI 94764A-0*	Reagent Vial	1 vial	25 vials
HI 93764-0	Nessler Reagent	4 drops (0.25 mL)	1 bottle

* *Reagent Vial identification: 01xx (xx represents the reagent lot code).*

Note: Store the unused vials in a cool and dark place.



REAGENT SET

HI 94764A-25 Reagents for 25 tests



For other accessories see page 70.

MEASUREMENT PROCEDURE





- Select the **Ammonia LR** method following one of the procedures described in the "Method Selection" section (see page 15).
- Remove the cap from Reagent Vial.
- Add exactly 5.0 mL of sample to the vial, while keeping the vial at a 45-degree angle.
- Replace the cap and mix by inverting the vial a couple of times. This is the blank.
- Place the vial into the holder and push it completely down.




- Press **Zero** and wait for vial identification. If that was successfully done, the display will show “-0.0-” when the meter is zeroed and ready for measurement.



- Remove the vial.
- Remove the cap and add 4 drops of HI 93764-0 Nessler Reagent.
- Replace the cap tightly and mix by inverting the vial a couple of times. This is the sample.
- Place the vial into the holder and push it completely down.
- Press **Timer** and the display will show the countdown prior to the measurement and the message “Reaction Time”. Alternatively, wait for 3 minutes and 30 seconds and press **Read**. Wait for vial identification. If that was successfully done, the instrument will perform the reading.



- The instrument displays concentration in **mg/L of ammonia nitrogen (NH₃-N)**. Press **▲** or **▼** to access the second level functions and then press the **Chem Frm** key to convert the result in mg/L ammonia (NH₃).



- Press **▲** or **▼** to return to the measurement screen.

INTERFERENCES

Interference may be caused by:

- Organic compounds like: chloramines, various aliphatic and aromatic amines, glycine or urea above 10 ppm. To eliminate these interferences distillation is required.
- Organic compounds like: aldehydes, alcohols (e.g. ethanol) or acetone above 0.1 %. To eliminate these interferences distillation is required.
- Sulfide: may cause turbidity.

AMMONIA HIGH RANGE

SPECIFICATIONS

Range	0 to 100 mg/L (as $\text{NH}_3\text{-N}$)
Resolution	1 mg/L
Accuracy	± 1 mg/L or ± 5 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 420 nm
Method	Adaptation of the <i>ASTM Manual of Water and Environmental Technology</i> , D1426-92, Nessler method. The reaction between ammonia and reagents causes a yellow tint in the sample.

REQUIRED REAGENTS

Code	Description	Q.ty/test	Q.ty/set
HI 94764B-0*	Reagent Vial	1 vial	25 vials
HI 93764-0	Nessler Reagent	4 drops (0.25 mL)	1 bottle

* *Reagent Vial identification: 02xx (xx represents the reagent lot code). GREEN RECTANGLE ON THE LABEL.*

Note: Store the unused vials in their container in a cool and dark place.



REAGENT SET

HI 94764B-25 Reagents for 25 tests



For other accessories see page 70.

MEASUREMENT PROCEDURE



- Select **Ammonia HR** method following one of the procedures described in the "Method Selection" section (see page 15).
- Remove the cap from Reagent Vial.
- Add exactly 1.0 mL of sample to the vial, while keeping the vial at a 45-degree angle.
- Replace the cap and mix by inverting the vial a couple of times. This is the blank.
- Place the vial into the holder and push it completely down.





- Press **Zero** and wait for vial identification. If that was successfully done, the display will show “-0.0-” when the meter is zeroed and ready for measurement.



- Remove the vial.
- Remove the cap and add 4 drops of HI 93764-0 Nessler Reagent.
- Replace the cap tightly and mix by inverting the vial a couple of times. This is the sample.
- Place the vial into the holder and push it completely down.
- Press **Timer** and the display will show the countdown prior to the measurement and the message “Reaction Time”. Alternatively, wait for 3 minutes and 30 seconds and press **Read**. Wait for vial identification. If that was successfully done, the instrument will perform the reading.



- The instrument displays concentration in **mg/L of ammonia nitrogen (NH₃-N)**. Press **▲** or **▼** to access the second level functions and then press the **Chem Frm** key to convert the result in mg/L ammonia (NH₃).



- Press **▲** or **▼** to return to the measurement screen.

INTERFERENCES

Interference may be caused by:

- Organic compounds like: chloramines, various aliphatic and aromatic amines, glycine or urea above 100 ppm. To eliminate these interferences distillation is required.
- Organic compounds like: aldehydes, alcohols (e.g. ethanol) or acetone above 1 %.
To eliminate these interferences distillation is required.
- Sulfide: may cause turbidity.

CHLORINE, FREE

SPECIFICATIONS

Range	0.00 to 5.00 mg/L Cl ₂
Resolution	0.01 mg/L
Accuracy	±0.03 mg/L or ±4 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 525 nm
Method	Adaptation of the EPA method 330.5 and <i>Standard Methods for the Examination of Water and Wastewater</i> , 20 th edition, 4500-Cl G, DPD method. The reaction between free chlorine and the DPD reagent causes a pink tint in the sample.

REQUIRED REAGENTS

Code	Description	Q.ty/test
HI 93701-0	DPD Powder Reagent	1 packet

REAGENT SETS


HI 93701-01 Reagents for 100 tests

HI 93701-03 Reagents for 300 tests




For other accessories see page 70.

MEASUREMENT PROCEDURE



- Select **Chlorine, Free** method following one of the procedures described in the “Method Selection” section (see page 15).
- Take an empty vial.
- Fill the vial with 10 mL of unreacted sample, then replace the cap. This is the blank.
- Place the vial into the holder and push it completely down.





- Press **Zero**. The display will show “-0.0-” when the meter is zeroed and ready for measurement.




- Remove the cap and add the content of one packet of HI 93701-0 Free Chlorine Reagent.
- Replace the cap and shake gently to mix for about 20 seconds. This is the sample.
- Place the vial into the holder and push it completely down.



- Press **Timer** and the display will show the countdown prior to the measurement and the message "Reaction Time". Alternatively, wait for 1 minute and press **Read**. The instrument will perform the reading.



- The instrument displays concentration in **mg/L of free chlorine**. Press ▲ or ▼ to access the second level functions and then press the **Log** key to store the reading.



- Press ▲ or ▼ to return to the measurement screen.

INTERFERENCES

Interferences may be caused by:

Bromine (Br_2)

Iodine (I_2)

Oxidized forms of Chromium and Manganese

Ozone (O_3)

Alkalinity above 250 mg/L CaCO_3 or acidity above 150 mg/L CaCO_3 will not reliably develop the full amount of color or it may rapidly fade. To resolve this, neutralize the sample with diluted HCl or NaOH.

In case of water with hardness greater than 500 mg/L CaCO_3 , shake the sample for approximately 2 minutes after adding the powder reagent.

CHLORINE, TOTAL

SPECIFICATIONS

Range	0.00 to 5.00 mg/L Cl ₂
Resolution	0.01 mg/L
Accuracy	±0.03 mg/L or ±4 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 525 nm
Method	Adaptation of the EPA method 330.5 and <i>Standard Methods for the Examination of Water and Wastewater</i> , 20 th edition, 4500-Cl G, DPD method. The reaction between chlorine and the DPD reagent causes a pink tint in the sample.

REQUIRED REAGENTS

Code	Description	Q.ty/test
HI 93711-0	DPD Powder Reagent	1 packet

REAGENT SETS


HI 93711-01 Reagents for 100 tests

HI 93711-03 Reagents for 300 tests



For other accessories see page 70.

MEASUREMENT PROCEDURE



- Select **Chlorine, Total** method following one of the procedures described in the “Method Selection” section (see page 15).
- Take an empty vial.
- Fill the vial with 10 mL of unreacted sample, then replace the cap. This is the blank.
- Place the vial into the holder and push it completely down.




- Press **Zero**. The instrument will perform a zero sequence and after a few seconds the display will show “-0.0-”. Now the meter is zeroed and ready for measurement.




- Remove the cap and add the content of one packet of HI 93711-0 Total Chlorine Reagent.
- Replace the cap and shake gently to mix for about 20 seconds. This is the sample.
- Place the vial into the holder and push it completely down.



- Press **Timer** and the display will show the countdown prior to the measurement and the message "Reaction Time". Alternatively, wait for 2 minutes and 30 seconds and press **Read**. The instrument will perform the reading.



- The instrument directly displays concentration in **mg/L of total chlorine**. Press **▲** or **▼** to access the second level functions and then press the **Log** key to store the reading.



- Press **▲** or **▼** to return to the measurement screen.

INTERFERENCES

Interferences may be caused by:

Bromine (Br_2)

Iodine (I_2)

Oxidized forms of Chromium and Manganese

Ozone (O_3)

Alkalinity above 250 mg/L CaCO_3 or acidity above 150 mg/L CaCO_3 will not reliably develop the full amount of color or it may rapidly fade. To resolve this, neutralize the sample with diluted HCl or NaOH.

In case of water with hardness greater than 500 mg/L CaCO_3 , shake the sample for approximately 2 minutes after adding the powder reagent.

NITRATE

SPECIFICATIONS

Range	0.0 to 30.0 mg/L NO ₃ ⁻ -N
Resolution	0.1 mg/L
Accuracy	±1.0 mg/L or ±5 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 420 nm
Method	Chromotropic acid method. The reaction between nitrate and the reagents causes a yellow tint in the sample.

REQUIRED REAGENTS

Code	Description	Q.ty/test	Q.ty/set
HI 94766V-0*	Reagent Vial	1 vial	50 vials
HI 93766-0	Nitrate Reagent	1 packet	50 packets

* *Reagent Vial identification: 05xx (xx represents the reagent lot code).*

Note: Store the unused vials in a cool and dark place.

REAGENT SET

HI 94766-50 Reagents for 50 tests

For other accessories see page 70.

MEASUREMENT PROCEDURE



Before using the reagent kit carefully read all the instructions and the Material Safety Data Sheet (MSDS). Pay particular attention to all warnings, cautions and notes. Failure to do so may result in serious injury to the operator.

- Select the **Nitrate** method following one of the procedures described in the "Method Selection" section (see page 15).
- Remove the cap from a Reagent Vial.



- Add exactly 1.0 mL of sample to the vial, while keeping the vial at a 45-degree angle.

- Replace the cap tightly and invert the vial 10 times. This is the blank.




Warning: the vial will become hot during mixing, be careful when handling it.

Note: The method is technique sensitive: to obtain reproducible results it is strongly recommended to follow carefully the "invert" procedure described on page 10.

- Place the vial into the holder and push it completely down.






- Press **Zero** and wait for vial identification. If that was successfully done, the display will show “-0.0-” when the meter is zeroed and ready for measurement.



- Remove the vial.
- Remove the cap and add the content of one packet of HI 93766-0 Nitrate Reagent.
- Replace the cap tightly and invert the vial 10 times. This is the reacted sample.

Note: the method is technique sensitive: to obtain most reproducible results it is strongly recommended to follow carefully the “invert” procedure described on page 10.

- Place the vial into the holder and push it completely down.
- Press **Timer** and the display will show the countdown prior to the measurement and the “Reaction Time” message. Alternatively, wait for 5 minutes and press **Read**. Wait for vial identification. If that was successfully done, the instrument will perform the reading.



- The instrument displays concentration in **mg/L** of nitrate-nitrogen ($\text{NO}_3\text{-N}$). Press \blacktriangle or \blacktriangledown to access the second level functions and then press the **Chem Frm** key to convert the result in mg/L of nitrate (NO_3^-).



- Press \blacktriangle or \blacktriangledown to return to the measurement screen.

INTERFERENCES

Interference may be caused by:

Barium (Ba^{2+}) above 1 mg/L

Chloride (Cl^-) above 1000 mg/L

Nitrite (NO_2^-) above 50 mg/L

Samples containing up to 100 mg/L nitrite may be measured after the following treatment: add 400 mg of urea to 10 mL of sample, mix until completely dissolved, then proceed with the usual measurement procedure.

NITROGEN, TOTAL LOW RANGE

SPECIFICATIONS

Range	0.0 to 25.0 mg/L N
Resolution	0.1 mg/L
Accuracy	± 1.0 mg/L or ± 5 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 420 nm
Method	Chromotropic acid method. A persulfate digestion converts all forms of nitrogen to nitrate. Then the reaction between nitrate and the reagents causes a yellow tint in the sample.

REQUIRED REAGENTS

<u>Code</u>	<u>Description</u>	<u>Q.ty/test</u>	<u>Q.ty/set</u>
HI 94767A-B*	Digestion Vial	1 vial	50 vials
DEIONIZED120	Deionized Water	2 mL	1 bottle
PERSULFATE/N	Potassium Persulfate	1 packet	50 packets
BISULFITE/N	Sodium Metabisulfite	1 packet	50 packets
HI 93767-0	Total Nitrogen Reagent	1 packet	50 packets
HI 94766V-OLR**	Reagent Vial	1 vial	50 vials

* *Digestion Vial identification: 16xx (xx represents the reagent lot code). GREEN RECTANGLE ON THE LABEL.*

** *Reagent Vial identification: 06xx (xx represents the reagent lot code). RED RECTANGLE ON THE LABEL.*

Notes:

- The instrument does not recognize the Digestion Vial (16xx).
- Store the unused vials in their container in a cool and dark place.

REAGENT SET

HI 94767A-50 Reagents for up to 49 tests. Contains:

Box 1: HI 94767A-50 Reagent Set

Box 2: HI 94767A&B-50 Reagent Set, for Nitrogen Total Low Range.

For other accessories see page 70.

REQUIRED ACCESSORIES

HI 839800-01	Hanna reactor (115 VAC)
HI 839800-02	Hanna reactor (230 VAC)
HI 740216	Test tube cooling rack (25 holes)
HI 740217	Laboratory bench safety shield

For other accessories see page 70.

MEASUREMENT PROCEDURE



Before using the reagent kit carefully read all the instructions and the Material Safety Data Sheet (MSDS). Pay particular attention to all warnings, cautions and notes. Failure to do so may result in serious injury to the operator.


Reagent Blank Correction: This method requires a reagent blank correction. A single blank vial may be used more than once; the blank vial is stable up to one week if stored in a dark place at room temperature. Always use the same lot of reagents for blank and samples. For most accurate measurement run a blank for each set of measurements.

- Preheat the Hanna Reactor HI 839800 to 105 °C (221 °F). For correct use of the reactor follow Reactor Instruction Manual.


The optional HI 740217 safety shield is strongly recommended.

DO NOT USE AN OVEN OR MICROWAVE samples may leak and generate a corrosive and possibly explosive atmosphere.


- Remove the cap from two Digestion Vials (16xx code vials).




- Add the content of one packet of Potassium Persulfate for Total Nitrogen analysis to each vial.



- Add exactly 2.0 mL of sample to one vial (sample vial), and 2.0 mL of deionized water to the other vial (blank vial), while keeping the vials at a 45-degree angle.



- Replace the cap tightly and shake the vials vigorously for about 30 seconds until all the powder is completely dissolved.



- Insert the vials into the reactor and heat them for 30 minutes at 105 °C.
Note: To obtain most accurate results, it is strongly recommended to remove the vials from the reactor after 30 minutes.



- At the end of the digestion period switch off the reactor, place the vials in the test tube rack and allow to cool to room temperature.


Warning: The vials are still hot, be careful in handling them.




- Select **Nitrogen, Total LR** method following one of the procedures described in the “Method Selection” section (see page 15).
- For this method the instrument provides 3 reaction timers which can be used throughout the procedure.
- Remove the cap from the vials and add the content of one packet of Sodium Metabisulfite for Total Nitrogen analysis to each vial. Replace the cap tightly and shake gently the vials for 15 seconds.




- Press **Timer** to start the 3 minutes timer or wait for 3 minutes (without shaking the vials) to allow the reaction to complete.




- Remove the cap from the vials and add the content of one packet of HI 93767-0 Total Nitrogen Reagent to each vial. Replace the cap tightly and shake gently the vials for 15 seconds.




- Press **Start** to activate the 2 minutes timer or wait for 2 minutes (without shaking the vials) to allow the reaction to complete.



- Remove the cap from two Reagent Vials (06xx code vials).




- Add exactly 2.0 mL of digested sample (sample vial) to one Reagent Vial (sample vial), and 2.0 mL of digested blank (blank vial) to the other Reagent Vial (blank vial), while keeping the vials at a 45-degree angle.
- Replace the cap tightly and invert the vials 10 times.



Warning: The vials will become hot during mixing, be careful when handling them.



Note: The method is technique sensitive; to obtain most reproducible results it is strongly recommended to follow carefully the "invert" procedure described on page 10.

- Place the blank vial into the holder and push it completely down.




Nitrogen, Total LR

- Press **Start** and the display will show the countdown prior to the measurement and the “Reaction Time” message. Alternatively, wait for 5 minutes and press **Zero**.






- Wait for vial identification. If that was successfully done, the display will show “-0.0-” when the meter is zeroed and ready for measurement.



- Remove the blank vial.




- Place the sample vial into the holder and push it completely down.

- Press **Read** and wait for vial identification. If that was successfully done, the instrument will perform the reading. The instrument displays concentration in **mg/L** of total nitrogen (N).



- Press **▲** or **▼** to access the second level functions.

- Press the **Chem Frm** key to convert the result in mg/L of ammonia (NH_3) and mg/L of nitrate (NO_3^-).



- Press **▲** or **▼** to return to the measurement screen.

The method detects all organic and inorganic forms of nitrogen present in the sample.

INTERFERENCES

Interference may be caused by:

Bromide (Br^-) above 60 mg/L

Chloride (Cl^-) above 1000 mg/L

Chromium (Cr^{3+}) above 0.5 mg/L

NITROGEN, TOTAL HIGH RANGE

SPECIFICATIONS

Range	10 to 150 mg/L N
Resolution	1 mg/L
Accuracy	± 3 mg/L or ± 4 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 420 nm
Method	Chromotropic acid method. A persulfate digestion converts all forms of nitrogen to nitrate. Then the reaction between nitrate and the reagents causes a yellow tint in the sample.

REQUIRED REAGENTS

<u>Code</u>	<u>Description</u>	<u>Q.ty/test</u>	<u>Q.ty/set</u>
HI 94767B-B*	Digestion Vial	1 vial	50 vials
DEIONIZED120	Deionized Water	0.5 mL	1 bottle
PERSULFATE/N	Potassium Persulfate	1 packet	50 packets
BISULFITE/N	Sodium Metabisulfite	1 packet	50 packets
HI 93767-0	Total Nitrogen Reagent	1 packet	50 packets
HI 94766V-OHR**	Reagent Vial	1 vial	50 vials

* *Digestion Vial identification: 17xx (xx represents the reagent lot code). RED RECTANGLE ON THE LABEL.*

** *Reagent Vial identification: 07xx (xx represents the reagent lot code). GREEN RECTANGLE ON THE LABEL.*

Notes: • The instrument does not recognize the Digestion Vial (17xx).

• Store the unused vials in their container in a cool and dark place.

REAGENT SET

HI 94767B-50 Reagents for up to 49 tests. Contains:

Box 1: HI 94767B-50 Reagent Set

Box 2: HI 94767A&B-50 Reagent Set, for Nitrogen Total High Range.

For other accessories see page 70.

REQUIRED ACCESSORIES

HI 839800-01 Hanna reactor (115 VAC)

HI 839800-02 Hanna reactor (230 VAC)

HI 740216 Test tube cooling rack (25 holes)

HI 740217 Laboratory bench safety shield

For other accessories see page 70.

MEASUREMENT PROCEDURE




Before using the reagent kit carefully read all the instructions and the Material Safety Data Sheet (MSDS). Pay particular attention to all warnings, cautions and notes. Failure to do so may result in serious injury to the operator.

Reagent Blank Correction: This method requires a reagent blank correction. A single blank vial may be used more than once; the blank vial is stable up to one week if stored in a dark place at room temperature. Always use the same lot of reagents for blank and samples. For most accurate measurement, run a blank for each set of measurements.


- Preheat the Hanna Reactor HI 839800 to 105 °C (221 °F). For correct use of the reactor follow Reactor Instruction Manual. Use of the optional HI 740217 safety shield is strongly recommended.

DO NOT USE AN OVEN OR MICROWAVE samples may leak and generate a corrosive and possibly explosive atmosphere.


- Remove the cap from two Digestion Vials (17xx cap vials).




- Add the content of one packet of Potassium Persulfate for Total Nitrogen analysis to each vial.



- Add exactly 0.5 mL of sample to one vial (sample vial), and 0.5 mL of deionized water to the other vial (blank vial), while keeping the vials at a 45-degree angle.



- Replace the cap tightly and shake vigorously the vials for about 30 seconds until all the powder is completely dissolved.




- Insert the vials into the reactor and heat them for 30 minutes at 105°C.
Note: to obtain most accurate results, it is strongly recommended to remove the vials from the reactor after 30 minutes.





- At the end of the digestion place the vials in the test tube rack and allow to cool to room temperature.
Warning: as the vials are still hot, be careful in handling them.




- Select **Nitrogen, Total HR** method following one of the procedures described in the "Method Selection" section (see page 15).
- For this method the instrument provides 3 reaction timers which can be used throughout the procedure.
- Remove the cap from the vials and add the content of one packet of Sodium Metabisulfite for Total Nitrogen analysis to each vial. Replace the cap tightly and shake gently the vials for 15 seconds.





- Press **Timer** to start the 3 minutes timer or wait for 3 minutes (without shaking the vials) to allow the reaction to complete.




- Remove the cap from the vials and add the content of one packet of HI 93767-0 Total Nitrogen Reagent to each vial. Replace the cap tightly and shake gently the vials for 15 seconds.




- Press **Start** to activate the 2 minutes timer or wait for 2 minutes (without shaking the vials) to allow the reaction to complete.



- Remove the cap from two Reagent Vials (07xx code vials).




- Add exactly 2.0 mL of digested sample (sample vial) to one Reagent Vial (sample vial), and 2.0 mL of digested blank (blank vial) to the other vial (blank vial), while keeping the vials at a 45-degree angle.
- Replace the cap tightly and invert the vials 10 times.






Warning: the vials will become hot during mixing, be careful when handling them.

Note: the method is technique sensitive: to obtain most reproducible results it is strongly recommended to follow carefully the "invert" procedure described on page 10.




- Place the blank vial into the holder and push it completely down.





- Press **Start** and the display will show the countdown prior to the measurement and the "Reaction Time" message. Alternatively, wait for 5 minutes and press **Zero**.



- Wait for vial identification. If that was successfully done, the instrument will perform a zero sequence and after a few seconds the display will show "-0.0-". Now the meter is zeroed and ready for measurement.
- Remove the blank vial and place the sample vial into the holder and push it completely down.
- Press **Read** and wait for vial identification. If that was successfully done, the instrument will perform the reading. The instrument displays the concentration in mg/L of total nitrogen (N).



- Press **▲** or **▼** to access the second level functions.
- Press the **Chem Frm** key to convert the result in mg/L of ammonia (NH_3) and mg/L of nitrate (NO_3^-).



- Press **▲** or **▼** to return to the measurement screen.

The method detects all organic and inorganic forms of nitrogen present in the sample.

INTERFERENCES

Interference may be caused by:

Bromide (Br^-) above 240 mg/L

Chloride (Cl^-) above 3000 mg/L

Chromium (Cr^{3+}) above 0.5 mg/L

OXYGEN DEMAND, CHEMICAL LOW RANGE

SPECIFICATIONS

Range	0 to 150 mg/L COD
Resolution	1 mg/L
Accuracy	± 5 mg/L or ± 5 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 420 nm
Method	Adaptation of the USEPA 410.4 approved method for the COD determination on surface waters and wastewaters. Oxidizable organic compounds reduce the dichromate ion (orange) to the chromic ion (green). The amount of remaining dichromate is determined.

REQUIRED REAGENTS

Code	Description	Q.ty/test	Q.ty/set
COD94LR*	Reagent Vial	1 vial	25 vials
DEIONIZED120	Deionized Water	2 mL	OPTIONAL

* *Reagent Vial identification: 12xx (xx represents the reagent lot code). RED RECTANGLE ON THE LABEL.*

Note: Store the unused vials in their container in a cool and dark place.

REAGENT SET

HI 94754A-25 Reagents for up to 24 tests

REQUIRED ACCESSORIES

- HI 839800-01 Hanna reactor (115 VAC)
- HI 839800-02 Hanna reactor (230 VAC)
- HI 740216 Test tube cooling rack (25 holes)
- HI 740217 Laboratory bench safety shield

For other accessories see page 70.


MEASUREMENT PROCEDURE




Before using the reagent kit carefully read all the instructions and the Material Safety Data Sheet (MSDS). Pay particular attention to all warnings, cautions and notes. Failure to do so may result in serious injury to the operator.

Reagent Blank Correction: This method requires a reagent blank correction. A single blank vial may be used more than once. The blank vial is stable for several months (room temperature). For most accurate measurement, run a blank for each set of measurements and always use the same lot of reagents for blank and samples.

- Choose a homogeneous sample. Samples containing settleable solids need to be homogenized with a blender.
- Preheat the Hanna Reactor HI 839800 to 150 °C (302 °F). For correct use of the reactor follow Reactor Instruction Manual. The optional HI 740217 safety shield is strongly recommended.
DO NOT USE AN OVEN OR MICROWAVE samples may leak and generate a corrosive and possibly explosive atmosphere.
- Remove the cap from two Reagent Vials.



- Add exactly 2.0 mL of sample to one vial (sample vial), and 2.0 mL of deionized water to the other vial (blank vial), while keeping the vials at a 45-degree angle. Replace the cap tightly and mix by inverting each vial a couple of times.




Warning: The vials will become hot during mixing, be careful when handling them.

- Insert the vials into the reactor and heat them for 2 hours at 150 °C.



- At the end of the digestion period switch off the reactor. Wait for twenty minutes to allow the vials to cool to about 120 °C.



- Invert each vial several times while still warm, then place them in the test tube rack.

Warning: The vials are still hot, be careful when handling them.

- Leave the vials in the tube rack to cool to room temperature. Do not shake or invert them, the samples may become turbid.






- Select **Oxygen Demand, (COD) LR** method following one of the procedures described in the “Method Selection” section (see page 15).




- Place the blank vial into the holder and push it completely down.



- Press **Zero** and wait for vial identification. If that was successfully done, the instrument will perform a zero sequence and after a few seconds the display will show “-0.0-”. Now the meter is zeroed and ready for measurement.



- Remove the blank vial.
- Place the sample vial into the holder and push it completely down.
- Press **Read** and wait for vial identification. If that was successfully done, the instrument will perform the reading.



- The instrument displays concentration in **mg/L** of oxygen demand.

INTERFERENCES

Interference may be caused by:

Chloride (Cl⁻) above 2000 mg/L.

Samples with higher chloride concentration should be diluted.

OXYGEN DEMAND, CHEMICAL MEDIUM RANGE

SPECIFICATIONS

Range	0 to 1500 mg/L COD
Resolution	1 mg/L
Accuracy	± 15 mg/L or ± 4 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 610 nm
Method	Adaptation of the USEPA 410.4 approved method for the COD determination on surface waters and wastewaters. Oxidizable organic compounds reduce the dichromate ion (orange) to the chromic ion (green). The amount of chromic ion formed is determined.

REQUIRED REAGENTS

Code	Description	Q.ty/test	Q.ty/set
COD94MR*	Reagent Vial	1 vial	25 vials
DEIONIZED120	Deionized Water	2 mL	OPTIONAL

* Reagent Vial identification: **13xx** (xx represents the reagent lot code).

Note: Store the unused vials in their container in a cool and dark place.

REAGENT SET

HI 94754B-25 Reagents for up to 24 tests

REQUIRED ACCESSORIES

HI 839800-01 Hanna reactor (115 VAC)

HI 839800-02 Hanna reactor (230 VAC)

HI 740216 Test tube cooling rack (25 holes)

HI 740217 Laboratory bench safety shield

For other accessories see page 70.

MEASUREMENT PROCEDURE



Before using the reagent kit carefully read all the instructions and the Manual Safety Data Sheet (MSDS). Pay particular attention to all warnings, cautions and notes. Failure to do so may result in serious injury to the operator.

Reagent Blank Correction: This method requires a reagent blank correction. A single blank vial may be used more than once. The blank vial is stable for several months (room temperature). For most accurate measurement, run a blank for each set of measurements and always use the same lot of reagents for blank and samples.

- Choose a homogeneous sample. Samples containing settleable solids need to be homogenized with a blender.
- Preheat the Hanna Reactor HI 839800 to 150 °C (302 °F). For correct use of the reactor follow Reactor Instruction Manual. Use of the optional HI 740217 safety shield is strongly recommended.

DO NOT USE AN OVEN OR MICROWAVE samples may leak and generate a corrosive and possibly explosive atmosphere.


- Remove the cap from two Reagent Vials.



- Add exactly 2.0 mL of sample to one vial (sample vial), and 2.0 mL of deionized water to the other vial (blank vial), while keeping the vials at a 45-degree angle. Replace the cap tightly and mix by inverting each vial a couple of times.

Warning: The vials will become very hot during mixing, be careful when handling them.

- Insert the vials into the reactor and heat them for 2 hours at 150 °C.




- At the end of the digestion period switch off the reactor. Wait for twenty minutes to allow the vials to cool to about 120 °C.

- Invert each vial several times while still warm, then place them in a test tube rack.


Warning: The vials are still hot, be careful when handling them.

- Leave the vials in the tube rack to cool to room temperature. Do not shake or invert them, the samples may become turbid.






- Select the **Oxygen Demand, (COD) MR** method following one of the procedures described in the “Method Selection” section (see page 15).




- Place the blank vial into the holder and push it completely down.



- Press **Zero** and wait for vial identification. If that was successfully done, the instrument will perform a zero sequence and after a few seconds the display will show “-0.0-”. Now the meter is zeroed and ready for measurement.



- Remove the blank vial.
- Place the sample vial into the holder and push it completely down.
- Press **Read** and wait for vial identification. If that was successfully done, the instrument will perform the reading.



- The instrument displays concentration in **mg/L** of oxygen demand on the LCD.

INTERFERENCES

Interference may be caused by:

Chloride (Cl^-) above 2000 mg/L.

Samples with higher chloride concentration should be diluted.

OXYGEN DEMAND, CHEMICAL HIGH RANGE

SPECIFICATIONS

Range	0 to 15000 mg/L COD
Resolution	10 mg/L
Accuracy	± 150 mg/L or ± 3 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 610 nm
Method	Adaptation of the USEPA 410.4 approved method for the COD determination on surface waters and wastewaters. Oxidizable organic compounds reduce the dichromate ion (orange) to the chromic ion (green). The amount of chromic ion formed is determined.

REQUIRED REAGENTS

Code	Description	Q.ty/test	Q.ty/set
COD94HR*	Reagent Vial	1 vial	25 vials
DEIONIZED120	Deionized Water	0.2 mL	OPTIONAL

* *Reagent Vial identification: 24xx (xx represents the reagent lot code). GREEN RECTANGLE ON THE LABEL.*

Note: Store the unused vials in their container in a cool and dark place.

REAGENT SET

HI 94754C-25 Reagents for up to 24 tests

REQUIRED ACCESSORIES

- HI 839800-01 Hanna reactor (115 VAC)
- HI 839800-02 Hanna reactor (230 VAC)
- HI 740216 Test tube cooling rack (25 holes)
- HI 740217 Laboratory bench safety shield

For other accessories see page 70.

MEASUREMENT PROCEDURE



Before using the reagent kit carefully read all the instructions and the Material Safety Data Sheet (MSDS). Pay particular attention to all warnings, cautions and notes. Failure to do so may result in serious injury to the operator.

Reagent Blank Correction: This method requires a reagent blank correction. A single blank vial may be used more than once. The blank vial is stable for several months (room temperature). For most accurate measurement, run a blank for each set of measurements and always use the same lot of reagents for blank and samples.

- Choose a homogeneous sample. Samples containing settleable solids need to be homogenized with a blender.
- Preheat the Hanna Reactor HI 839800 to 150 °C (302 °F). For correct use of the reactor follow Reactor Instruction Manual. Use of the optional HI 740217 safety shield is strongly recommended.

DO NOT USE AN OVEN OR MICROWAVE samples may leak and generate a corrosive and possibly explosive atmosphere.


- Remove the cap from two Reagent Vials.



- Add exactly 0.2 mL of sample to one vial (sample vial), and 0.2 mL of deionized water to the other vial (blank vial), while keeping the vials at a 45-degree angle. Replace the cap tightly and mix by inverting each vial a couple of times.

Warning: The vials will become very hot during mixing, be careful when handling them.

- Insert the vials into the reactor and heat them for 2 hours at 150 °C.




- At the end of the digestion period switch off the reactor. Wait for twenty minutes to allow the vials to cool to about 120 °C.

- Invert each vial several times while still warm, then place them in the test tube rack.


Warning: The vials are still hot, be careful when handling them.

- Leave the vials in the tube rack to cool to room temperature. Do not shake or invert them anymore otherwise the samples may become turbid.






- Select the *Oxygen Demand, (COD) HR* method following one of the procedures described in the “Method Selection” section (see page 15).




- Place the blank vial into the holder and push it completely down.



- Press **Zero** and wait for vial bar code identification. If that was successfully done, the instrument will perform a zero sequence and after a few seconds the display will show “-0.0-”. Now the meter is zeroed and ready for measurement.
- Remove the blank vial.



- Place the sample vial into the holder and push it completely down.
- Press **Read** and wait for vial identification. If that was successfully done, the instrument will perform the reading.
- The instrument directly displays concentration in **mg/L** of oxygen demand.



INTERFERENCES

Interference may be caused by:

Chloride (Cl^-) above 20000 mg/L.

Samples with higher chloride concentration should be diluted.

PHOSPHORUS, REACTIVE

SPECIFICATIONS

Range	0.00 to 1.60 mg/L P
Resolution	0.01 mg/L
Accuracy	± 0.05 mg/L or ± 5 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 610 nm
Method	Adaptation of the EPA method 365.2 and <i>Standard Methods for the Examination of Water and Wastewater</i> , 20 th edition, 4500-P E, ascorbic acid method. The reaction between orthophosphate and the reagent causes a blue tint in the sample.

REQUIRED REAGENTS

Code	Description	Q.ty/test	Q.ty/set
HI 94758A-0*	Reagent Vial	1 vial	50 vials
HI 93758-0	Phosphorus Reagent	1 packet	50 packets

* Reagent Vial identification: **30xx** (xx represents the reagent lot code). RED RECTANGLE ON THE LABEL.





REAGENT SET

HI 94758A-50 Reagents for 50 tests

For other accessories see page 70.



MEASUREMENT PROCEDURE

- Select the **Phosphorus, Reactive** method following one of the procedures described in the "Method Selection" section (see page 15).
- Remove the cap from a Reagent Vial.
- Add exactly 5.0 mL of sample to the vial, while keeping the vial at a 45-degree angle.
- Replace the cap and mix by inverting the vial a couple of times. This is the blank.
- Place the vial into the holder and push it completely down.
- Press **Zero** and wait for vial identification. If that was successfully done, the instrument will perform a zero sequence and after a few seconds the display will show "**-0.0-**". Now the meter is zeroed and ready for measurement.




- Remove the vial.


- Remove the cap and add the content of one packet of HI 93758-0 Phosphorus Reagent.
- Replace the cap tightly and shake gently to mix for 2 minutes until most of the powder is dissolved. This is the reacted sample.
- Place the vial into the holder and push it completely down.
- Press **Timer** and the display will show the countdown prior to the measurement and the "Reaction Time" message. Alternatively, wait for 3 minutes and press **Read**.



- Wait for vial identification. If that was successfully done, the instrument will perform the reading. The instrument displays concentration in mg/L of phosphorus (P).



- Press \blacktriangle or \blacktriangledown to access the second level functions and then press the **Chem Frm** key to convert the result in mg/L of phosphate (PO_4^{3-}) and mg/L of P_2O_5 .



- Press \blacktriangle or \blacktriangledown to return to the measurement screen.

Note: for accurate measurements

- 1) wash glassware only with phosphate-free detergents
- 2) clean all glassware with 1 : 1 hydrochloric acid solution and rinse with deionized water.

INTERFERENCES

Interference may be caused by:

Arsenate at any level

Silica above 50 mg/L

Sulfide above 6 mg/L.

To eliminate sulfide: add Bromine Water drop-wise until a pale yellow color develops; remove Bromine Water excess by adding Phenol solution drop-wise.

Turbidity and suspended matter in large amounts may cause interference because the reaction conditions may dissolve suspended matter or cause desorption of phosphates from particles. Turbidity or suspended matter should be removed before measurement by treatment with active carbon and by prior filtration.

PHOSPHORUS, ACID HYDROLYZABLE

SPECIFICATIONS

Range	0.00 to 1.60 mg/L P
Resolution	0.01 mg/L
Accuracy	± 0.05 mg/L or ± 5 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 610 nm
Method	Adaptation of the EPA method 365.2 and <i>Standard Methods for the Examination of Water and Wastewater</i> , 20 th edition, 4500-P E, ascorbic acid method. A mild acid digestion converts condensed inorganic forms of phosphates to orthophosphate. Then the reaction between orthophosphate and the reagents causes a blue tint in the sample.

REQUIRED REAGENTS

Code	Description	Q.ty/test	Q.ty/set
HI 94758V-OAH*	Reagent Vial	1 vial	50 vials
HI 93758B-0	NaOH Solution 1.20 N 2 mL		1 bottle
HI 93758-0	Phosphorus Reagent	1 packet	50 packets

* Reagent Vial identification: 31xx (xx represents the reagent lot code).

Note: Store the unused vials in their container in a cool and dark place.

REAGENT SET

HI 94758B-50 Reagents for 50 tests

For other accessories see page 70.

REQUIRED ACCESSORIES

HI 839800-01 Hanna reactor (115 VAC)

HI 839800-02 Hanna reactor (230 VAC)

HI 740216 test tube cooling rack (25 holes)

HI 740217 Laboratory bench safety shield

For other accessories see page 70.

MEASUREMENT PROCEDURE



Before using the reagent kit carefully read all the instructions and the Material Safety Data Sheet (MSDS). Pay particular attention to all warnings, cautions and notes. Failure to do so may result in serious injury to the operator.

- Preheat the Hanna Reactor HI 839800 to 150 °C (302 °F). For correct use of the reactor follow Reactor Instruction Manual.

The optional HI 740217 safety shield is strongly recommended.

DO NOT USE AN OVEN OR MICROWAVE samples may leak and generate a corrosive and possibly explosive atmosphere.

- Remove the cap from a Reagent Vial.
- Add exactly 5.0 mL of sample to the vial, while keeping the vial at a 45-degree angle.
- Replace the cap and mix by inverting the vial a couple of times.
- Insert the vial into the reactor and heat it for 30 minutes at 150 °C.





- At the end of the digestion place the vials carefully in the test tube rack and allow to cool to room temperature.






Warning: The vials are still hot, be careful when handling them.



- Select the **Phosphorus, Acid Hydrolyzable** method following one of the procedures described in the "Method Selection" section (see page 15).
- Remove the cap from the vial and add exactly 2.0 mL of Sodium Hydroxide (NaOH) Solution 1.20 N, while keeping the vial at a 45-degree angle.
- Replace the cap tightly and mix by inverting the vial a couple of times. This is the blank.
- Place the vial into the holder and push it completely down.




- Press **Zero** and wait for vial identification. If that was successfully done, the instrument will perform a zero sequence and after a few seconds the display will show “-0.0-”. Now the meter is zeroed and ready for measurement.



- Remove the vial.
- Remove the cap and add the content of one packet of HI 93758-0 Phosphorus Reagent.
- Replace the cap tightly and shake gently to mix for 2 minutes until most of the powder is dissolved. This is the sample.
- Place the vial into the holder and push it completely down.
- Press **Timer** and the display will show the countdown prior to the measurement and the “Reaction Time” message. Alternatively, wait for 3 minutes and press **Read**.





- Wait for vial identification. If that was successfully done, the instrument will perform the reading. The instrument displays concentration in **mg/L of phosphorus (P)**.



The method detects free (orthophosphate) and condensed inorganic forms (meta-, pyro- and other polyphosphates) of phosphates present in the sample.

- Press **▲** or **▼** to access the second level functions and then press the **Chem Frm** key to convert the result in mg/L of phosphate (PO_4^{3-}) and mg/L of P_2O_5 .



- Press **▲** or **▼** to return to the measurement screen.

Note: for accurate measurements

1) wash glassware only with phosphate-free detergents

2) clean all glassware with 1 : 1 hydrochloric acid solution and rinse with deionized water.

INTERFERENCES

Interference may be caused by:

Arsenate at any level

Silica above 50 mg/L

Sulfide above 9 mg/L.

To eliminate sulfide: add Bromine Water drop-wise until a pale yellow color develops; remove Bromine Water excess by adding Phenol solution drop-wise.

Turbidity and suspended matter in large amounts may cause interference because the strongly acidic reaction conditions may dissolve suspended matter or cause desorption of phosphates from particles. Before measurement, turbidity or suspended matter should be removed by treatment with active carbon and by prior filtration.

PHOSPHORUS, TOTAL LOW RANGE

SPECIFICATIONS

Range	0.00 to 1.15 mg/L P
Resolution	0.01 mg/L
Accuracy	± 0.05 mg/L or ± 6 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 610 nm
Method	Adaptation of the EPA method 365.2 and <i>Standard Methods for the Examination of Water and Wastewater</i> , 20 th edition, 4500-P E, ascorbic acid method. A persulfate digestion converts organic and condensed inorganic forms of phosphates to orthophosphate. Then the reaction between orthophosphate and the reagents causes a blue tint in the sample.

REQUIRED REAGENTS

Code	Description	Q.ty/test	Q.ty/set
HI 94758V-0*	Reagent Vial	1 vial	50 vials
HI 93758C-0	NaOH Solution 1.54 N	2 mL	1 bottle
HI 93758-0	Phosphorus Reagent	1 packet	50 packets
PERSULFATE/P	Potassium Persulfate	1 packet	50 packets

* Reagent Vial identification: **32xx** (xx represents the reagent lot code). RED RECTANGLE ON THE LABEL.

Note: Store the unused vials in their container in a cool and dark place.

REAGENT SET

HI 94758C-50 Reagents for 50 tests

For other accessories see page 70.

REQUIRED ACCESSORIES

HI 839800-01 Hanna reactor (115 VAC)

HI 839800-02 Hanna reactor (230 VAC)

HI 740216 test tube cooling rack (25 holes)

HI 740217 Laboratory bench safety shield

For other accessories see page 70.

MEASUREMENT PROCEDURE



Before using the reagent kit carefully read all the instructions and the Material Safety Data Sheet (MSDS). Pay particular attention to all warnings, cautions and notes. Failure to do so may result in serious injury to the operator.

- Preheat the Hanna Reactor HI 839800 to 150 °C (302 °F). For correct use of the reactor follow Reactor Instruction Manual.

Use of the optional HI 740217 safety shield is strongly recommended.

DO NOT USE AN OVEN OR MICROWAVE samples may leak and generate a corrosive and possibly explosive atmosphere.

- Remove the cap from a Reagent Vial.
- Add exactly 5.0 mL of sample to the vial, while keeping the vial at a 45-degree angle.
- Add the content of one packet of Potassium Persulfate for Phosphorus analysis. Replace the cap and shake gently the vial until all the powder is completely dissolved.




- Insert the vial into the reactor and heat it for 30 minutes at 150 °C.
- At the end of the digestion place the vials carefully in the test tube rack and allow to cool to room temperature.




Warning: The vials are still hot, be careful when handling them.





- Select the ***Phosphorus, Total*** method following one of the procedures described in the "Method Selection" section (see page 15).
- Remove the cap from the vial and add exactly 2.0 mL of Sodium Hydroxide (NaOH) Solution 1.54 N, while keeping the vial at a 45-degree angle.
- Replace the cap tightly and mix by inverting the vial a couple of times. This is the blank.
- Place the vial into the holder and push it completely down.





- Press **Zero** and wait for vial identification. If that was successfully done, the instrument will perform a zero sequence and after a few seconds the display will show “-0.0-”. Now the meter is zeroed and ready for measurement.



- Remove the vial.
- Remove the cap and add the content of one packet of HI 93758-0 Phosphorus Reagent.
- Replace the cap tightly and shake gently to mix for about 2 minutes until all the powder is completely dissolved. This is the sample.
- Place the vial into the holder and push it completely down.
- Press **Timer** and the display will show the countdown prior to the measurement and the “Reaction Time” message. Alternatively, wait for 3 minutes and press **Read**.





- Wait for vial identification. If that was successfully done, the instrument will perform the reading. The instrument displays concentration in **mg/L of phosphorus (P)**.



The method detects free (orthophosphate) and condensed inorganic forms (meta-, pyro- and other polyphosphates) of phosphates present in the sample.

- Press **▲** or **▼** to access the second level functions and then press the **Chem Frm** key to convert the result in mg/L of phosphate (PO_4^{3-}) and mg/L of P_2O_5 .



- Press **▲** or **▼** to return to the measurement screen.

Note: for accurate measurements

1) wash glassware only with phosphate-free detergents

2) clean all glassware with 1 : 1 hydrochloric acid solution and rinse with deionized water.

INTERFERENCES

Interference may be caused by:

Arsenate at any level

Silica above 50 mg/L

Sulfide above 90 mg/L.

Turbidity and suspended matter in large amounts may cause interference because the strongly acidic reaction conditions may dissolve suspended matter or cause desorption of phosphates from particles. Before measurement, turbidity or suspended matter should be removed by treatment with active carbon and by prior filtration.

PHOSPHORUS, REACTIVE HIGH RANGE

SPECIFICATIONS

Range	0.0 to 32.6 mg/L P
Resolution	0.1 mg/L
Accuracy	± 0.5 mg/L or ± 5 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 420 nm
Method	Adaptation of the <i>Standard Methods for the Examination of Water and Wastewater</i> , 20 th edition, 4500-P C, vanadomolybdophosphoric acid method. The reaction between orthophosphate and the reagents causes a yellow tint in the sample.

REQUIRED REAGENTS

Code	Description	Q.ty/test	Q.ty/set
HI 94763A-0*	Reagent Vial	1 vial	50 vials
DEIONIZED120	Deionized Water	5 mL	1 bottle

* *Reagent Vial identification: 33xx (xx represents the reagent lot code). GREEN RECTANGLE ON THE LABEL.*

Note: Store the unused vials in their container in a cool and dark place.

REAGENT SET

HI 94763A-50 Reagents for up to 49 tests

For other accessories see page 70.

MEASUREMENT PROCEDURE


Reagent Blank Correction: This method requires a reagent blank correction. A single blank vial may be used more than once; the blank vial is stable up to two weeks (room temperature). For most accurate measurement, run a blank for each set of measurements and always use the same lot of reagents for blank and samples.

- Select the *Phosphorus, Reactive High Range* method following one of the procedures described in the "Method Selection" section (see page 15).
- Remove the cap from two Reagent Vials.




- Add exactly 5.0 mL of sample to one vial (sample vial), and 5.0 mL of deionized water to the other vial (blank vial), while keeping the vials at a 45-degree angle.

- Replace the cap and mix by inverting each vial a couple of times.





- Place the blank vial into the holder and push it completely down.




- Press **Timer** and the display will show the countdown prior to the measurement and the "Reaction Time" message. Alternatively, wait for 7 minutes and press **Zero**.






- Wait for vial identification. If that was successfully done, the instrument will perform a zero sequence and after a few seconds the display will show "-0.0-". Now the meter is zeroed and ready for measurement.



- Remove the blank vial.
- Place the sample vial into the holder and push it completely down.
- Press **Read** and wait for vial recognition. If that was successfully done, the instrument will perform the reading. The instrument displays concentration in **mg/L of phosphorus (P)**.



- Press ▲ or ▼ to access the second level functions and then press the **Chem Frm** key to convert the result in mg/L of phosphate (PO_4^{3-}) and mg/L of P_2O_5 .



- Press ▲ or ▼ to return to the measurement screen.

Note: for accurate measurements

- 1) wash glassware only with phosphate-free detergents
- 2) clean all glassware with 1 : 1 hydrochloric acid solution and rinse with deionized water.

INTERFERENCES

Interference may be caused by:

Bismuth

Fluoride

pH: the sample should have a neutral pH

Sulfide: to eliminate sulfide add Bromine Water drop-wise until a pale yellow color develops; remove Bromine Water excess by adding Phenol solution drop-wise.

Temperature: the method is temperature sensitive.

It is recommended to run measurements at $T = 20$ to $25\text{ }^{\circ}\text{C}$:

$T < 20\text{ }^{\circ}\text{C}$ causes a negative error

$T > 25\text{ }^{\circ}\text{C}$ causes a positive error

Turbidity and suspended matter in large amounts may cause interference because the strongly acidic reaction conditions may dissolve suspended matter or cause desorption of phosphates from particles. Before measurement, turbidity or suspended matter should be removed by treatment with active carbon and by prior filtration.

PHOSPHORUS, TOTAL HIGH RANGE

SPECIFICATIONS

Range	0.0 to 32.6 mg/L (P)
Resolution	0.1 mg/L
Accuracy	± 0.5 mg/L or ± 5 % of reading @ 25 °C, whichever is greater
Light Source	Tungsten lamp with narrow band interference filter @ 420 nm
Method	Adaptation of the <i>Standard Methods for the Examination of Water and Wastewater</i> , 20 th edition, 4500-P C, vanadomolybdophosphoric acid method. A persulfate digestion converts organic and condensed inorganic forms of phosphates to orthophosphate. Then the reaction between orthophosphate and the reagents causes a yellow tint in the sample.

REQUIRED REAGENTS

Code	Description	Q.ty/test	Q.ty/set
HI 94758V-OHR*	Reagent Vial	1 vial	50 vials
DEIONIZED120	Deionized Water	5 mL	1 bottle
HI 93758C-0	NaOH Solution 1.54 N	2 mL	1 bottle
HI 93763B-0	Molybdovanadate Reagent	0.5 mL	1 bottle
PERSULFATE/P	Potassium Persulfate	1 packet	50 packets

* *Reagent Vial identification: 34xx (xx represents the reagent lot code). GREEN RECTANGLE ON THE LABEL.*

Note: Store the unused vials in their container in a cool and dark place.

REAGENT SET

HI 94763B-50 Reagents for up to 49 tests

For other accessories see page 70.

REQUIRED ACCESSORIES

HI 839800-01 Hanna reactor (115 VAC)

HI 839800-02 Hanna reactor (230 VAC)

HI 740216 test tube cooling rack (25 holes)

HI 740217 Laboratory bench safety shield

For other accessories see page 70.

MEASUREMENT PROCEDURE




Before using the reagent kit carefully read all the instructions and the Material Safety Data Sheet (MSDS). Pay particular attention to all warnings, cautions and notes. Failure to do so may result in serious injury to the operator.

Reagent Blank Correction: This method requires a reagent blank correction. A single blank vial may be used more than once; the blank vial is stable up to one day (room temperature). For most accurate measurement, run a blank for each set of measurements and always use the same lot of reagents for blank and samples.


- Preheat the Hanna Reactor HI 839800 to 150 °C (302 °F). For correct use of the reactor follow Reactor Instruction Manual. The optional HI 740217 safety shield is strongly recommended.

DO NOT USE AN OVEN OR MICROWAVE samples may leak and generate a corrosive and possibly explosive atmosphere.


- Remove the cap from two Reagent Vials.



- Add exactly 5.0 mL of sample to one vial (sample vial), and 5.0 mL of deionized water to the other vial (blank vial), while keeping the vials at a 45-degree angle.



- Add the content of one packet of Potassium Persulfate for Phosphorus analysis to each vial. Replace the cap tightly and shake gently the vials until all the powder is completely dissolved.



- Insert the vials into the reactor and heat them for 30 minutes at 150 °C.




- At the end of the digestion place the vials carefully in the test tube rack and allow to cool to room temperature.

Warning: The vials are still hot, be careful when handling them.




Phosphorus, Total HR


- Select the **Phosphorus, Total High Range** method following one of the procedures described in the “Method Selection” section (see page 15).
- Remove the cap from the vials and add exactly 2.0 mL of Sodium Hydroxide (NaOH) Solution 1,54 N to each vial, while keeping the vials at a 45-degree angle. Replace the cap tightly and mix by inverting the vials a couple of times.





- Remove the cap from the vials and add exactly 0.5 mL of HI 93763B-0 Molybdovanadate Reagent to each vial, while keeping the vial at a 45-degree angle. Replace the cap tightly and mix by inverting the vials a couple of times.






- Place the blank vial into the holder and push it completely down.
- Press **Timer** and the display will show the countdown prior to the measurement and the “Reaction Time” message. Alternatively, wait for 7 minutes and press **Zero**.



- Wait for vial identification. If that was successfully done, the instrument will perform a zero sequence and after a few seconds the display will show “-0.0-”. Now the meter is zeroed and ready for measurement.





- Remove the blank vial.
- Place the sample vial into the holder and push it completely down.
- Press **Read** and wait for vial identification. If that was successfully done, the instrument will perform the reading. The instrument displays concentration in **mg/L of phosphorus (P)**.



The method detects free (orthophosphate), condensed inorganic forms (meta-, pyro- and other polyphosphates) and organic forms of phosphates present in the sample.

- Press **▲** or **▼** to access the second level functions and then press the **Chem Frm** key to convert the result in mg/L of phosphate (PO_4^{3-}) and mg/L of P_2O_5 .



- Press **▲** or **▼** to return to the measurement screen.

Note: For accurate measurements

- 1) wash glassware only with phosphate-free detergents
- 2) clean all glassware with 1 : 1 hydrochloric acid solution and rinse with deionized water.

INTERFERENCES

Arsenate

pH: the sample should have a neutral pH

Temperature: the method is temperature sensitive.

It is recommended to add the Molybdovanadate Reagent and to run measurements at $T = 20$ to $25\text{ }^{\circ}\text{C}$:


$T < 20\text{ }^{\circ}\text{C}$ causes a negative error

$T > 25\text{ }^{\circ}\text{C}$ causes a positive error

Turbidity and suspended matter in large amounts may cause interference because the strongly acidic reaction conditions may dissolve suspended matter or cause desorption of phosphates from particles. Before measurement, turbidity or suspended matter should be removed by treatment with active carbon and by prior filtration.

ERRORS AND WARNINGS

The instrument shows clear messages when erroneous conditions appear and when measured values are outside the expected range. These messages are described below.



Method identification was unsuccessful: The instrument cannot identify the vial's bar code or the vial has no code.



Wrong vial: The vial's bar code is not the expected one.




No Light: The light source is not functioning properly.




Light Leak: There is too much light when dark current is measured.



Inverted vials: The sample and the zero vials are inverted.




Wrong vial: The current method doesn't correspond to the current vial.



Different reagent lot: The reagent lot doesn't correspond to the current lot code.

DATA MANAGEMENT

The analyzed data can be managed using Hanna's product HI92000, Windows® Compatible Software.



STANDARD METHODS

Description	Range	Method
Ammonia LR	0.00 to 3.00 mg/L	Nessler
Ammonia HR	0 to 100 mg/L	Nessler
Chlorine, Free	0.00 to 5.00 mg/L	DPD
Chlorine, Total	0.00 to 5.00 mg/L	DPD
Nitrate	0.0 to 30.0 mg/L	Chromotropic Acid
Nitrogen, Total LR	0.0 to 25.0 mg/L	Chromotropic Acid
Nitrogen, Total HR	10 to 150 mg/L	Chromotropic Acid
COD LR	0 to 150 mg/L	Dichromate, Mercuric Sulfate
COD MR	0 to 1500 mg/L	Dichromate, Mercuric Sulfate
COD HR	0 to 15000 mg/L	Dichromate, Mercuric Sulfate
Phosphorus, Reactive	0.00 to 1.60 mg/L	Ascorbic Acid
Phosphorus, Acid Hydrolyzable	0.00 to 1.60 mg/L	Ascorbic Acid
Phosphorus, Total LR	0.00 to 1.15 mg/L	Ascorbic Acid
Phosphorus, Reactive HR	0.0 to 32.6 mg/L	Vanadomolybdophosphoric Acid
Phosphorus, Total HR	0.0 to 32.6 mg/L	Vanadomolybdophosphoric Acid

Windows® is registered Trademark of "Microsoft Co."

ACCESSORIES

REAGENT SETS

HI 93701-01	100 free chlorine tests
HI 93701-03	300 free chlorine tests
HI 93711-01	100 total chlorine tests
HI 93711-03	300 total chlorine tests
HI 94754A-25	25 COD LR tests
HI 94754B-25	25 COD MR tests
HI 94754C-25	25 COD HR tests
HI 94758A-50	50 reactive phosphorus tests
HI 94758B-50	50 acid hydrolyzable phosphorus tests
HI 94758C-50	50 total phosphorus LR tests
HI 94763A-50	50 reactive phosphorus HR tests
HI 94763B-50	50 total phosphorus HR tests
HI 94764A-25	25 ammonia LR tests
HI 94764B-25	25 ammonia HR tests
HI 94766-50	50 nitrate tests
HI 94767A-50	50 total nitrogen LR tests
HI 94767B-50	50 total nitrogen HR tests

OTHER ACCESSORIES

HI 839800-01	Hanna Reactor (115 VAC)
HI 839800-02	Hanna Reactor (230 VAC)
HI 731318	Cloth for wiping vials (4 pcs.)
HI 731340	200 µL automatic pipette
HI 731341	1000 µL automatic pipette
HI 731342	2000 µL automatic pipette
HI 731350	tips for 200 µL automatic pipette (25 pcs.)
HI 731351	tips for 1000 µL automatic pipette (25 pcs.)
HI 731352	tips for 2000 µL automatic pipette (4 pcs.)
HI 740142	1 mL graduated syringe
HI 740143	1 mL graduated syringe (6 pcs.)
HI 740226	5 mL graduated syringe
HI 740144	Pipette tip (6 pcs.)
HI 740157	Plastic refilling pipette (20 pcs.)
HI 740216	Test tube cooling rack (25 holes)
HI 740217	Laboratory bench safety shield
HI 92000	Windows® Compatible Software
HI 920013	PC Connection Cable
HI 93703-50	Vial cleaning solution (230 mL)