



Summit

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English – page 3

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Contents

Warnings and Safety Precautions	4
Getting Started	5
Installation	6
Below-Balance Weighing	10
Leveling the Balance	11
Operation	12
Display and Operating Elements	12
Basic Weighing Function	13
Calibration/Adjustment	15
Configuration (Operating Menu)	18
Functions of the Keys during Configuration	18
Menu Navigation; Example: Setting the Language	19
Parameter Settings: Menu	20
Parameter Settings: Overview	21
Input: ID Number, Date and Time	25
Application Programs	27
Counting	28
Weighing in Percent	30
Calculation	32
Animal Weighing/Averaging	34
Net-total Formulation	36
Totalizing	38
Density Determination	40
Mass Unit Conversion	42
ISO/GLP-compliant Printout/Record	45
Data Interface	47
Troubleshooting Guide	48
Care and Maintenance	49
Recycling	50
Overview	51
Specifications	52
Accessories	54
CE Marking	55

Warnings and Safety Precautions

Safety

- To prevent damage to the equipment, please read these operating instructions carefully before using the balance.
 -  Do not use this equipment in hazardous areas.
 -  The balance may be opened only by trained service technicians.
 -  Disconnect the balance from power before connecting or disconnecting peripheral devices.
 -  If you operate the balance under ambient conditions subject to higher safety standards, you must comply with the applicable installation regulations.
 -  Exposure to excessive electromagnetic interference can cause the readout value to change. Once the disturbance has ceased, the instrument can be used again in accordance with its intended purpose.
Make sure that no liquid enters the equipment housing; use only a slightly moistened cloth to clean the balance.
 -  Make sure the voltage rating printed on the power supply is identical to your local line voltage.
 - Proceed with extreme caution when using pre-wired RS-232 connecting cables, as the pin assignments may not be compatible with Denver equipment.
Before connecting the cable, check all pin assignments against the cabling diagrams and disconnect any lines that are assigned differently.
-  If there is visible damage to the equipment or power cord, disconnect the equipment from power and lock it in a secure place to ensure that it cannot be used for the time being.
- Connect only Denver accessories, as these are optimally designed for use with your Summit balance. The operator shall be solely responsible for installation and testing of any modifications to Denver equipment, including connection of cables or equipment not supplied by Denver. On request, Denver will be happy to provide information on operating specifications (in accordance with the Standards for defined immunity to interference).
- Do not open the balance housing. If the seal is broken, this will void all claims under the manufacturer's warranty.
- If you have any problems with your balance, contact your local Denver customer service center.
- ## Symbols
- The following symbols are used in these instructions:
- indicates required steps
 - indicates steps required only under certain conditions
 - > describes what happens after you have performed a particular step
 - indicates an item in a list
 -  indicates a hazard

Getting Started

Storage and Shipping Conditions

- Do not expose the balance to extreme temperatures, moisture, shocks, blows or vibration.

Unpacking the Equipment

- After unpacking the equipment, please check it immediately for any external damage.
- If you detect any damage, proceed as directed in the chapter entitled "Care and Maintenance," under "Safety Inspection."
- Save the box and all parts of the packaging for any future transport. Disconnect all cables before packing the balance for shipping.

Equipment Supplied

- Balance
- Weighing pan
- Pan support (only for models with a round weighing pan)
- AC adapter

Additional equipment supplied with models S/SI-234, -114, -64:

- Sliding-door draft shield
- Drip/breeze ring
- Draft shield base plate

Additional equipment supplied with models S/SI-603, -403, -203:

- Round glass draft shield with cover

Installation

Choose a location that is not subject to the following negative influences:

- Heat (heater or direct sunlight)
- Drafts from open windows and doors
- Excessive vibration during weighing
- Excessive moisture

Conditioning the Balance


Moisture in the air can condense on cold surfaces whenever the equipment is moved to a substantially warmer place. To avoid the effects of condensation, allow the balance to sit for 2 hours, at room temperature, before plugging into AC power.

Seal on Balances Verified for Use in Legal Metrology in the EU*:

EU legislation requires that a control seal be affixed to verified balances. The control seal consists of a sticker. If the seal is broken, the verification becomes null and void and the balance must be re-verified.

* Including the Signatories of the Agreement on the European Economic Area


Installation

Setting Up the Balance


Instruments with a sliding-door draft shield:

- Place components inside the chamber in the following order:
 - Draft shield base plate
 - Drip/breeze ring
 - Pan support
 - Weighing pan

Instruments with a round glass draft ring:

- Position the components listed below in the order given:
 - Place the lower lid on the balance with the raised edge facing upwards and turn it until it is firmly in position
 - Pan support
 - Weighing pan
 - Glass draft ring
 - Place the upper lid on the draft shield ring with the raised edge facing downwards


Instruments with a round weighing pan

- Position the components listed below in the order given:
 - Pan support
 - Weighing pan

Instruments with a rectangular weighing pan:

- Place the weighing pan on the balance


Connecting the Balance to AC Power/Safety Precautions

Use only original Denver AC adapters.

For use within

- Europe: part no. 6971412
- U.S./Canada: part no. 6971413
- Connect the angle plug to the balance
- Connect the AC adapter to the wall outlet (mains)

AC Adapter with Country-specific Power Cord

Some models come with separate country-specific power cords for the AC adapter.

- Connect the angle plug to the balance
- Select the power cord for your area and connect it to the AC adapter
- Plug the power cord into the wall outlet (mains)
- Use an original Denver AC adapter with a wide input voltage range (100 to 240 V~), order no. 6971966, and replaceable power cord:

6900900 (Europe)	6900905 (AUS)
6900901 (US/CDN)	6900902 (ZA)
6971945 (UK)	6971977 (Argentina)
6971973 (India)	6971978 (China)
6971980 (Denmark)	6971975 (Israel)
6971776 (Italy)	

Safety Precautions

Plug-in AC Adapter 6971412/6971413:

The AC adapter rated to Class 2 can be plugged into any wall outlet without additional safety precautions.

Benchtop AC Adapter 6971966:

The AC adapter rated to Class 1 can be plugged into any wall outlet without additional safety precautions.

The ground terminal is connected to the balance housing, which can be additionally grounded for operation. The data interface is also electrically connected to the balance housing (ground).

NOTE: This equipment has been tested and found to comply with the limits pursuant to part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with these instructions, may cause harmful interference to radio communications.

For information on the specific limits and class of this equipment, please refer to the Declaration of Conformity. Depending on the particular class, you are either required or requested to correct the interference.


If you have a Class A digital device, you need to comply with the FCC statement as follows: "Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense."

If you have a Class B digital device, please read and follow the FCC information given below:

"However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:


- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help."

Before you operate this equipment, check which FCC class (Class A or Class B) it has according to the Declaration of Conformity included. Be sure to observe the information of this Declaration.

Connecting Electronic Peripheral Devices

- Make sure to unplug the balance from AC power before you connect or disconnect a peripheral device (printer or computer) to or from the interface port.

Warmup Time


To ensure accurate results, the balance must warm up for 30 minutes before operation. Only after this time will the instrument have reached the required operating temperature.

Using Verified Balances in Legal Metrology in the EU*:

- Make sure to allow the equipment to warm up for at least 24 hours after initial connection to AC power or after a relatively long power outage.

Operation Outside the Temperature Range

⚠ Operating the balance beyond the temperature range of +10...+30°C (50°...86°F). Differences from the specifications listed in the chapter on "Specifications" are possible.


Anti-theft Locking Device

To attach an anti-theft locking device (chain or similar), use the lug for attachment on the back of the balance.

- Secure the balance at the place of installation by fastening a chain or lock.

* Including the Signatories of the Agreement on the European Economic Area


Below-Balance Weighing

A port for a below-balance weighing hanger is located on the bottom of the balance.

- Below-balance weighing is not permitted in legal metrology.

- Open cover plate on the bottom of the balance.


Important: set the balance on its side to access the cover plate. DO NOT turn the balance upside-down.

- Using the built-in hook 1: Attach the sample (e.g., using a suspension wire) to the hanger.

- Bore hole 2 (not on models S/SI-203, -602, -402, -8001, -4001): Carefully fasten the special hanger, or order a hanger directly from Denver.

- Install a shield for protection against drafts if necessary.


Leveling the Balance

Purpose:

- To compensate for unevenness at the place of installation

Always level the balance again any time after it has been moved to a different location. Only the 2 front feet are adjusted to level the balance.

- Retract the two rear feet (only on models with a rectangular weighing pan).
 - Turn the 2 front feet as shown in the diagram until the air bubble is centered within the circle of the level indicator.
- > In most cases this will require several adjustment steps.
- On models with a rectangular weighing pan: Lower the 2 rear feet until they touch the surface on which the balance rests.


Anti-theft Locking Device

- To secure the balance at the place of installation, fasten a chain or a lock to the lug located on the rear panel of the balance.

Operation

Overview of Display Elements

Pos. Designation


- | | |
|----|--|
| 1 | Weight unit |
| 2 | Menu level indicator |
| 3 | Symbol: "GLP printing mode active" |
| 4 | Symbol: "Printing mode active" |
| 5 | Symbol: "Application program active" |
| 6 | Calculated-value indicator
(i.e., not a weight value) |
| 7 | Symbol: Gross or net value |
| 8 | Symbols for active application
(Δ , ∇ , %, \odot , \pm , A, C) |
| 9 | Symbol: Calibration/adjustment function |
| 10 | Symbols for zero range |

Pos. Designation

- | | |
|----|--|
| 11 | Busy symbol: command is being processed
(for example, "Wait for stability icon);
after you turn on the power, \diamond will be
displayed until you press a key |
| 12 | Weight value displayed in selected weight
unit

Symbols:
<< Save settings and exit the operating menu
< One menu level higher
V Scroll through menu items
> Next item on current menu level
↓ Select a parameter setting |

Overview of Operating Elements

Pos. Designation

- | | |
|---|--|
| 1 | Taring |
| 2 | Start an application program |
| 3 | Select an application program Open the
operating menu |
| 4 | Clear
This key is generally used to cancel functions:
– Quit application program
– Cancel calibration/adjustment routine
Exit the operating menu |

Pos. Designation

- | | |
|---|--|
| 5 | Start calibration/adjustment routine |
| 6 | Data output:
Press this key to send readout values to
the built-in data interface. |
| 7 | On/off |
| 8 | Level indicator |

Basic Weighing Function

Preparation

- Switch on the balance: Press 
 - Tare the balance, if necessary:
Press 
 - If necessary, change the configuration settings:
See the chapter entitled
“Configuration”
 - If desired, load the factory settings:
See the chapter entitled
“Configuration”
- Additional Functions
- Switching off the balance:
Press 
 - Balance in standby mode:
the current time is displayed

Features

- Taring the balance
- Printing weights

Using Verified Balances as Legal Measuring Instruments in the EU*:


The type-approval certificate for verification applies only to non-automatic weighing instruments. For automatic operation with or without auxiliary measuring devices, you must comply with the regulations applicable to the place of installation.

- Before using the balance as a legal measuring instrument, calibrate and adjust it at the place of use using the built-in motorized calibration weight; for details, see “Calibration/Adjustment” in this chapter.
 - The temperature range (°C) indicated on the verification label must not be exceeded during operation.
- Example:
BD ED 200
 +10°C to +30°C

* Including the Signatories of the Agreement on the European Economic Area

Example

Simple Weighing

Step	Key (or instruction)	Display/Printout
Balance in standby mode		11:54 AM
1. Switch on the balance Self-test is performed, followed by automatic initial tare function.	①	0.0 g
2. Place container on weighing pan (in this example: 11.5 g).		+ 11.5 g
3. Tare the balance	 TARE	0.0 g
4. Place sample in container (in this example: 132 g).		+ 132.0 g
5. Print weight.	 PRINT	N + 132.0 g

Calibration and Adjustment

Purpose

Calibration is the determination of any difference between the measured value displayed and the true weight (mass) of a sample. Adjustment is the correction of this difference, or its reduction to an allowable level within maximum permissible error limits.

Using Verified Balances as Legal Measuring Instruments in the EU*: Before using your balance as a legal measuring instrument, internal calibration must be performed at the place of installation.

Features

Calibration/adjustment can be performed only when:

- there is no load on the balance
- the balance is tared
- the internal signal is stable
- for external calibration, the value displayed for the calibration weight on the balance does not differ from the nominal weight value by more than 2%

If these conditions are not met, an error message is displayed ("ERR 02").

You can use any of the following weight units in calibration/adjustment: **CAL.UNIT: GRAMS, KILOGR. or POUNDS** (not for verified models)


External Calibration in Verified Balances

When the balance is used in legal metrology, external calibration is blocked as follows:

- The setting of the menu access switch is locked (see "To block calibration/adjustment")
- The cap over the menu access switch is sealed

To block calibration/adjustment:

- Select **CAL.-ADJ.: BLOCKED** in the menu
- Close the menu access switch on the back of the balance

For details on generating an ISO/GLP-compliant printout of calibration/adjustment results, see page 45.

Following calibration/adjustment, the application program is cleared.

Internal Calibration/Adjustment

In the operating menu, select **CAL.-ADJ.: CAL.INT.** before beginning.

The built-in motorized calibration weight is applied and removed automatically for internal calibration.

- Select calibration/adjustment:

- Press **CAL**
- > The built-in weight is applied automatically
- > The balance is adjusted
- > The built-in calibration weight is removed.

* Including the Signatories of the Agreement on the European Economic Area

Internal Calibration/Adjustment (Only on SI models)

Set the following parameters:

SETUP: BAL.SCAL.: CAL.-ADJ.: CAL.INT. (menu code 1.1.9.4)

The built-in motorized calibration weight is applied and removed automatically for internal calibration.



Step	Key (or instruction)	Display
1. Tare the balance	(TARE)	0.0 g
2. Start calibration	(CAL)	CAL.INT.
	The built-in weight is applied automatically	CAL.RUN.
3. Calibration/adjustment executed		CAL.END
4. The built-in weight is removed		0.0 g

External Calibration

Parameters (changes in factory settings):

SETUP: BAL.SCAL.: CAL.-ADJ.: CAL.EXT. (menu code 1.1.9. 1)

The required calibration weight is configured at the factory (see "Specifications")

Step	Key (or instruction)	Display
1. Tare the balance	(TARE)	0.0 g
2. Start calibration.	(CAL)	CAL.EXT.
	Once you store the zero point the required calibration weight is prompted (flashing display)	- 5000.0 g
3. Apply the prompted calibration weight (in this example: 5000 g) Weight too light: a minus sign “-” is shown Weight too heavy: a plus sign “+” is shown		5000.0 g
	The display stops flashing as soon as the weight value is within the defined limit.	
4. Calibration/adjustment executed; then the calibration weight is displayed		CAL.END
		+ 5000.0 g
5. Remove the calibration weight		0.0 g

Configuration (Operating Menu)

You can configure the balance; i.e., adapt it to individual requirements.

Functions of the Keys during Configuration

Symbol	Key	Function
V	Press and hold	Scroll through menu items
>		One menu level lower
↓		Confirm menu item
	Press and hold	Save settings and exit menu from any position
<<		Save settings and exit menu
<		One menu level higher
[****]		Indicates menu level

Menu Navigation

Example: Setting the Language

Step	Key (or instruction)	Display
1. Open the menu: In weighing mode: first menu item is shown	(SELECT) (hold)	APPL IC.
2. Scroll upward within the menu level; after the last menu code, the first code is displayed again	Repeatedly: (SELECT)	INPUT ... LANGUAGE.
3. Select menu level (scrolls to the right)	(ENTER)	ENGLISH °
4. Change setting: Scroll until the desired setting is shown	(SELECT)	ESPAÑOL
5. Confirm the menu code; “o” indicates the active setting	(ENTER)	ESPAÑOL °
6. Return to the next higher menu level (from the second level) <input type="radio"/> Set other menu items as desired	(CLEAR)	LENGUA
7. Save settings and exit menu	Repeatedly: (CLEAR)	
or		
<input type="radio"/> Exit menu without saving changes	(POWER)	
> Restart your application		0.0 g

Parameter Settings: Menu

Level 1 [•]	Level 2 [••]	Level 3 [•••]	Menu code
SETUP	<ul style="list-style-type: none"> BAL.SCAL. balance parameters INTERF. Interface PRNT.OUT Settings for print function EXTRAS Additional functions RESET 	<ul style="list-style-type: none"> AMBIENT Ambient conditions APP.FILT. Application filter STAB.RNG. Stability range TARING Taring1) AUTOZER. Auto zero WT.UNIT Basic weight unit DISPLAY Display accuracy1) CAL./ADJ. Function of the key CAL.UNIT Weight unit for calibration1) BAUD Baud rate PARITY Parity STOPBIT Number of stop bits HANDSHK Handshake mode DATBITS Number of data bits DAT.REC. Output: SBI (ASCII) or printout PRINT (manual/automatic) STOPAUT. Stop automatic printing TAR./PRT. Tare balance after ind. print PRT.INIT. Printout of appl. parameters FORMAT Line format for printout GLP ISO/GLP-compliant printout TIME: 12/24 h DATE: Format MENU SIGNAL Acoustic signal (beep) KEYS Keypad EXT.KEY External switch function ON/MODE Power-on mode BACRLIT Display backlighting MENU Factory settings 	1.1.1. 1.1.2. 1.1.3. 1.1.5. 1.1.6. 1.1.7. 1.1.8. 1.1.9. 1.1.11. 1.5.1. 1.5.2. 1.5.3. 1.5.4. 1.5.5. 1.5.6. 1.6.1. 1.6.2. 1.6.4. 1.6.5. 1.6.6. 1.6.7. 1.6.8. 1.6.9. 1.8.1. 1.8.2. 1.8.3. 1.8.4. 1.8.5. 1.8.6. 1.9.1.
APPLIC. Application programs	<ul style="list-style-type: none"> WEIGH UNIT Toggle wt. unit COUNT Counting PERCENT Weighing in percent NET-TOT Net-total formulation TOTAL Totalizing ANIMALW. Animal weighing CALC. Calculation DENSITY Density determination 	<ul style="list-style-type: none"> DISP.DIG. Display accuracy1) RESOLUT. Resolution REF.UPDATE Autom. ref. sample updating DEC.PLCS Decimal places COMP.PRT. Printout of components COMP.PRT. Printout of components ACTIVITY. Animal activity START METHOD (operator) DEC.PLCS Decimal places DEC.PLCS Decimal places 	2.1. 2.2.2. 2.3.1. 2.3.2. 2.4.1. 2.5.1. 2.6.1. 2.7.1. 2.7.2. 2.8.1. 2.8.2. 2.9.1.
INPUT Input	ID NO., DATE, TIME	Input: ID no., date, time	3.1/2/3.
INFO Information	VERSION, SER.NO., MODEL	Display software ver., serial no., model	4.1/2/3.
LANGUAG.	<ul style="list-style-type: none"> ENGLISH (factory setting) DEUTSCH German FRANCE French ITAL. Italian ESPAÑOL Spanish CODES Menu shows codes (not texts) 		5.1. 5.2. 5.3. 5.4. 5.5. 5.6.

¹⁾ Setting cannot be changed on verified balances

Parameter Settings: Overview

o = Factory setting √ = User-defined setting

Level 1 [•]	Level 2 [••]	Level 3 [•••]	Level 4 [••••]	Menu code
SETUP —	BAL.SCAL. balance parameters	AMBIENT — Ambient conditions (Filter adaptation)	V-STABLE Very stable o STABLE UNSTABL.	1.1.1.1 1.1.1.2 1.1.1.3 1.1.1.4
		APP.FILT. Application filter	o FINAL RD Final readout mode FILLING Filling mode	1.1.2.1 1.1.2.2
		STAB.RNG. Stability range	1/4 DIG. (digit) 1/2 DIG. 1 - DIGIT o 2 - DIGIT 4 - DIGIT 8 - DIGIT	1.1.3.1 1.1.3.2 1.1.3.3 1.1.3.4 1.1.3.5 1.1.3.6
		TARING ¹⁾ Taring	W/O STAB W/o stability o W/ STAB After stability	1.1.5.1 1.1.5.2
		AUTOZER. Auto zero	OFF o ON	1.1.6.1 1.1.6.2
		WT.UNIT Basic weight unit	For list of units, see "Toggling between Weight Units"	1.1.7.1 1.1.7.23
		DISP.DIG. ¹⁾ Display accuracy	o ALL MINUS 1	1.1.8.1 1.1.8.2
		CAL./ADJ. Function of the CAL key	o CAL.EXT. External cal./adj. ¹⁾ CAL. INT. Internal cal./adj. ²⁾ BLOCKED key blocked	1.1.9.1 1.1.9.2 1.1.3.3
		CAL.UNIT. Unit ¹⁾ for calibration weight	o GRAMS KILOGR. Kilograms POUNDS	1.1.11.1 1.1.11.2 1.1.11.3

¹⁾) Setting cannot be changed on verified balances

²⁾) Only on SI models

Level 1 [•]	Level 2 [••]	Level 3 [•••]	Level 4 [••••]	Menu code
SETUP	INTERF. Interface	BAUD rate	600 1200 2400 4800 9600 19200	1.5.1.3 1.5.1.4 1.5.1.5 1.5.1.6 1.5.1.7 1.5.1.8
		PARITY	o ODD EVEN NONE	1.5.2.3 1.5.2.4 1.5.2.5
		STOPBIT	o 1 BIT 2 BITS	1.5.3.1 1.5.3.2
		HANDSHK.	SFTWARE HARDWARE NONE	1.5.4.1 1.5.4.2 1.5.4.3
		DATABIT	o 7 BITS 8 BITS	1.5.5.1 1.5.5.2
		DAT.REC. Com- munication mode	o SBI (ASCII) ¹⁾ PRINTER (GLP-printout)	1.5.6.1 1.5.6.2
	PRNT.OUT Printing fct.	PRINT	MAN. W/O W/o stability MAN.WITH W/ stability AUT.W/O Autom. w/o stability AUT.WITH. Autom. w/ stability	1.6.1.1 1.6.1.2 1.6.1.3 1.6.1.4
		STOPAUT. Stop automatic printing	o OFF Not possible ON Use print key PRINT	1.6.2.1 1.6.2.2
		TAR./PRT. Tare the balance after individual printout	o OFF ON	1.6.4.1 1.6.4.2

- 1) Note concerning verified balances as legal measuring instruments in the EU*:
 In the setting "SBI", the non-verified display digit is not automatically identified.
 Please take the corresponding measures or adjust the settings on the peripheral device.

* Including the signatories of the Agreement on the European Economic Area.

Level 1 [•]	Level 2 [••]	Level 3 [•••]	Level 4 [••••]	Menu code
SETUP —	PRNT.OUT Printing	PRT.INIT. Printing application parameters	OFF o ALL All parameters MAINPAR. Main parameters	1.6.5.1 1.6.5.2 1.6.5.2
		FORMAT Line format for printout	16CHAR. 16 characters (w/o ID) 22CHAR. 22 characters (w/ ID) 2CHAR. 2 characters (w/ ID) (Date/Time and weight value)	1.6.6.1 1.6.6.2 1.6.6.3
		GLP Printout as ISO/GLP-compliant	OFF CAL.-ADJ. Only for calib./adj. ALWAYS All printouts	1.6.7.1 1.6.7.2 1.6.7.3
		TIME	24H 24-hour format 12H 12-hour format "AM/PM"	1.6.8.1 1.6.8.2
		DATE	DD/MMM/YY Day/month/year MMM/DD/YY Month/day/year	1.6.9.1 1.6.9.2
	EXTRAS Additional functions	MENU	CAN EDIT Can change settings R.D. ONLY Read only	1.8.1.1 1.8.1.2
		SIGNAL Acoustic signal	OFF ON	1.8.2.1 1.8.2.2
		KEYS Keypad	FREE LOCKED	1.8.3.1 1.8.3.2
		EXT.KEY Function of the external switch	PRINT (PRINT) Z/TARE (TARE) CAL. (CAL) SELECT (SELECT) CF (CLEAR) ENTER (ENTER) BLOCKED: External key blocked	1.8.4.1 1.8.4.2 1.8.4.3 1.8.4.4 1.8.4.5 1.8.4.6 1.8.4.9
		ON MODE Power-on mode	OFF/ON Off/on/standby STANDBY On/standby AUTO ON Auto on	1.8.5.1 1.8.5.2 1.8.5.3
		BACKLIT Display backlighting	OFF ON	1.8.6.1 1.8.6.2
RESET Reset menu	MENU Factory settings	YES Restore defaults NO Do not restore defaults		1.9.1.1 1.9.1.2


Level 1 [•]	Level 2 [••]	Level 3 [•••]	Level 4 [••••]	Menu code
APPL.IC. Applic. programs	WEIGH			2.1.
	UNIT	DISP.DIG.1)	o ALL	2.2. 2. 1
	Toggle units	Display accuracy	MINUS 1	2.2. 2. 2
	COUNTING	RESOLUT.	o DISP.ACC. Display accuracy	2.3. 1. 1
		Resolution	10-FOLD 10 times > disp.	2.3. 1. 2
		REF.UPDT.	o OFF	2.3. 2. 1
		Autom. reference updating	AUTO	2.3. 2. 2
	PERCENT	DEC.PLCS	NONE No dec. places	2.4. 1. 1
	Weighing in percent	Decimal places	1 DEC.PL. 1 decimal place	2.4. 1. 2
			2 DEC.PL. 2 decimal places	2.4. 1. 3
			3 DEC.PL. 3 decimal places	2.4. 1. 4
NET-TOT	COMP.PRT.		OFF	2.5. 1. 1
	Net-total	Component printout	o ON	2.5. 1. 2
TOTAL	COMP.PRT.		OFF	2.6. 1. 1
	Totalizing	Component printout	o ON	2.6. 1. 2
ANIMALW. Animal weighing	ACTIVITY.		CALM (fluct.: 2% of test obj.)	2.7. 1. 1
		Animal activity	o ACTIVE (fluct.: 5% of test obj.)	2.7. 1. 2
			V.ACTIVE (fluct.: 20% of test obj.)	2.7. 1. 3
	START		MANUAL	2.7. 2. 1
			o AUTO. Automatic	2.7. 2. 2
CALC. Calculation	METHOD		o MUL. Multiplier	2.8. 1. 1
	(operator)		o DIV. Divisor	2.8. 1. 2
	DEC.PLCS		NONE No dec. places	2.8. 2. 1
	Decimal places	o 1 DEC.PL. 1 decimal place	2.8. 2. 2	
			2 DEC.PL. 2 decimal places	2.8. 2. 3
			3 DEC.PL. 3 decimal places	2.8. 2. 4
DENSITY	DEC.PLCS		NONE No dec. places	2. 9. 1. 1
	Density determination	Decimal places	o 1 DEC.PL. 1 decimal place	2. 9. 1. 2

¹⁾ Setting cannot be changed on verified balances


Input: ID number, Date and Time

	Level 1 [•]	Level 2 [••]	Level 3 [•••]	Menu Code
INPUT	ID NO.		ID input for ISO/GLP compliant data record; 7 characters max. Permitted characters: 0 to 9; A to Z; dash/hyphen; space	3. 1.
Input	DATE		Menu item for setting the date	3. 2.
	TIME		Menu item for setting the time	3. 3.

Depending on the setting you select for the menu item "SETUP PRNT.OUT DATE," the date will be displayed in the following format:

Format	Display: Date
DD.MMM.YY	
MMM.DD.YY	

Depending on the setting you select for the menu item "SETUP PRNT.OUT TIME," the time will be displayed in the following format:

Time	Display: Time
24-hour format	
12-hour format	

Example: ID No., Date and Time

Step	Key (or instruction)	Display
1. Open the menu: In weighing mode; first menu item is displayed	(SELECT) hold	APPLIC.
2. Select "Input"	(SELECT)	INPUT
3. Select input for ID no.	(ENTER) twice	ID NO.
4. Set or change the ID no. – hold down key to automatically change the digit(s):	(SELECT)	3-----
5. Scroll within the 7-digit ID no.	(SELECT) hold (ENTER) or (CLEAR)	3-ABC 12
6. Save input when you have reached the last digit of the ID no.	(ENTER)	ID NO.
7. Select "Date"	(SELECT), (ENTER)	FEB. 08.
8. Change setting – hold down key to automatically change the digit(s):	(SELECT)	FEB. 10.
9. Toggle between Day/Month/Year positions	(SELECT) hold (ENTER) or (CLEAR)	FEB. 10.
10. Save setting when you reach the "YEAR" position	(ENTER)	DATE
11. Select "Time"	(SELECT), (ENTER)	A 10.46.23
12. Change setting – hold down key to automatically change the digits:	(SELECT) (SELECT) hold	A 11.46.23
13. Toggle between Hour/Minute/Second positions	(ENTER) or (CLEAR)	A 11.46.32
14. Set seconds to zero	(SELECT)	A 11.47.00
15. Save setting when you have reached the "second" position	(ENTER)	TIME
16. Save all settings and exit the menu	(CLEAR) repeatedly	
> Restart your application		0.0 g

Device Information

Level 1 [•]	Level 2 [••]	Level 3 [•••]	Example	Menu code
INFO Information	VERSION	Show software version	REL.36.0 1	4. 1.
	SER. NO.	Show serial number (To toggle focus between upper and lower display sections, press (SELECT))	1080 1234	4. 2.
	MODEL	Show model designation (to change focus from upper to middle to lower display section and back, press (SELECT))	SI-800 1	4. 3.

Display of Menu Items: Text or Codes

LANGUAGE	ENGLISH (factory setting)	5. 1.
	DEUTSCH German	5. 2.
	FRANC. French	5. 3.
	ITAL. Italian	5. 4.
	ESPAÑOL Spanish	5. 5.
	CODES Menu shows codes (not texts)	5. 6.

Application Programs

Using Verified Balances as Legal Measuring Instruments in the EU*:
All application programs can be selected on balances used as legal measuring instruments.

Calculated values are alternately indicated with the following symbols:

- Percent = %
- Piece count (Counting) = pcs
- Computed value = o,Δ

* Including the Signatories of the Agreement on the European Economic Area

Counting

Display symbol: 


Purpose

With the Counting program you can determine the number of parts that each have approximately equal weight. To do this, a known number of parts (the reference sample quantity) is weighed first, and the individual piece weight (reference weight) is calculated from this result. Thus the number of parts subsequently placed on the balance can be determined from their weight.

Preparation

- Select the Counting application in the menu:
see "Configuration."
- Set the following parameters:

APPLIC. Application program

o = Factory setting

Changing the Reference Sample Quantity

Activate function:

Press the  key

Select the desired reference sample quantity (1 to 100):

In increments of 1: Press the  key briefly

In increments of 10:

Press and hold the  key.

The quantity is stored in battery-backed memory.

Reference Sample Updating

Automatic reference sample updating optimizes the counting accuracy. You can activate or deactivate this function in the menu.

Automatic reference sample updating is performed when the requirements, including the specified stability criterion, have been met.





The abbreviation  for "optimizing", is displayed briefly with the new reference sample quantity.

Printout: Counting

nRef	+	10	: Reference sample quantity
wRef	+	21.14 g	: Reference weight
Qnt	+	500 pcs	: Calculated quantity

Example : Counting parts of equal weight

Parameter: APPLIC. - COUNT. (menu code 2. 3.)

Step	Key (or instruction)	Display/Data output
1. Place empty container on the balance		+ 22.6 g
2. Tare the balance		0.0 g
3. Add reference sample quantity to container (in this example: 20 pcs)		
4. Changing the reference sample quantity:		REF 10pcs
5. Select reference sample quantity: In increments of 1 (1, 2, 3, etc. to 100) In increments of 10 (10, 20, etc. to 100)	Repeatedly:  Press briefly  press and hold	REF 20pcs
6. Confirm selected reference sample quantity and start application. The current reference weight remains stored until a new reference weight remains stored until a new reference is set or the power supply is interrupted		+ 20pcs nRef 20 pcs wRef 1.07 g
7. Add desired number of pieces.		+ 500pcs
8. If desired, print quantity		Qnt + 500 pcs
9. Toggle display between mean piece weight, weight, quantity	Repeatedly: 	+ 1.07 g Δ * + 535.0 g * + 500pcs *
10. Unload the balance		- 2 1pcs *
11. Repeat as needed, starting from Step 7		
12. Delete reference value		0.0 g

Weighing in Percent

Display symbol: %

Purpose

This application program allows you to obtain weight readouts in percent which are in proportion to a reference weight.

Preparation

- Select the Weighing in percent application in the menu: see "Configuration."
- Set the following parameters:

APPLIC. Application program

└ PERCENT Weighing in percent

- └ DEC.PLCS. Decimal places
 - └ NONE No decimal places
 - 1 DEC. PL. 1 decimal place
 - 2 DEC. PL. 2 decimal places
 - 3 DEC. PL. 3 decimal places

o = Factory setting

Changing the Reference Percentage

Activate function:

Press the  key

Select the desired reference (1 to 100):

In increments of 1: Press the  key briefly

In increments of 10: Press and hold the  key.

The percentage is stored in battery-backed memory.




Printout: Weighing in percent

pRef	100	: Reference percentage
Wxx%	111.6 g	: Reference weight net xx% for selected reference percentage
Prc	+ 94.9 %	: Calculated reference percentage

Example: Determining residual weight in percent

Parameter settings: APPLIC.: PERCENT (menu code 2. 4.)

Reference percentage: REF 100%

Step	Key (or instruction)	Display/Data output
1. Tare the balance	TARE	0.0 g
2. Place sample equal to 100% on the balance (in this example: 111.6 g)		
3. Information: Enter reference percentage (Changing the reference: see the previous page)	SELECT	REF 100 %
4. Initialize the balance The current reference weight remains stored until a new reference is set or the power supply is interrupted	ENTER	+ 100.0 % * pRef 100 % Wxx% 111.6 g
5. Remove sample (e.g., for drying)		
6. Place unknown weight on balance (in this example: 322.5 g)		+ 94.9 % *
7. If desired, print percentage	PRINT	Prc + 94.9 %
8. Toggle display between weight and percentage	Repeatedly: SELECT	+ 105.9 g * + 94.9 % *
9. Clear display of residual weight and reference percentage	CLEAR	+ 105.9 g
10. If desired, print net residual weight	PRINT	N + 105.9 g

Calculation

Display symbol: C

Purpose

With this application program you can calculate weight value using a multiplier or divisor. This can be used, for example, to determine the weight per unit area, or "gsm" weight (grams per square meter), of paper.

Preparation

- Select the Calculation application in the menu:
see "Configuration."
- Set the following parameters:

APPLIC. Application program

└ CALC. Calculation

└ METHOD Operator

└ o MUL. Multiplier
└ DIV. Divisor

DEC.PLCS. Decimal places

└ NONE No decimal places
└ o 1 DEC. PL. 1 decimal place
└ 2 DEC. PL. 2 decimal places
└ 3 DEC. PL. 3 decimal places

o = Factory setting

Setting the Factor or Divisor

Activate function:

Press the **SELECT** key

Select a number of up to 7 digits and, if needed, one decimal point (0.000001 to 9999999):

In increments of 1: Press the **SELECT** key briefly

To increase the value without pressing repeatedly:

Press and hold the **SELECT** key.

The selected operator is stored in battery-backed memory.

Printout: Calculation

Mul	+	1.2634	:	Multiplier
Div	+	0.6237	:	Divisor
Res	+	79.7	o :	Result

Example:

Calculating the weight per unit area of paper: An A4 sheet of paper is used in this example, with surface dimensions of $0.210 \text{ m} \times 0.297 \text{ m} = 0.06237 \text{ m}^2$. To determine the weight per unit area, the total weight is divided by the surface.

Parameter settings:

APPLIC.: CALC...: METHOD: DIV. (menu code 2. 8. 1. 2)

Step	Key (or instruction)	Display/Data output
1. Tare the balance		0.00 g
2. Activate divisor input		-----0.
3. Set the divisor (in this example: Position the decimal point, Enter numerals Repeatedly	 5 × 2 ×	0.06237:0000006000 or press and hold: , etc.
4. Store the divisor and initialize the balance The current divisor remains stored in battery-backed memory until the setting is changed		+ 0.0 ° Div 0.6237
5. Weight per unit area: Place an A4 sheet of paper on the balance	 	+ 79.7 ° *
6. If desired, print result		Res + 79.7 °
7. Toggle display between weight and calculated value	Repeatedly: 	+ 4.97 g *
8. Unload the balance	 	+ 0.0 ° *
9. Repeat as needed, starting from Step 5		

Animal Weighing/Averaging


Display symbol:

Purpose

Use this program to determine the weights of unstable samples (e.g., live animals) or to determine weights under unstable ambient conditions. With this program, the balance calculates the weight as the average of a defined number of individual weighing operations (also referred to as "sub-weighing operations").

Preparation

- Select the Animal weighing application in the menu: see "Configuration."
- Set the following parameters:

Changing the Number of Subweighing Operations

Activate function:

Press the key

Select the desired number of measurement (1 to 100):

In increments of 1:

Press the key briefly

In increments of 10:

Press and hold the key.




The selected number of measurements is stored in battery-backed memory.

Printout: Animal weighing

mDef	20	: Number of sub-weighing operations
x-Net +	410.1 g	: Calculated average

Example: Determining animal weight with automatic start and 20 subweighing operations (measurements)

Parameter settings: APPLIC.: ANIMALW. (menu code 2. 7.)

Step	Key (or instruction)	Display/Data output
1. Place animal weighing bowl on the balance		22.6 g
2. Tare the balance		0.0 g
3. Change the number of subweighing operations:		REF 30
4. Set number of measurements: In increments of 1 (1, 2, 3, etc. to 100)	Repeatedly:  Press briefly	
In increments of 10 (10, 20, etc. to 100)	 press and hold	REF 20
5. Confirm number of measurements and start automatic animal weighing The number of measurements remains stored in battery-backed memory until the setting is changed		+ 0.0 g *
6. Place first animal in bowl. The balance delays the start of measurements until the difference between 2 measurements meets the criterion		000 20 19 1
<input type="radio"/> If you selected the manual start mode, press		
7. Read off the result The result is displayed with the "*" symbol (= calculated value) and remains displayed until the sample (animal) is removed from the load plate (bowl)		+ 410.1 g Δ*
		mDef 20 x-Net + 410.1 g
8. Unload the balance		+ 0.0 g *
9. Weigh next animal (if desired) Next weighing series begins automatically		

Net-total Formulation


Display symbol: ↴

Purpose

With this application program you can weigh individual components either by their individual weight or by the total weight. You can print out both the total weight and the individual weights of the components.

Preparation

- Select the Net-total application in the menu:
see "Configuration."
- Set the following parameters:

o = Factory setting

Features

- Weigh up to 99 components from "0" to a defined total component weight.
- Store component weights (printout shows Comp xx), with
 - display zeroed automatically after value is stored, and
 - automatic printout
- Clear component memory following cancellation of the weighing sequence (by pressing **CLEAR**) and printout of the total weight.
- Toggling between component weight and total weight by pressing and holding **SELECT** (< 2 sec).
- Printout of the total of the individual component weights (T - Comp)




Printout: Net-total formulation

Comp 2+	278.1 g	: Second component
T-Comp+	2117.5 g	: Sum of components

Example: Counting parts into a container

Parameter settings:

APPLIC.: NET-TOT (menu code 2.5.)

Step	Key (or instruction)	Display/Data output
1. Place empty container on the balance.		65.0 g
2. Tare the balance		0.0 g
3. Add first component		+ 120.5 g
4. Store component data		+ 0.0 g * NET Comp 1+ 120.5 g
5. Add next component		+ 70.5 g * NET
6. Store component data		+ 0.0 g * NET Comp 2+ 70.5 g
7. Weigh in further components as desired	Repeat steps 5 and 6	
8. Fill to desired final value view the current total weight value:		+ 191.0 g *
9. Print total weight and clear the component memory		+ 2117.5 g T-Comp+ 2117.5 g

Totalizing

Display symbol: ↓


Purpose

With this application program you can add values from successive, mutually independent weight values to a total that exceeds the capacity of the balance.

Preparation

- Select the Totalizing application in the menu:
see "Configuration."
- Set the following parameters:

APPLIC. Application program

o = Factory setting

Features

- Totalizing memory for up to 99 values
- Store component weights (print-out shows Comp xx), with automatic printout
- Toggle display between the current individual weight value and the value in totalizing memory by pressing **SELECT**
- Printout of the total of the individual component weights (**S - Comp**)
- To close the application program and print the total weight: press **CLEAR**




Printout: Totalizing

Comp 2+	278.1 g	: Second component
S-Comp+	2117.5 g	: Totalizing memory

Example: Totalizing weight values

Parameter settings:

APPLIC.: TOTAL: COMP.PRT: ON (menu code 2. 6. 1. 2)

Step	Key (or instruction)	Display/Data output
1. Tare the balance	TARE	0.0 g
2. Place sample balance (in this example: 380 g)		+ 380.0 g
3. Store value in memory	ENTER	+ 380.0 g * Comp 1+ 380.0 g
4. Remove sample		+ 0.0 g *
5. Place the next sample on the balance (in this example, 575 g)		+ 575.0 g *
6. Store value in memory	ENTER	+ 955.0 g * + 575.0 g * Comp 2+ 575.0 g
7. View the value in totalizing memory	SELECT	+ 955.0 g Δ*
8. Weigh in further components as desired	Repeat steps 5 and 6	
9. Print total weight and clear the totalizing memory	CLEAR	0.0 g S-Comp+ 2117.5 g

Density Determination

Display symbol: ΔΔ

Purpose

This application program lets you determine the density of solid substances using the buoyancy method. You can have results displayed with one decimal place, or no decimal places: see "Configuration." Note: the sample holder and suspension wire used in the example below are not included with the balance.

Printout for Density Determination

Wa	+	20.0	g	: Weight in air
WfL	+	15.0	g	: Weight in liquid
Rho		4.0	o	: Result: density of the sample

Preparation

- Select the Density Determination application in the menu:
see "Configuration."
- Set the following parameters:



APPLIC. Application program
└ DENSITY Density determination
 └ DEC.PLCS. Decimal places
 └ NONE No decimal places
 o DEC.PL. 1 decimal place

o = Factory setting

Example: Determining the density of a solid sample.

Parameter settings:

APPLIC.: DENSITY: DEC.PLCS: DEC.PL. (menu code 2. 9. 1. 2)

Step	Key (or instruction)	Display/Data output
1. Attach sample holder to suspension wire		
2. Tare the balance	(TARE)	0.0 g
3. Start application program	(ENTER)	
4. Confirm "AIR" display	(ENTER)	AIR ?
5. Determine the weight of the sample in air: Place sample on the balance		+ 20.0 g *
6. Store value for weight in air	(ENTER)	
7. Remove sample from the balance		WATER ?
8. Determine weight in liquid: place sample in holder		
9. Confirm "WATER" display	(ENTER)	0.0 g *
10. Immerse sample in liquid		+ 15.0 g *
11. Store value for weight in liquid, view result, and print	(ENTER)	+ 4.0 o * Wa + 20.0 g WfL + 15.0 g Rho 4.0 o
12. Delete result	(CLEAR)	
13. Repeat as desired, starting from Step 3.		

Mass Unit Conversion


Purpose

With this application program you can change the weight value displayed from the basic weight unit to any of 4 application weight units (see table on next page).

Preparation

- Select the Unit application for toggling weight units:
see chapter on "Configuration" (Parameter Settings)
- Set the following parameters:

APPLICATION program


o = Factory setting

Printout for Wt. Unit Toggling

+ 100.0 g	: Weight with 16-character data output format
N + 0.22046 lb	: Weight with 22-character data output format
13-Jan-2005 08:35	: Data output format for two-line printout:
N + 3.5275 oz t	: Date/time and weight

Example: Change display from the basic unit (in this example, grams [g]) to pounds [lb] and then to Troy ounces [ozt].

Set the following parameters: *APPLIC.:UNIT* (code 2. 2.)

Step	Key (or instruction)	Display
Preparation:		
1. Begin selection of an application weight unit	(SELECT)	NONE ° [•]
2. Select an application unit; in this example, pounds (see table on next page)	Repeatedly: (SELECT)	POUNDS
3. Confirm the weight unit (pounds)	(ENTER)	POUNDS °
4. Select the next application weight unit; in this example: Troy ounces (see table on next page)	(ENTER), Repeatedly: (SELECT)	NONE ° [••] TROY OZ.
5. Confirm weight unit (Troy ounces)	(ENTER)	TROY OZ. °
6. Select other application units if desired (max. 4 total) (otherwise, confirm <i>NONE</i> by pressing (ENTER))		[•••]
7. Save selection	(CLEAR)	0.00 g
Conversion:		
8. Place sample on balance		+ 100.00 g
9. Toggle unit for weight value	Repeatedly: (ENTER)	+ 0.22046 lb + 3.5275 ozt

The following weight units are available in your Summit balance (in legal metrology, only units permitted by national law are available):

Menu item	Unit	Conversion factor	Display symbol
1) USERDEF. ¹⁾	Grams	1,00000000000	o
2) GRAMS (factory setting)	Grams	1.00000000000	g
3) KILOGR.	Kilograms	0.00100000000	kg
4) CARATS	Carats	5.00000000000	o
5) POUNDS	Pounds	0.00220462260	lb
6) OUNCES	Ounces	0.03527396200	oz
7) TROY OZ.	Troy ounces	0.03215074700	ozt
8) HKTAEL	Hong Kong taels	0.02671725000	tl
9) SNG.TAEL.	Singapore taels	0.02645544638	tl
10) TWN.TAEL	Taiwanese taels	0.026666666000	tl
11) GRAINS	Grains	15.4323583500	GN
12) PENNY.WT.	Pennyweights	0.64301493100	dwt
13) MILLIGR.	Milligrams	1000.00000000	mg
14) PT.P.LB.	Parts per pound	1.12876677120	o
15) CHN.TAEL	Chinese taels	0.02645547175	tl
16) MOMMES	mommes	0.26670000000	m
17) AUSTR.CT.	Austrian carats	5.00000000000	Kt
18) TOLA	Tola	0.08573333810	o
19) BAHT	Baht	0.06578947436	b
20) MESGHAL	Mesghal	0.21700000000	o
21) TONS	Tons	0.00000100000	t
22) LB / OZ ²⁾	Pounds : ounces	0.03527396200	lb oz
23) NEWTON	Newton	0.00980665000	N

¹⁾ = User-defined weight unit; can be loaded in the balance over an optional RS-232 or USB interface using a computer program.

²⁾ = The format for display of pounds/ounces cannot be changed: xx:yy.yy
x=lb, y=oz

 Some weight units may be blocked from use in legal metrology, depending on national verification laws.

ISO/GLP-compliant Printout/Record

Features

You can have device information, ID texts and date and time printed before (GLP header) and after (GLP footer) the values of a weighing series. These parameters include:

GLP header:

- Date
- Time at beginning of measurement
- Balance manufacturer
- Balance model
- Balance serial number
- Software version number
- Identification number of the current sampling operation

GLP footer:

- Date
- Time at end of measurement
- Field for operator signature

△ Operating the Balance with Printer:

Select the following settings on the balance and on the printer:

- Software handshake:
`SETUP: INTERF.: HANDSHK.: SOFTWARE.`
(menu code 1. 5. 4. 1)

Configuration

- Setting menu codes for the printout (see "Configuration"):
 - ISO/GLP-compliant printout or record only for calibration/adjustment:
`SETUP: PRNT.OUT: GLP: CAL.-ADJ.`
(menu code 1. 6. 7. 2) or
ISO/GLP-compliant printout or record always on:
`SETUP: PRNT.OUT: GLP: ALWAYS ON`
(menu code 1. 6. 7. 3)

- Line format for printout: include data ID codes (22 characters; factory setting):

`SETUP: PRNT.OUT: FORMAT: 22 CHAR.`
(menu code 1. 6. 6. 2)

- Formats for time:

`SETUP: PRNT.OUT: TIME: 24 H`
(menu code 1. 6. 8. 1) or

`SETUP: PRNT.OUT: TIME: 12 H`
(menu code 1. 6. 8. 2)

- Formats for date:

`SETUP: PRNT.OUT: DATE: DD.MMM.YY`
(menu code 1. 6. 9. 1) or

`SETUP: PRNT.OUT: DATE: MMM.DD.YY`
(menu code 1. 6. 9. 2)

△ No ISO/GLP-compliant record is output if any of the following settings are configured:

`SETUP: PRNT.OUT: PRINT: AUT.W/O` or
`AUT.WITH` (menu code 1. 6. 1. 3, 1. 6.
1. 4,) or `FORMAT: 16 CHAR.` (menu
code 1. 6. 6. 1)

Function Keys

Transfer header and first measured value: press **PRINT**

> The header is included with the first printout/data record.

To output header and reference data automatically when an application program is active: press **ENTER**

Exit the application:

- 1) To send the GLP footer:
press **CLEAR**
- 2) Quit application program:
press **CLEAR** again

The ISO/GLP-compliant printout can contain the following lines:

-----	Dotted line
17-Aug-2005 10:15	Date/time (beginning of measurement)
DENVER	Balance manufacturer
Mod. SI-8001	Model
Ser. no. 10105355	Balance serial number
Ver. no. 00-36-01	Software version
ID 2690 923	ID
-----	Dotted line
L ID	Measurement series no.
nRef 10 pcs	Counting: reference sample quantity
wRef 21.14 g	Counting: reference weight
Qnt + 567 pcs	Counting result
-----	Dotted line
17-Aug-2005 10:20	Date/time (end of measurement)
Name:	Field for operator signature
-----	Blank line
-----	Dotted line

ISO/GLP-compliant printout for external calibration/adjustment:

-----	Dotted line
17-Aug-2005 10:30	Date/time (beginning of measurement)
DENVER	Balance manufacturer
Mod. SI-8001	Model
Ser. no. 10105352	Balance serial number
Ver. no. 00-36-01	Software version
ID 2690 923	ID
-----	Dotted line
Cal. Ext. Test	Calibration/adjustment mode
Set + 5000.0 g	Calibration weight
Diff. + 0.2 g	Difference determined in calibration
Cal. Ext. Complete	Confirmation of completed calibration procedure
Diff. 0.0 g	Difference from target following adjustment
-----	Dotted line
17-Aug-2005 10:32	Date/time (end of measurement)
Name:	Field for operator signature
-----	Blank line
-----	Dotted line

Data Interface

Purpose

Your balance comes equipped with an interface port for connection to a computer or other peripheral device. You can use a computer to change, start and/or monitor the functions of the balance and the application programs.

Preparation

You can set these parameters for other devices in the Setup menu (see the chapter entitled "Configuring the Balance").

You will also find a detailed description of the available data interface commands in the file "Data Interface Descriptions for S and SI Models", which you can download from the Denver website (www.denverinstrument.com).

The many and versatile properties of these balances can be fully utilized for printing out records of the results when you connect your balance to a Denver data printer. The recording capability for printouts makes it easy for you to work in compliance with ISO/GLP.

Hardware Handshake

With a 4-conductor interface, 1 or 2 characters can be transmitted after CTS.

These connections must be made when the balance is connected through the RS-232C port.

Pin Assignment Chart for Data Interface

Summit Series

Balance 9-pin port

RxD	2
TxD	3
DTR	4
Signal ground	5
CTS	8

Standard RS-232

9-pin connector

3	TxD
2	RxD
4	DTR
5	Signal ground
8	CTS

Summit Series

Balance 9-pin port

RxD	2
TxD	3
DTR	4
Signal ground	5
CTS	8

Standard RS-232

25-pin connector

2	TxD
3	RxD
20	DTR
7	Signal ground
5	CTS

Troubleshooting Guide

Error codes are shown on the main display for approx. 2 seconds. The program then returns automatically to the previous mode.

Display	Cause	Solution
No segments appear on the display	No AC power is available	Check the AC power supply
	The power supply is not plugged in	Plug in the power supply
HIGH	The load exceeds the balance capacity	Unload the balance
LOW or ERR 54	Something is touching the weighing pan	Move the object that is touching the weighing pan
ERR 54, typical	Weighing system defect	Contact Denver Service Center
APP.ERR.	Cannot store data: Load on weighing pan too light or no sample on pan while application is active	Increase load
BIS.ERR.	Data output not compatible with output format	Change the configuration in the operating menu
PRT.ERR.	Interface port for printer output is blocked	Reset the menu factory settings, or Contact your local Denver Service Center
ERR 02	Calibration parameter not met; e.g.: – balance not tared – load on weighing pan	Calibrate only when zero is displayed – Press (TARE) to tare the balance – Unload the balance
ERR 10	The (TARE) key is blocked when there is data in the second tare memory (net-total); only 1 tare function can be used at a time	Press (CLEAR) to clear the tare memory and release the tare key
ERR 11	Tare memory not allowed	Press (TARE)
The weight readout changes constantly	Unstable ambient conditions (excessive vibration or draft) at the place of installation A foreign object is caught between weighing pan and balance housing	Set up the balance in another area Remove the foreign object
The weight readout is obviously wrong	The balance was not calibrated/adjusted Balance not tared before weighing	Calibrate/adjust the balance Tare or zero the balance before weighing

If any other errors occur, contact your local Denver Service Center.

Contact information: <http://www.denverinstrument.com>

Care and Maintenance

Service


On request, Denver Instrument can provide you with individual service information.

Repairs

Repair work must be performed by trained service technicians. Any attempt by untrained persons to perform repairs may result in considerable hazards for the user.

Cleaning

- Unplug the AC adapter from the wall outlet (mains supply). If you have an interface cable connected to the balance port, unplug it from the port.
- ⚠ Do not used solvents or aggressive cleaning agents on the following parts: keypad overlay, power jack, data interface
- After cleaning, wipe down the balance with a soft, dry cloth.

On analytical balances remove and clean the weighing pan as follows:

- Reach beneath the drip/breeze and lift it carefully, together with the pan support, to avoid damaging the weighing system.

- ⚠ Make sure that no liquid enters the balance housing.

Cleaning Stainless Steel Surfaces

Clean all stainless steel parts regularly. Remove the stainless steel weighing pan and thoroughly clean it separately. Use a damp cloth or sponge to clean stainless steel parts on the balance. You can use any household cleaning agent that is suitable for use on stainless steel. Clean stainless steel surfaces only by wiping them down. Then rinse the equipment thoroughly, making sure to remove all residues. Afterwards, allow the equipment to dry. If desired, you can apply oil to the cleaned surfaces as additional protection.

Recycling

Safety Inspection

If there is any indication that safe operation of the balance is no longer warranted:

- Turn off the power and disconnect the equipment from AC power immediately.

- > Lock the equipment in a secure place to ensure that it cannot be used for the time being.

Notify your nearest Denver Service Center. Repair work must be performed by trained service technicians.

We recommend having the power supply inspected by a certified electrician at regular intervals, according to the following checklist:

- Insulating resistance: > 7 megaohms measured with a constant voltage of at least 500 volts at a 500 K-ohm load
- Leakage current: < 0.05 mA measured with a properly calibrated multimeter

Information and Instructions on Disposal and Repairs

Packaging that is no longer required must be disposed of at the local waste disposal facility. The packaging is made of environmentally friendly materials that can be used as secondary raw materials.

The equipment, including accessories and batteries, does not belong in your regular household waste. The EU legislation requires its Member States to collect electrical and electronic equipment and dispose of it separately from other unsorted municipal waste with the aim of recycling it.

In Germany and many other countries, Denver Instrument takes care of the return and legally compliant disposal of its electrical and electronic equipment on its own. These products may not be placed with the household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other Member States of the European Economic Area (EEA), please contact our service technicians on location or our Service Center in Goettingen, Germany:

Denver Instrument GmbH
Robert-Bosch-Breite 10
37079 Goettingen, Germany

In countries that are not members of the European Economic Area (EEA) or where no Denver affiliates, subsidiaries, dealers or distributors are located, please contact your local authorities or a commercial disposal operator.

Prior to disposal and/or scrapping of the equipment, any batteries should be removed and disposed of in local collection boxes.

Denver, its affiliates, subsidiaries, dealers and distributors will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the accompanying leaflet/manual or visit our Internet website (www.denverinstrument.com) for comprehensive information that includes our service addresses to contact if you plan to send your equipment in for repairs or proper disposal.

Overview

Specifications

Specifications

Built-in motorized calibration weight	All SI models	
AC power source/power requirements, voltage, frequency	AC adapter 230 V or 115 25 V, +15% to -20%, 48–60 Hz	
Power consumption	VA	maximum 16; typical 8 (STNG6)
Approx. hours of operation with the YRB05Z rechargeable battery pack (backlighting on)	h	35

Ambient Conditions

The specifications given here are ensured under the following ambient conditions:

Operating temperature range	+10 to +30°C (273 to 303 K, 50 to 86°F)
-----------------------------	---

Specifications for Individual Models

Model		SI-234 IS-234	SI-114 IS-114	SI-64 IS-64
Weighing capacity	g	230	110	60
Readability	g	0.0001	0.0001	0.0001
Tare range (subtractive)	g	230	110	60
Repeatability (std. deviation)	≤± g	0.0001	0.0001	0.0001
Linearity	≤± g	0.0002	0.0002	0.0002
Response time (average)	s	2.5		
Sensitivity drift within +10 to +30°C	≤±/K	$2 \cdot 10^{-6}$		
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels; display update: 0.1–0.4 (depends on filter level selected)		
External calibration weight	g	200 (E2)	100 (E2)	50 (E2)
Net weight, approx.:		4.8 4.4 kg, 10.6 9.7 lb		
Weighing pan size		90 Ø mm, 3.54 Ø in.		
Weighing chamber height		230 mm, 9.1 in.		
Dimensions (W × D × H)		230 × 303 × 330 mm, 9.1 × 11.9 × 13.0 in.		

Model		SI-603 IS-603	SI-403 IS-403	SI-203 IS-203
Weighing capacity	g	600	400	200
Readability	g	0.001	0.001	0.001
Tare range (subtractive)	g	600	400	200
Repeatability (std. deviation)	≤± g	0.001	0.001	0.001
Linearity	≤± g	0.002	0.002	0.002
Response time (average)	s	1		
Sensitivity drift within +10 to +30°C	≤±/K	$2 \cdot 10^{-6}$		
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels; display update: 0.05–0.4 (depends on filter level selected)		
External calibration weight (of at least accuracy class...)	g	500 (E2)	200 (E2)	200 (F1)
Net weight, approx.:		3.6 3.2 kg, 7.9 7.0 lb		
Weighing pan size		115 Ø mm, 4.52 Ø in.		
Dimensions (W × D × H)		230 × 303 × 136 mm, 9.1 × 11.9 × 5.4 in.		

Specifications for Individual Models

Model		SI-6002 I S-6002	SI-4002 I S-4002	SI-2002 I S-2002	SI-602 I S-602	SI-402 I S-402
Weighing capacity	g	6000	4000	2000	600	400
Readability	g	0.01	0.01	0.01	0.01	0.01
Tare range (subtractive)	g	6000	4000	2000	600	400
Repeatability (std. deviation)	<± g	0.01	0.01	0.01	0.01	0.01
Linearity	<± g	0.02	0.02	0.02	0.02	0.02
Response time (average)	s	1.1	1.1	1.1	1.1	1.0
Sensitivity drift within +10 to +30°C	≤±/K	$2 \cdot 10^{-6}$	$2 \cdot 10^{-6}$	$2 \cdot 10^{-6}$	$5 \cdot 10^{-6}$	$5 \cdot 10^{-6}$
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels; display update: 0.1–0.4 (depends on filter level selected)				
External calibration weight (of at least accuracy class...)	g	5000 (E2)	2000 (E2)	2000 (F2)	500 (F2)	200 (F2)
Net weight, approx.:		3.5 3.1 kg, 7.7 6.8 lb		2.6 2 kg, 5.7 4.4 lb		
Weighing pan size		180 × 180 mm, 7.1 × 7.1 in.			150 Ø mm, 5.91 Ø in.	
Dimensions (W × D × H)		230 × 303 × 91 mm, 9.1 × 11.9 × 3.58 in.			230 × 303 × 87 mm, 9.1 × 11.9 × 3.4 in.	
Model		SI-8001 I S-8001	SI-4001 I S-4001			
Weighing capacity	g	8000	4000			
Readability	g	0.1	0.1			
Tare range (subtractive)	g	8000	4000			
Repeatability (std. deviation)	≤± g	0.1	0.1			
Linearity	≤± g	0.1	0.1			
Response time (average)	s	1	1			
Sensitivity drift within +10 to +30°C	≤±/K	$10 \cdot 10^{-6}$				
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels; display update: 0.05–0.4 (depends on filter level selected)				
External calibration weight (of at least accuracy class...)	g	5000 (F2)	2000 (F2)			
Net weight, approx.:		3.5 2.7 kg, 7.7 6.0 lb				
Weighing pan size		180 × 180 mm, 7.1 × 7.1 in.				
Dimensions (W × D × H)		230 × 303 × 91 mm, 9.1 × 11.9 × 3.6 in.				

Accessories

Draft ring for round pan units

In-use cover

Printer with connecting cable

Cable with 9-pin connector for balance

Weigh below hanger*

Calibration weights: contact your distributor for a complete list

* = not allowed in applications involving legal metrology
(such as legal for trade)

CE Marking

The balance complies with the following EC Directives and European Standards:

**Council Directive 89/336/EEC
“Electromagnetic compatibility
(EMC)”**

Applicable European Standards:

Limitation of emissions:

In accordance with product standard
EN 61326-1 Class B (residential area)

Defined immunity to interference:
in accordance with product standard
EN 61326-1
(minimum test requirements,
non-continuous operation)

Important Note:

The operator shall be responsible for
any modifications to Denver equip-
ment and must check and, if neces-
sary, correct these modifications.

On request, Denver will provide infor-
mation on the minimum operating
specifications (in accordance with the
Standards listed above for defined
immunity to interference).

**73/23/EEC “Electrical equipment
designed for use within certain
voltage limits”**

Applicable European Standards:

EN 60950

Safety of information technology
equipment including electrical busi-
ness equipment

EN 61010

Safety requirements for electrical
equipment for measurement, control
and laboratory use Part 1: General
requirements

If you use electrical equipment
in installations and under ambient
conditions requiring higher safety
standards, you must comply with the
provisions as specified in the applica-
ble regulations for installation in your
country.

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Copyright by Denver Instrument GmbH, Goettingen, Germany.
All rights reserved. No part of this publication may be reprinted
or translated in any form or by any means without the prior
written permission of Denver Instrument GmbH.

The status of the information, specifications and illustrations in
this manual is indicated by the date given below. Denver Instrument
GmbH reserves the right to make changes to the technology, fea-
tures, specifications, and design of the equipment without notice.
Status: April 2006, Denver Instrument GmbH, Goettingen,
Germany.

Europe, Asia and Australia:

Denver Instrument GmbH
Robert-Bosch-Breite 10
37079 Goettingen, Germany
Tel: +49-551-20977-30
Fax: +49-551-20977-39

Internet:

www.denverinstrument.com

U.K. and Ireland:

Denver Instrument Company
Denver House, Sovereign Way
Trafalgar Business Park
Downham Market
Norfolk PE38 9SW England
Tel: +44-136-63862-42
Fax: +44-136-63862-04

North and South America:

Denver Instrument Company
1401 17th Street, Suite 750
Denver, Colorado 80202
1-800-321-1135
Tel: +1-303-431-7255
Fax: +1-303-423-4831

98648-013-93