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
Abb./Fig. 1

1.1 Range of Application for the Eppendorf Comforpette 4700

- Piston stroke pipette with ejection mechanism and fixed volume
- Pipette made from resistant materials (e. g., polypropylene, dyflor, ceramics, stainless steel)
- Pipetting of liquids in the microliter range
- Polypropylene pipette tips
- Single volume pipettes: 1 – 5 – 10 – 20 – 25 – 30 – 40 – 50 – 100 – 200 – 250 – 300 – 400 – 500 – 1,000 µl
- Tri-volume pipettes, adjusted by turning the control button: 10 – 20 – 25 µl
20 – 25 – 50 µl
50 – 75 – 100 µl
100 – 200 – 250 µl
200 – 300 – 500 µl
500 – 750 – 1,000 µl
- Pipette Stand 4707 for pipettes and tips (Fig. 2).

Pipette with attached tip to functional unit.

IMPORTANT

Only use tips which are designed to fit the nose cone (Fig. 1), otherwise:
– a perfect seal
– tip ejection
are not guaranteed.
E. g. Eppendorf Comfortips, Standartips and Comfortips cristal.

With the Eppendorf Comforpette 4700, you have a pipette which provides reproducible results. It fits comfortably into your hand. Mess or damage caused by chemicals is reduced as a result of the cone shape.

The piston stroke is divided into three sections: measuring, blow-out and ejection. The piston stroke is subdivided by two intermediate stops:

- First stop: End of measuring displacement, beginning of the blow-out
- Second stop: End of blow-out, liquid completely dispensed
- Final stop: Ejection of pipette tip

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Abb./Fig. 4

2. Pipette tips

- Ready-to-use
- Storage in Pipette Stand 4707 (Fig. 2)
- Pipette stand capable of being expanded by simple attachment (Fig. 2)
- Pipette tips (1–10 μ l) cristal, silver-colored pushbutton (10–100 μ l) yellow, yellow pushbutton (200–1,000 μ l) blue, blue pushbutton
- Available as Standartips in bags or as Comfortips in tip storage boxes
- Under sterile conditions Eppendorf Steriltips
- Standartips – tip storage box with adjustable partition (Fig. 4)
- Comfortips – without cover in Pipette Stand 4707
- Eppendorf Comfortips cristal (Fig. 3), autoclavable up to 120°C
- Tips easily taken up by the pipette, attachment by hand no longer necessary
- When not in use, can be protected by closing lid.

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The Eppendorf Comforpette 4700 is simple to use. Nevertheless, the following instructions on how to operate should be observed:

- Carry out a few pipettings **without** liquid
- Press button up to the first stop: measuring volume is defined
- Second stop: dispensing of the remaining liquid
- Third stop: tip ejection
- Tips must always fit firmly onto nose cone.

Pipettes 1–10 µl:
silver-colored pushbutton,
cristal tips

Pipettes 10–100 µl:
yellow pushbutton,
yellow tips

Pipettes 200–1,000 µl:
blue pushbutton,
blue tips.

3.6 Volume Selection (Fig. 5)

The volume which is displayed opposite the red mark is pipetted by the tri-volume pipette.

Volume Selection:

- Press the button down to the mark in the volume ring (e. g. 500 : 750) and turn to the right or left
- Must lock!
- Mark must be visible again
- Pipette ready for operation.

3.7 Filling (Fig. 5)

IMPORTANT: Never pipette liquid **without tip**.

- Keep pipette almost vertical
- Press button up to first stop
- Immerse tips 2–3 mm into the liquid **Slowly** allow the button to glide back (never let it snap back)
- Slide tip out along the inside of the vessel
- Wipe off any droplets with lint-free tissue.

IMPORTANT: Never put down pipettes with filled tips!

When pipetting viscous liquids or liquids with a low surface tension (e. g. serum, organic solvents):

- Rinse tip two or three times with the liquid to be pipetted
- Wipe off tips with lint-free tissue.

To rinse the button, press and release to the first stop.

There are two possible ways emptying the tip:

1. Vessel
2. Chromatographic paper
- Hold tip against the inside the vessel or directly to the surface
- Press button down to the first stop and wait 1–3 sec.
- Press button on to the second stop
- Slide pipette out along the inside of the vessel.

When pipetting small volume into a larger amount of liquid, remaining sample should be rinsed out (see 3.3). Then blow-out to empty remaining liquid.

Note:

When pipetting different liquids, change the tips at the start of each run. Rinse them well and then dab tissue.

When pipetting serum, rinse tip each time. If small bubbles collect in the tip, the tip **must** be changed.

Differences in temperature of $\pm 5^\circ\text{C}$ from the room temperature have an effect on the precision of the pipetting volume.

Rinse tips several times with liquid.

Variations in accuracy must be tolerated.

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4. Calibration

4.1 Testing

For checking the pipetting volume, there are two controls of interest:

- Determining the precision (reproducibility)
- Checking the accuracy (absolute volume)

To determine the absolute error of the pipette, 5 measurements are usually necessary. If the precision is also to be tested, 10, preferably 20 pipettings, at least, are necessary in order to determine the standard deviation and coefficient variation of the series.

4.2 Photometric Testing

Production of exact dilutions with officially calibrated, degreased glass pipettes or burettes. Phenol red (phenolsulfonphthalein) has proved itself well suited for this purpose.

Phenol red, e. g., Merck
Typ No.: 7241 solid matter
Typ No.: P 1735 ampules

Original solution:

defined amount phenol red
+ 0.1 % soda solution
(1 g NaOH/l)
Hg 546 nm $A=1.0-1.5$
measured against 0.1 % NaOH

Diluting solution:

(Rinse cuvettes)
Hg 546 nm $A=0.1$
measured against 0.1 % NaOH

Gravimetric Test

A micro- or semi-micro balance with the appropriate sensitivity is required for this test.

After determining the corresponding empty weights, there are two possibilities:

- Weighing of a pipette tip filled with water (only for pipettes 1–10 μ l and 10–100 μ l)
- Weighing the water dispensed in sealable micro test tubes (pipettes 200–1,000 μ l).

For routine measurements, milligrams can be equated with microliters; for exact determinations, the density of the water must be taken into consideration.

Fig. 6 + 7

- Fill bottle 3/4 full with water or colored solution
- Slide the tubing into the neck of the bottle and seal tightly
- Insert red adapter for the pipette into slot
- Before testing, set mark at level of the base
- Close air vent with index finger.
- Press bottle with index finger and thumb of the other hand until the distance of the liquid level is 2–4 cm below the tubing end
- Open air vent
- Press pipette 10–1,000 μ l with tip into adapter and remove together from the slot
- To test the pipette with crystal tip, remove the adapter and press tip directly into tubing
- Pull upwards
- Push the marker approximately 5 mm below the meniscus
- Press pushbutton approx. 20 times.

Each time the button reaches the upper position, the meniscus should be in its original position. If it is constantly lower, the pipette or the tip seal is leaking.

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5.1 Physical Limitations

- Temperature of the liquid $\pm 5^{\circ}\text{C}$ above or below the room temperature
- Density greater than 2.0 g/cm^3 (20°C)
- Liquids with high vapour pressure (1.293 g/cm^3), e. g., n-butane
- organic liquids with temperatures greater than $+60^{\circ}\text{C}$.

5.2 Chemical Limitations


If used properly, damage caused by chemicals is not possible.

Improper Use:

Liquid in the pipette can cause crystal formation. The cone is then blocked. Suction is affected. Organic solvents can cause swelling of the sealing rings. Piston moves jerkily.

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If the pipette is properly used there are no dangers whatsoever. To avoid infection, when working with infectious substances the pipette tip should not be touched, this is possible by actuating the ejection mechanism.

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Droplets on the inside of the tip	Wettability of the plastic	- Use a new tip
Pipette drips, volume too small	Contaminated piston Tip is loose Piston loosened Wrong tip Worn seal Pipette nose cone damaged	- Clean piston (9.2) - Press tip on tightly - Tighten piston - Use Original Eppendorf - Exchange lower part (pipettes 1–10 μ l, 10–100 μ l) and 200–1,000 μ l) - Exchange nose cone (pipettes 10–1,000 μ l) Exchange lower part (pipettes 1–10 μ l)
Jerky return of the Pipette-pushbutton	Too much resistance on the piston from the seal (too much friction) Swelling of the sealing ring caused by improper pipetting of organic solvents or caustic liquids	- Clean piston (9.2) - Exchange lower part (pipettes 1–10 μ l, 10–100 μ l and 200–1,000 μ l)
Pipette does not pick up	Nose cone clogged (Pipettes 10–100 μ l) Piston bent (pipettes 1–10 μ l)	- Remove and clean cone - Exchange nose cone (9.2) - Return pipette

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8. Short instructions

These brief instructions may only be used if the operator is thoroughly acquainted with the detailed instructions, otherwise errors may not be avoidable.

- Take up the pipette from tip storage box
- Press button to first stop
- Immerse tip a few millimeters into the solution
- Allow button to slide back (never let it snap back)
- When removing the filled tip, slide along inside of vessel, possibly wiping off the tip with lint-free tissue
- When pipetting, hold tip against the inside of the vessel or (in the case of small volumes) to the surface of the liquid
- Press button slowly to the first stop and wait 1-3 seconds
- Press button to the second stop, pulling the tip from this position along the inside of the vessel or (in the case of small volumes) pulling from the liquid upper limit in the vessel
- To eject the tip, press the button to the final stop
- Keep pipette in the pipette stand.

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
Abb./Fig. 9

The pipette seals do not require maintenance. The piston must be greased.

- Clean externally with Isopropanol.
- Nose cone clogged:
 - Unscrew (Fig. 9)
 - Wipe piston with damp cloth
 - Connect to water jet clean with distilled water
 - Slightly grease O-Ring
 - Screw on.

The individual components shown in Fig. 10 are available as replacement parts.

Caution: After installation of spare parts or similar operations, check that the pipette works properly.

If the problem cannot be solved with the aid of the suggestions mentioned here, we suggest the pipette be returned to the Eppendorf distributor for repair or replacement.

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Liquid: Water
Temperature: approx. 20°C
No. of measurements: 10

Measurement: in accordance
with DIN 12 650

When non-aqueous or
moistening liquids (serum) are
used, these values may double

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1 Pipette with attached tip
1 Instruction Manual

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Technical specification subject
to change!

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Netheler + Hinz GmbH