
**DB-2D, DB-2DH, DB-3D
DB-3DL, DB-4D, DB-2TC
Dri-Block®**

OPERATOR'S MANUAL

DB-2D, DB-2DH, DB-3D, DB-3DL, DB-4D and DB-2TC**OPERATOR'S MANUAL**

CONTENTS	PAGE
SAFETY and INSTALLATION	2
English	2
Français	5
Deutsch	8
Español	11
Italia	14
THE DRI-BLOCK®	17
SPECIFICATION	18
DB-2D and DB-2DH	18
DB-3D and DB-3DL	19
DB-4D	20
DB-TC	21
Working environment	22
Unpacking	22
OPERATION	23
Installation	23
The overtemperature cutout	24
Power up screen	24
The front panel controls	25
Setting the operating temperature	27
After use	28
ADDITIONAL INFORMATION	29
Accessories	29
Operator maintenance	29
Calibration	30
Insert blocks and replacement parts	32

SAFETY AND INSTALLATION

Please read all the information in this booklet before using the unit.

WARNING

HIGH TEMPERATURES ARE DANGEROUS: they can cause serious burns to operators and ignite combustible material.

Techne have taken great care in the design of these units to protect operators from hazards, but operators should pay attention to the following points:

- USE CARE AND WEAR PROTECTIVE GLOVES TO PROTECT HANDS.
- DO NOT put hot objects on or near combustible objects.
- DO NOT operate the unit close to inflammable liquids or gases.
- DO NOT place any liquid directly in your unit.

OPERATOR SAFETY

All operators of Techne equipment must have available the relevant literature needed to ensure their safety. It is important that only suitably trained personnel operate this equipment, in accordance with the instructions contained in this manual and with general safety standards and procedures. If the equipment is used in a manner not specified by Techne the protection provided by the equipment to the operator may be impaired.

All Techne units have been designed to conform to international safety requirements and are fitted with an overtemperature cutout. On some models, the cutout is adjustable and should be set to suit the application. On all other models the cutout is preset to protect the unit.

If a safety problem should be encountered, switch off at the mains socket and remove the plug from the supply.

INSTALLATION

1. All Techne units are supplied with a power cable. This may be integral or plug-in.
2. Before connecting the mains supply, check the voltage against the rating plate. The rating plate is on the rear of the unit. Connect the mains cable to a suitable plug according to the table below.

Note that the unit must be earthed to ensure proper electrical safety.

Live	Brown	Black
Neutral	Blue	White
Earth	Green/yellow	Green

The fused plug supplied with the mains lead for use in the UK is fitted with the following value fuse to protect the cable: 5 amp.

The fuse in the unit protects the unit and the operator.

Note that units marked 230V on the rating plate work at 220V; units marked 120V work at 110V. In both cases, however, the heating rate will degrade by approximately 8%.

3. Plug the mains cable into the socket on the rear of the unit.
4. Place the unit on a suitable bench or flat workspace, or in a fume cupboard if required, ensuring that the air inlet vents on the underside are free from obstruction.
5. Note that the following symbols may be next to the indicator lamps on the front panel of the units and have the following meanings:
 - ~ the power indicator
 - ~~ the heater indicator
 - ! the overtemperature indicator
6. Symbols on or near the power switch of the unit have the following meanings:
 - I mains switch On
 - O mains switch Off

AFTER USE

When you have finished heating samples, remember that parts of the unit – the tubes, blocks and associated accessories – may be very hot. Take the precautions listed earlier.

GUARANTEE

The unit is guaranteed against any defect in material or workmanship for the period specified on the enclosed guarantee card. This period is from the date of purchase, and within this period all defective parts will be replaced free of charge provided that the defect is not the result of misuse, accident or negligence. Servicing under this guarantee should be obtained from the supplier.

Notwithstanding the description and specification(s) of the units contained in the operator's manual, Techne hereby reserves the right to make such changes as it sees fit to the units or to any component of the units.

This manual has been prepared solely for the convenience of Techne customers and nothing in this instruction book shall be taken as a warranty, condition or representation concerning the description, merchantability, fitness for purpose or otherwise of the units or components.

OPERATOR MAINTENANCE

NOTE: THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL.
REMOVING THE SIDE, FRONT OR REAR PANELS EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES.
THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.

In the unlikely event that you experience any problems with your unit which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault

observed and remember to return the unit in its original packing. Techne accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier. See the decontamination certificate supplied with your unit.

1. Cleaning

Before cleaning your unit ALWAYS disconnect it from the power supply and allow it to cool below 50°C.

Your unit can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

2. Overtemperature cutout

In the event of no heater power, check the mains plug and lead. Repeated operation of the cutout indicates a serious fault: you may need to return the unit to your supplier for repair.

3. Fuses

The unit is protected by one or two fuses. These should only be changed by suitably qualified personnel. If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

CONTACT INFORMATION

For technical, sales or servicing information, contact your local Techne dealer or

UK and Rest of World

Bibby Scientific Ltd
Beacon Road, Stone,
Staffordshire, ST15 0SA, UK
Tel: +44 (0)1785 812121
Fax: +44 (0)1785 813748
e-mail: techne@bibby-scientific.com
www.bibby-scientific.com

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THE DRI-BLOCK®

Please read all the information in this manual before using the unit.

Techne's Dri-Block® heaters provide a safe, dry, constant temperature source in the laboratory. The unit can be used for incubation, boiling, inactivation, wet ashing, sample concentration, enzyme analysis and many other clinical and industrial purposes.

The units cover the temperature range from 25°C to 200°C using machined aluminium alloy blocks as the heat transfer medium: temperature control circuits are built into the unit. Independent indicator lights indicate mains power, heater action and overtemperature cutout. An ON/OFF switch is mounted on the rear panel, together with the mains connector and a double fuse holder.

The units are constructed of polyphenylene sulphide plastic, which is strong yet lightweight and can resist sustained temperatures as high as 260°C. Even though the units heat up rapidly, highly efficient insulation ensures that cases remains cool enough to handle even at maximum operating temperatures.


The operating instructions in this Operator's Manual covers six of a range of Dri-Blocks®:

- DB-2D • DB-2DH • DB-3D • DB-3DL • DB-4D • DB-2TC

The features of individual models are outlined in the 'Specification' section of this manual.

Units not covered by this manual:

- DB-2A • DB-3 • DB-3A


Typical rear connections

SPECIFICATION

THE DB-2D AND DB-2DH DRI-BLOCKS®

Each unit has a recessed chamber which can hold one or two interchangeable insert blocks which accept test tubes or other sample containers. The temperature is set by 3 front panel buttons under the digital display.

The DB-2D and DB-2DH are supplied with an extraction tool for removing the blocks from the unit and two plastic half blocks for reducing heat losses when only one aluminium block is used.

Specification

	DB-2D	DB-2DH
Working temperature range	Ambient 0°C to 100°C	Ambient 0°C to 200°C
Temperature range	0.0°C to 105.0°C	0.0°C to 200.0°C
Temperature stability	±0.05°C at 40°C	±0.1°C at 40°C
Temperature stability	±0.1°C at 100°C	±0.15°C at 100°C
Temperature setting	3 push buttons	3 push buttons
Temperature display	4 digit LED	4 digit LED
Set point resolution	0.1°C	0.1°C
Set point accuracy	±1°C	±1°C
Max temperature variation between blocks	0.2°C at 40°C	0.2°C at 40°C
Heat up time 30°C to 37°C	8 minutes	11 minutes
Heat up time 30°C to 56°C	9 minutes	15 minutes
Heat up time 30°C to max	12 minutes	25 minutes
Number of blocks	2	2
Overall size L x W x H (mm)	202 x 260 x 105	202 x 260 x 105


Electrical supply

Voltage	Frequency	Power
230V	50Hz-60Hz	300W
120V	50Hz-60Hz	300W

THE DB-3D AND DB-3DL DRI-BLOCKS®

Each unit has a recessed chamber which can hold up to three inter-changeable insert blocks which accept test tubes or other sample containers. The temperature is set by 3 front panel buttons with digital display.

The DB-3D and DB-3DL are supplied with an extraction tool for removing the blocks from the unit.

**Specification**

	DB-3D	DB-3DL
Working temperature range	Ambient 0°C to 200°C	Ambient 0°C to 100°C
Temperature range	0.0°C to 200.0°C	0.0°C to 105.0°C
Temperature stability	±0.1°C at 40°C	±0.1°C at 40°C
Temperature stability	±0.15°C at 100°C	±0.15°C at 100°C
Temperature setting	3 push buttons	3 push buttons
Temperature display	4 digit LED	4 digit LED
Set point resolution	0.1°C	0.1°C
Set point to accuracy	±1°C	±1°C
Max temperature variation between blocks	0.2°C at 40°C	0.2°C at 40°C
Heat up time 30°C to 37°C	11 minutes	8 minutes
Heat up time 30°C to 56°C	15 minutes	9 minutes
Heat up time 30°C to max	25 minutes	12 minutes
Number of blocks	3	3
Overall size L x W x H (mm)	279 x 260 x 105	279 x 260 x 105


Electrical supply

Voltage	Frequency	Power
230V	50Hz-60Hz	450W
110V-120V	50Hz-60Hz	450W

THE DB-4D DRI-BLOCK®

This unit has a recessed chamber which can hold four interchangeable insert blocks which accept test tubes or other sample containers. The temperature is set by 3 front panel buttons under the digital display.

The DB-4D is supplied with an extraction tool for removing the blocks from the unit and two plastic half blocks for reducing heat losses when only one aluminium block is used.

**Specification****DB-4D**

Working temperature range	Ambient 0°C to 100°C
Temperature range	0.0°C to 100.0°C
Temperature stability	±0.1°C at 40°C
Temperature stability	±0.1°C at 100°C
Temperature setting	Push button
Temperature display	4 digit LED
Set point resolution	0.1°C
Set point to accuracy	±1°C
Max temperature variation between blocks	0.2°C
Heat up time 30°C to 37°C	13 minutes
Heat up time 30°C to 56°C	11 minutes
Heat up time 30°C to 100°C	15 minutes
Number of blocks	4
Overall size L x W x H (mm)	356 x 260 x 105


Electrical supply

Voltage	Frequency	Power
230V	50Hz-60Hz	600W
120V	50Hz-60Hz	600W

THE DB-2TC DRI-BLOCK®

The unit has a recessed chamber which can hold up to two interchangeable insert blocks which accept test tubes or other sample containers. The temperature is set by 3 front panel buttons with digital display.

The DB-TC is supplied with an extraction tool for removing the insert blocks from the unit.

**Specification**

	DB-2TC
Working temperature range	Ambient 0°C to 100°C
Temperature range	0°C to 100°C
Temperature stability	±0.1°C at 40°C
Temperature setting	Push button
Temperature scale graduation	No graduation
Temperature display	LED 4 digit
Uniformity within block at 37°C	±0.1°C
Uniformity within block at 100°C	±0.1°C
Block at maximum temperature	1°C
Display resolution	0.1°C
Set point accuracy	±1°C
Maximum number of blocks	2
Heater power	2 x 150W
Heat up time 30°C to 37°C	6 minutes
Heat up time 30°C to 56°C	14 minutes
Heat up time 30°C to 100°C	19 minutes
Overall size L x W x H (mm)	279 x 260 x 105
Net weight (without blocks)	2.9 kg

Electrical supply

Voltage	Frequency	Power
230V	50Hz - 60Hz	300W
110V - 120V	50Hz - 60Hz	300W

WORKING ENVIRONMENT (ALL UNITS)

The Dri-Block® units are designed to work safely under the following conditions:

Ambient temperature range 5°C to 40°C

Humidity Up to 95% relative humidity, non-condensing

Note: The control specifications are quoted for an ambient temperature of 20°C for units fitted with 2, 3 or 4 insert blocks type F3506. The specification will be closely held if the ambient temperature is in the range 10°C to 30°C. Outside this range the quoted figures may deteriorate but the unit will still work safely.

Radio frequency interference tested and passed to EN50081-1.

Immunity Tested and passed to EN50082-1.

UNPACKING

When unpacking the unit, check that the following have been removed from the packing: The unit; a guarantee card; a de-contamination certificate.


Within the guarantee period, shown on the guarantee card, we undertake to supply replacements free of charge for parts which may on examination prove to be defective, provided that the defect is not the result of misuse, accident or negligence. Any instrument requiring service under this guarantee should be sent to the supplier through whom it was purchased, or, in the case of difficulty, it should be carefully packed in its original packing and consigned, carriage paid, to Techne. Techne takes no responsibility for returned goods damaged in transit.

Returned goods will not be processed without a Returns Authorisation Number. Call the service department on +(44) (0)1785 812121 for a number. On all correspondence, please quote the serial number in full and/or the sales order number. **Please write the Returns Number on the outside of any packing.**

OPERATION

INSTALLATION

1. Select the insert aluminium alloy block(s) appropriate for the application. A list of available blocks is given on page 29 of this manual. Ensure that both the underside of the block(s) and the top of the hotplate are clean; efficient heat conduction between these two surfaces is essential.

Top view showing three insert blocks in place

2. Place the insert blocks onto the hotplate in the well of the unit and place the tubes containing the sample liquid in the blocks.
3. The heater design, temperature sensor and control circuit give good temperature control and uniformity, but it is important that there is a close fit of the tubes in the block to allow efficient heat transfer.
4. If less than the maximum number of insert blocks are fitted it is recommended that plastic half blocks are fitted to the space either side of the aluminium block. The dummy blocks have the effect of reducing heat loss, thus improving temperature stability.
5. Plug the mains cable into the socket in the rear of the unit. Connect to the mains electricity supply with the plug provided or one wired correctly for your supply. Switch the power ON from the switch located next to the mains input lead. The mains indicator on the front panel will light up (as well as the fuse light on the mains input unit).

SETTING THE OVERTEMPERATURE CUT-OUT

DB-2D, DB-2DH, DB-3DL, DB-4D

Since there is no user-adjustable cut-out, in the event of a fault occurring which could cause the heater to remain on, power is removed from the heater when the block temperature reaches approximately 110°C.

DB-3D

An adjustable overtemperature cut-out is fitted to these units. As supplied it is set to remove power from the heater should the block temperature exceed approximately 210°C. If you need to prevent the block temperature exceeding a lower value, the overtemperature cut-out may be re-adjusted downwards. Note that the cut-out should be set to trip at least 5°C above the maximum desired operating temperature.

To set the overtemperature cut-out, turn the reset button (located on the back panel of the unit) fully clockwise then:

- either:
 - Heat the block to the desired cut-out temperature and turn the reset button anti-clockwise until the unit just cuts out.
- or:
 - Heat the block to the required maximum operating temperature and turn the reset button anti-clockwise until the unit just cuts out. Now turn the reset button an additional small amount clockwise.

When the cut-out operates, the unit stops working and the overtemperature cut-out indicator lights up.

Remove the block from the unit by screwing the extractor tool into the thermometer hole and lifting the block vertically. Remember that high temperatures are dangerous; wear insulating gloves and take extreme care. Wait for the temperature to fall, perhaps by as much as 40°C, and press the reset button. The unit will now work as usual.


DB-2TC

The model DB-2TC accommodates 2 aluminium blocks with independent temperature controls so that each block may be set to a different temperature. The layout and operation of the controls for each block is as described above.

If using this model with only one block in place the side with no block should have its set temperature adjusted to 30°C.

POWER UP SCREEN

When you first switch on, the display will show the edition of the software which your unit uses. For example software issue "I" would be shown as follows:


It will display this for 1 second, then the actual temperature of the plate will be indicated.

THE FRONT PANEL CONTROLS

The front panel controls consist of three buttons for controlling the display, a four digit LED display and three indicators.


The "set" temperature

Button

The "set" temperature button displays the set temperature when pressed.


The UP ARROW Button

When the "set" temperature button is held down and the UP ARROW button is pressed, the set temperature is increased.

The DOWN ARROW


Button


When the "set" temperature button is held down and the DOWN ARROW button is pressed, the set temperature is decreased.

Speed of Change of Set Temperature


Each press of the UP ARROW or DOWN ARROW buttons will increase or decrease the set temperature by 0.1°C. If the buttons are held down the temperature change will accelerate to 5°C per second.

Power Indicator

The top indicator shows that there is power to the unit.

Heater Indicator

The next indicator shows when the heater is heating the unit. When the temperature is being set, and the new set temperature is higher than the temperature already in the unit, the heater indicator will be lit up as the unit tries to follow the set temperature. If the light is on continuously the heater is getting constant power. The only exception is described under over-temperature indicator. As the temperature approaches the set temperature the heater indicator will flash. When set temperature is reached the indicator will stay on for shorter periods. If the Dri-Block® is above the set temperature then the indicator will be off, as the heater is not getting any power.

Over-Temperature Indicator

If the unit should, for any reason, exceed the temperature set for the over-temperature cutout (see page 24), the over-temperature indicator will be lit. The heater will have been switched off and the unit will begin to cool even if the heater light is on (the light staying on or not depends on which circuit has sensed an over-temperature).

SETTING THE OPERATING TEMPERATURE


1. To display the set temperature on the digital display, press and hold the set temperature button.

To adjust the set temperature, press the set temperature button and hold it while pressing the up or down buttons. When the set temperature button is released, the measured temperature is displayed (in degrees Celsius).

2. The heater (and heater indicator) comes on if the set temperature is higher than the current block temperature.
3. When the measured temperature approaches the set temperature, the heater indicator will begin to flash. As the measured temperature stabilises the indicator will stay on for shorter periods.
4. Due to variations in heat losses with different designs of insert block, the actual temperature may vary. The units are calibrated using model F3506 insert blocks. If you need to control the temperature to a greater accuracy than the instrument's set point, place a thermometer in a designated hole in one of the blocks. For greater accuracy still, you can place the thermometer in the sample liquid in one of the test tubes. In either of these cases, it may be necessary to readjust the set temperature to achieve the precise temperature required. Allow the temperature of the unit to stabilise after each adjustment.
5. There will be a time lag between the heater plate and the insert block achieving the set temperature due to thermal contact between them. This is particularly true of the DB-3D where the display may reach temperature in 2 or 3 minutes while the insert block may take 12 minutes or more to reach temperature. The guide times are shown in the specifications.

AFTER USE

1. When you have finished heating samples, remember that parts of the unit – the tubes, blocks and associated accessories – may be very hot. Take the precautions listed earlier. We recommend that the blocks should be allowed to cool to 70°C before being removed from the heating unit. They will still have to be handled with care.
2. If you need to remove an insert block while it is hot, screw the extractor tool into the thermometer hole and lift the block out vertically. Never leave the extractor tool in the block while it is being used in the Dri-Block®.

Removing an insert block with the extractor tool.

ADDITIONAL INFORMATION

ACCESSORIES

Insert blocks are made of aluminium alloy and must be ordered separately from the heater units. The following blocks are available from Techne:

Part	Tube Size	Number of Holes	Size
F3501	Plain block	None	95 x 76 x 51
F3502	6mm	30	95 x 76 x 51
F3503	10mm	20	95 x 76 x 51
F3504	12mm	20	95 x 76 x 51
F3505	13mm	20	95 x 76 x 51
F3506	15mm	12	95 x 76 x 51
F3507	16mm	12	95 x 76 x 51
F3508	19mm	8	95 x 76 x 51
F3509	25mm	6	95 x 76 x 51
F3510	10mm cuvettes	2 channels	95 x 76 x 51
F4460	Plain block	Thermometer hole only	95 x 76 x 51
F4461	9mm/7mm	10/20	95 x 76 x 51
F4462	24mm	6	95 x 76 x 51
F4463	26mm	6	95 x 76 x 51
F4464	1.5ml Eppendorf	20	95 x 76 x 51
F4465	0.5ml	30	95 x 76 x 51
F4466	Plastic half block	Plain	95 x 37 x 50
F4470	2 ml Eppendorf	20	95 x 76 x 51
F4471	0.2ml	72	95 x 76 x 51

1. Overtemperature cutout, DB-3D and DB-2TC

The overtemperature cut-out is a sensitive mechanical device and mechanical shock can cause it to trip.

- In the event of no heater power, check the mains plug and lead, then reset the cut-out control (located at the rear of the unit).
- Repeated operation of the cut-out indicates a serious fault: you may need to return the unit to your supplier for repair.

2. Fuses

Your unit is protected by two fuses. These should only be changed by suitably qualified personnel. If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

DB-2D	230V	2 x F2A	120/100V	2 x F4A
DB-3D	230V	2 x F4A	120/100V	2 x F6.3A
DB-4D	230V	2 x F4A	120/100V	2 x F6.3A
DB-2DH	230V	2 x F3.15A	120/100V	2 x F6.3A
DB-3DL	230V	2 x F3.15A	120/100V	2 x F6.3A
DB-2TC	230V	2 x F2A	120/100V	2 x F4A

CALIBRATION

Remember that if you change the calibration from that set at the factory you may change the calibration at all temperatures. You may get different calibration with different blocks.

Make sure that you have good thermal contact between the blocks and the heater plate in the Dri-Block®. In order to ensure that the calibration you are setting is correct, you will need to use an independently calibrated probe or thermometer.

Measure the actual temperature of the particular block you want to use in the Dri-Block® using the calibrated probe or thermometer. If the calibration is not correct then you can follow this procedure. Make sure that the block and the unit are below 35°C.

- 1 Set the temperature display to 0.0°C. This will ensure that the block is at ambient temperature for the next procedure.
- 2 Hold down the Up and Down buttons and then press the set button at the same time for 5 seconds. The display will change from the set temperature, 0.0, to the block temperature to indicate the start of self-calibration.
- 3 After a period of 17 minutes for example with the DB-2D or 25.5 minutes for the DB-3D the display will change from the block temperature (eg 40.0°C) to "EEEE". Measure the actual temperature of the block using the thermometer or probe. Press the 'SET' button and either the 'UP' or the 'DOWN' button to adjust the display to the same temperature as the measured value.
- 4 Press the 'UP' and 'DOWN' buttons together to confirm the value. The display will return to the block temperature and the heater will come full on.
- 5 After a period of 29.75 minutes for the DB-2D or 34 minutes for the DB-3D the display will change from the block temperature (eg 95.0°C for the DB-2D or 195°C for the DB-3D) to "EEEE".
- 6 Measure the actual temperature of the block using the thermometer. Press the 'SET' button and either the 'UP' or 'DOWN' button to adjust the display to the same temperature as the measured value.
- 7 Press the 'UP' and the 'DOWN' buttons together to confirm the value. The display will return to the block temperature and the unit will control with the new calibration parameters.
- 8 Set the temperature to the particular temperature you require. Measure the temperature when the block temperature has stabilised. Repeat the procedure as necessary; you may need to put in an off-set.

CALIBRATION - TWIN CONTROL (DB-2TC)

These units are calibrated during manufacture using blocks designed for 16mm tubes. The different mass of aluminium in blocks designed for other sizes or types of container can cause small discrepancies between the displayed temperature and the actual temperature of the blocks.

If the block heater is to be used with other blocks and very accurate temperature setting is required then the unit should be re-calibrated for the blocks to be used.

You will need an accurately calibrated thermometer for this procedure to be worthwhile. The block heater can only be as accurate as the thermometer used to calibrate it.

Follow the procedure outlined on the next page.

1. Begin with items that are cool.
2. The blocks should be pushed firmly to the right hand side of the recess with no air gap between blocks. The right hand block will be used as the calibration block.
3. Switch on the unit and set the temperature display to 30°C.
4. Press all three buttons together and hold for 5 seconds until the amber "heat" light stops flashing and remains permanently on.
5. Self-calibration will now begin and the unit will begin to heat.
6. After approximately 30 minutes the unit will have stabilised at approximately 40°C and the display will change from showing temperature to read "EEEE".
7. Accurately measure the block temperature within the thermometer hole in the right hand block.
8. Press the set button and release.
9. Press and hold the set button and use either the up or down buttons to adjust the display temperature to equal the actual temperature measured.
10. Press the up and down buttons together to confirm the value. The amber "heat" light will stop flashing and remain permanently on.
11. Self-calibration will now continue and the unit will ramp up to approximately 5°C below maximum, i.e. either 95°C, or 195°C, depending on the model.
12. After approximately 30 minutes the unit will have stabilised and the display will change from showing temperature to read "EEEE".
13. Accurately measure the block temperature within the thermometer hole in the right hand block.
14. Press the set button and release.
15. Press and hold the set button and use either the up or down buttons to adjust the display temperature to equal the actual temperature measured.
16. Press the up and down buttons together to confirm the value.
17. The unit will now begin to control with the new calibration parameters.
18. Allow 30 minutes for the temperature to stabilise.
19. Check that both the display temperature and the thermometer temperature are now between 99°C and 101°C.
20. Press the set button and release to confirm the calibration. The set temperature will revert to 30°C.
21. Reset the set temperature to maximum (100°C) and allow a few minutes to confirm that the unit actually achieves that temperature.
22. Calibration is now complete and the unit may be used normally.

INSERT BLOCKS AND REPLACEMENT PARTS

The following parts may be obtained from your Techne® dealer if replacements or alternatives are required:

Part Number	Description
6101308	Extraction tool
FCABLEUK	Mains cable with 13amp UK plug (5 amp fuse)
FCABLEEU	Mains cable with R/A Schuko plug
FCABLEUS	Mains cable with 3 pin US plug
FMW11	Techne HI-TEMP 96 TM well plate, pack of 25
FMW12	HI-TEMP 96 TM well cover, pack of 50
F4466	Plastic half block

**UK**

Bibby Scientific Ltd.
Beacon Road, Stone
Staffordshire ST15 0SA
United Kingdom
Tel: +44 (0)1785 812121
Fax: +44 (0)1785 810405
e-mail: sales@bibby-scientific.com
www.bibby-scientific.com

France

Bibby Scientific Limited
Bâtiment Le Deltaparc Parc Silic PN2
7 rue du Canal
BP 55437 VILLEPINTE
95944 ROISSY Charles de Gaulle
France
Tel: +33(0)148 63 78 03
Fax: +33(0)148 63 78 01
e-mail: ventes@bibby-scientific.com
www.bibby-scientific.com

Middle East

Bibby Scientific Middle East Ltd.
PO Box 27842, Engomi 2433
Nicosia
Cyprus
Tel: +357 22 660 423
Fax: +357 22 660 424
e-mail: sales@bibby-scientificme.com

North & South America

Bibby Scientific US Inc.
t/a Techne Inc.
3 Terri Lane, Suite 10
Burlington, NJ 08016 USA
Toll Free (in NA): 800-225-9243
Tel: +1 609 589 2560
Fax: +1 609-589-2571
e-mail: labproducts@techneusa.com
www.techneusa.com

Asia

Bibby Scientific - Singapore
Prudential Tower, Level 26
30 Cecil Street
Singapore 049712
Tel: +65 6631 2976
Fax: +44 (0)1785 810405
e-mail: info@bibby-scientific.com
www.bibby-scientific.com