






Accessories for 6.5406.000

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5406.000

Qty.	Order no.	Description	
1 pcs	6.0341.100	Double Pt-wire electrode for coulometry	
		Indicator electrode used for coulometric Karl Fischer titration.	
1 pcs	6.0344.100	Generator electrode with diaphragm	
		Generator electrode with diaphragm for coulometric Karl Fischer titrations, SGJ 29/22. This generator electrode is well suited for: <ul style="list-style-type: none"> • samples with very low absolute water contents < 50 ppm • coulometric water content determination in oil samples • applications where ketone reagents are used • applications where low-conducting reagents are used (addition of chloroform, xylene >10%) 	
1 pcs	6.1403.030	KF adsorber tube for coulometer cell	
		For KF coulometers.	
1 pcs	6.1437.000	SGJ stopper SGJ 14	

1 pcs	6.1446.060	Stopper / B-14/15 / M10	
2 pcs	6.1448.020 Set of 5 pcs.	Septum 16 mm, 5 pieces	
1 pcs	6.1464.320	KF titration vessel / 80-250 mL / coulometric	
1 pcs	6.1903.030	Stirring bar / 25 mm Stirring bar with magnetic core, PTFE covering, length 25 mm.	
1 pcs	6.2047.020	Titration vessel holder Titration vessel holder for coulometric cells. Used with 756, 831 KF Coulometer	

1 pcs 6.2104.020 Electrode cable / 1 m / F
 For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).



1 pcs 6.2104.120 Electrode cable / 1 m / H
 With plug H. For connecting generator electrode - KF Coulometer.



1 pcs 6.2701.040 GL 18 screw cap with hole
 Made of PBT. With a hole for the 6.1448.020 septum. For the 6.5405.000 cell without diaphragm and the 6.1464.32X and 6.1465.320 titration vessels.



2 pcs 6.2713.000 Sleeve
 For the fat-free ground-joint connection with the standard ground-joint SGJ 14.



1 pcs 6.2713.010 Sleeve
 Sleeve for the fat-free ground-joint connection with the standard ground-joint SGJ 29.



1 pcs 6.2713.020 Sleeve
 For the fat-free ground-joint connection with the standard ground-joint SGJ 19.



1 pcs**6.2730.030****Stopper**

Complete with nipple and O-ring. For 6.1414.030 titration vessel lids. For KF instruments, VA measurement stands, VA Computrace and IC Dilution.

