

Note: *Indicates models provided with 50 disposable blades. (S-35)

		No.130A	No.160	No.240S	No.240D
	Body unit (L)x(W)x(T)mm	130 x 40 x 8	160 x 45 x 12	240 x 45 x 12	240 x 45 x 12
	Blade fixing plate (L)x(W)x(T)mm	50 x 20 x 7.6	80 x 20 x 5.8	80 x 20 x 5.8	80 x 20 x 5.8
	Weight (g)	184	375	545	586
Material	Body unit	SKD-11	SKD-11	SKD-11	SKD-11
	Knob parts	Stainless steel	Stainless steel	Stainless steel	Stainless steel
Finish	Body unit	Ultra-precision ground finish	Ultra-precision ground finish	Ultra-precision ground finish	Ultra-precision ground finish
Accessories	Plastic box	1	1	1	1
	Hexagonal screwdriver	1	1	1	1
	Spare hexagonal screws	2	3	3	6
	Disposable blades (S-35 Type)	1	1	1	1

Note: *Indicates models provided with 50 disposable blades. (S-35)
Other models are provided without blades.

		No.130AE	No.160E	No.240E
	Body unit (L)x(W)x(T)mm	130 x 40 x 9.2	160 x 45 x 12	240 x 47.5 x 12.0
	Blade fixing plate (L)x(W)x(T)mm	50 x 22 x 6	80 x 20 x 7.2	160 x 22 x 7.5
	Knob length x diameter(W)mm	20 x 8	22 x 10	38 x 17
	Weight (g)	216	383	680
Material	Body unit	SKD-11	SKD-11	SKD-11
	Knob parts	Stainless steel	Stainless steel	Stainless steel
Finish	Body unit	Ultra-precision ground finish	Ultra-precision ground finish	Ultra-precision ground finish
Accessories	Plastic box	1	1	1
	Fixing plate bolt	1	1	
	Spare spring	4	4	3
	L-shape wrench (for adjustment)			1
	Disposable blades (S-35 Type)	1	1	

Note: **Indicates models provided with 20 disposable blades. (C-35)
Other models are provided without blades.

		No.130C-A	No.130C-S	No.130C-AO	No.160C
	Body unit (L)x(W)x(T)mm	130 x 35 x 6	130 x 35 x 6	130 x 33.5 x 9.5	160 x 38 x 10.2
	Blade fixing plate (L)x(W)x(T)mm	50 x 29 x 4.8	50 x 29 x 4.8	50 x 29 x 4.8	60 x 29 x 4.8
	Weight (g)	128	153	188	280
Material	Body unit	SKD-11	SKD-11	SKD-11	SKD-11
	Knob parts	Stainless steel	Stainless steel	Stainless steel	Stainless steel
Finish	Body unit	Ultra-precision ground finish	Ultra-precision ground finish	Ultra-precision ground finish	Ultra-precision ground finish
	Fixing plate	Teflon coat	Teflon coat	Teflon coat	Teflon coat
Accessories	Plastic box	1	1	1	1
	Hexagonal screwdriver	1	1	1	1
	Spare hexagonal screws	2	2	2	3
	Disposable blades (C-35 Type)	1	1	1	