

4609-5 with 4610-5



Mistral Multi-Mixers®

- ♦ Eleven specialized, interchangeable attachments allow for the performance of a wide variety of laboratory mixing tasks
- ♦ Compact size conserves valuable bench top space
- ♦ Speed is adjustable from 50 to 2400 rpm with 0.120" (0.3 cm) circular motion
- ♦ Heavy duty DC motor
- ♦ Counterbalanced drive mechanism with flexible suspension-mounted on rubber feet to increase stability
- ♦ Rugged steel body

construction with baked enamel finish

- ♦ Mixer attachments are constructed of corrosion-resistant anodized aluminum and can be changed in seconds
- ♦ CE and CSA certified

ORDERING INFORMATION

Model	Volts*	Watts	Dimensions W" x L" x H" (cm)	Net Wt. lbs. (kg)	Ship Wt. lbs. (kg)
4600	120	120	8 $\frac{1}{8}$ x 10 $\frac{1}{2}$ x 7 (20.6 x 26.7 x 17.8)	18 (8.2)	22 (10)
4600-1	240				

*Both models rated 50/60 Hz.

ACCESSORIES

Model	Description
4601	Utility Tray, 6 $\frac{1}{2}$ "L x 12"W (16.5 x 30.5 cm) with rubber nonslip mat
4602	Universal Holder Attachment for flask, beakers and bottles up to 500 ml capacity
4603	Separatory Funnel Attachment for 1 funnel up to 500 ml capacity
4604	Adjustable Bottle Holder Attachment, horizontally holds bottles up to 500 ml capacity
4605	Four-Place Holder Attachment, tri-grip clamps hold round bottom flasks or bottles up to 250 ml capacity
4606	Test Tube Holder Attachment holds 36 13 mm test tubes
4607	Round Bottom Flask Holder Attachment, vertically holds 1 flask or bottle up to 500 ml capacity
4608	Rubber Pad Attachment for hand-held test tubes 4 $\frac{5}{8}$ " x 5 $\frac{1}{2}$ " (11.7 x 14 cm)
4609	Holder Attachment holds Model 4610 Sieve Set
4609-5	Holder Attachment holds Model 4610-5 Sieve Set
4610	Sieve Set, 2" (5 cm) diameter, 4 stackable holders, #25, 35, 45, 60, 80, 120, 170 and 230 mesh screens
4610-5	Sieve Set, 3" (7.6 cm) diameter, 4 stackable holders, #18, 35, 60 and 120 mesh screens
4611	Microtube platform holds 64 0.5 ml and 81 1.5 ml microtubes
4612	Replacement Foam Pad for 4611 platform