

PART OF THE **DURAN®** BOTTLE SYSTEM

TRADITIONAL


DURAN® REAGENT BOTTLE



**DURAN
WHEATON
KIMBLE**

Excellence in your hands

DURAN® REAGENT BOTTLE NARROW NECK



With precision ground neck. All glass components, therefore also suitable for storage of aggressive media, which could attack plastic parts.

Typical application: storage of liquids.

Cat. No.	Capacity (mL)	Neck	d (OD) (mm)	h (mm)	Remark	Pack Unit
Neck with standard ground joint and glass flat-head stopper						
21 165 08 09	10	10/19	28	52	Non ISO size.	10
21 165 14 02	25	12/21	36	64		10
21 165 17 02	50	14/15	42	80		10
21 165 24 07	100	14/15	52	96		10
21 165 36 09	250	19/26	70	130		10
21 165 44 08	500	24/29	86	164		10
21 165 54 04	1000	29/32	107	200		10
21 165 63 06	2000	29/32	134	248		10
21 165 73 02	5000	45/40	182	323		1
21 165 86 07	10000	60/46	227	398		1
21 165 91 06	20000	60/46	288	492		1
Neck with standard ground joint						
21 164 08 08	10	10/19	28	52		10
21 164 14 01	25	12/21	36	64		10
21 164 17 01	50	14/15	42	80		10
21 164 24 06	100	14/15	52	96		10
21 164 36 08	250	19/26	70	130		10
21 164 44 07	500	24/29	86	164		10
21 164 54 03	1000	29/32	107	200		10
21 164 63 05	2000	29/32	134	248		10
21 164 73 01	5000	45/40	182	323		1
21 164 86 06	10000	60/46	227	398		1
21 164 91 05	20000	60/46	288	492		1

