

Qty	Part No	Description	List Price	Distributor Discount	Your Price
1	ECO-555-X3-12	ECOCHYLL X3 12 LITER HIGH SPEED EVAPORATOR COMPLETE SYSTEM, TO INCLUDE:	\$39,999.00	\$9,999.75	\$29,999.25
1	ECO-555-X3	Ecodyst EcoChyll X3, 115V dual condenser double loop Halar coated stainless steel coils			
1	ECO-777-133	Condenser Body, Left Side, 120mm Flat O-Ring Flange, Vacuum Port 1" OD Pipe			
1	ECO-777-134	Condenser Body, Right Side, 120mm Flat O-Ring Flange, 1" OD Pipe			
2	ECO-777-199	120mm Quick Release Clamp			
2	ECO-777-197	120mm Viton O-ring			
1	ECO-777-005	1" PFA Elbow Compression Fitting			
1	ECO-777-015	1" PFA Tee Compression Fitting			
1	ECO-777-010	1" PFA Straight Compression Fitting			
1	ECO-777-007	Vacuum tubing, 5ft length			
1	ECO-777-008	Advanced Self Cleaning Dry Vacuum with Digital Controller			
2	ECO-777-010	1" to 1" PFA Compression Fitting			
1	ECO-777-011	Inlet feed tubing 5ft length			
1	ECO-777-115	12 Liter Evaporator Flask, 150mm Flat O-Ring Flange, 1.5" Bottom Drain Valve			
1	ECO-777-203	150mm 4 Neck Stainless Steel Evaporator Lid, 45/50 joint CN, 24/40 joint SN, 2-Elbow 1.5" TC Pipes			
1	ECO-777-034	Zero Dead Space Drain Valve, 2" Beaded Pipe, 1" Beaded Pipe Sidearm, GL-40 Thread			
1	ECO-777-003	150mm Quick Release Clamp			
1	ECO-777-004	150mm Viton O-Ring			
1	ECO-777-037	2" Coupling for Zero Dead Space Drain Valve			
1	ECO-777-016	19mm Polished Stir Shaft, overall length 30.6"			
1	ECO-777-018	5.5" O.D. Teflon Agitator, High Viscosity			
1	ECO-777-019	19mm Stirrer Bearing, 45/50			
1	ECO-777-020	Stirrer Shaft Coupling, 19mm, 26" OAL			
1	ECO-777-108	Compact Overhead Stirrer with mounting rod			
1	ECO-777-233	12 Liter Heating Mantle, 4" Diameter Hole in Bottom, 2 Circuits Rated @100 Watts, 230V, 2000 watts total			
1	ECO-777-023	High Power Digital Temperature Controller, 230V, 15 amps, 3450 watts			
1	ECO-777-024	Coiled cord type T type plug			
1	ECO-777-106	Teflon coated Temperature probe 1/4x22			
1	ECO-777-107	PTFE 1/4" Thermocouple Adapter, 24/40			
1	ECO-777-129	5 Gallon Keg Receiving System			
1	ECO-777-128	Beam Torque Wrench with 7/16" deep socket			
1	ECO-777-145	Stainless Steel Tube Fitting, Union, 1" Tube OD			
1	ECO-777-124	Positive Displacement Pump for emptying keg			
		Estimated Shipping and Handling	\$1,500.00	\$0.00	\$1,500.00

Sales Amount: **\$31,499.25**

Total Discount: **\$9,999.75**

Total Amount: **\$31,499.25**

Qty	Part No.	Optional Items Description	List Price	Distributor Discount	Your Price
1	ECO-777-104	Insulating Mantle Top, for 12 Liter RBF, 6.5" Diameter Hole, Jacket is an Insulator Only	\$600.00	\$90.00	\$510.00