

Accessories and Options

A range of different containers and lids is available for the GRINDOMIX GM 200 and GM 300 for optimum adaptation to a particular application. These include:

- **Patented gravity lid**
Automatically adjusts the grinding chamber volume to the changing sample volume.
- **Gravity lid with overflow channels**
Ideally suited to homogenize samples with a high water content.
- **Stainless steel container**
Minimum wear when hard sample materials are processed.
- **Reduction lid**
Reduces the chamber volume of the GM 200.
- **Serrated blade knife**
Used for particularly tough samples such as fatty, streaky meat.
- **Accessories for cryogenic grinding**
Applications with dry ice are carried out in the GM 300 with a full metal knife and a special lid.



Accessories
GRINDOMIX GM 300

