

**Thermo IEC Centra®
GP8 Centrifuges**

GSA

Quiet Centra GP8 centrifuges give you the same capacity, separating power, and versatility of floor model centrifuges, but at a fraction of the space and cost!

Select spin times of 1 second to 360 minutes with exact precision. Set speeds from 500 to 6,000 rpm in 10 rpm increments, or set by required g-force directly. LED display gives both actual and set values with the touch of a key. "At Speed" timer mode enhances the accuracy of time-sensitive runs.

Membrane key-pad offers protection from splashes and offers 5 acceleration rates, 5 deceleration rates, and a coast position—choose fast for rapid processing, gentle for delicate samples. Microprocessor control stores 35 of your most frequently used parameters.

Fail-safe cover interlock system disables centrifuge when the cover is open; releases only at the end of a run. Rotors and adapters sold separately.

GP8R

**Two Year
Warranty**Zero rotation safety
interlock prevent lid
openingExclusive "At Speed"
timer only begin
countdown after
program speed has
been reachedStore up to 35
programs in the
memorySelect from
5 deceleration
and 5 acceleration
modes

GP8

Catalog Number	Description	List Each	Discount Each
TX23550	Centra GP8, 115 V	\$6,189.00	\$5,523.00
TX23551	Centra GP8R, Refrigerated, 115 V	7,849.00	7,173.00

For 230 V add -1 to catalog number (i.e. TX23500-1)

SPECIFICATIONS Thermo IEC Centra® GP8 Centrifuges

Maximum Speed :	5,100/3,400 rpm (Fixed Angle/Swing Bucket)
Maximum RCF:	4,190/2,550 x g (Fixed Angle/Swing Bucket)
Speed Range:	500 to 6,000 rpm in 10 rpm increments
Timer:	1 sec 360 min by 1 min + hold mode
Temperature Range:	-5° to 40°C (GP8R only)
Dimensions (W x H x D"):	GP8/GP8R Cover Closed: 17 x 24 x 23 Cover Open: 36 x 24 x 23
Shpg Weight (lb):	245/315(GP8R)

Rotors for Thermo IEC Centra GP8 Centrifuges


Use with Centrac Adapter below.

**Microplate
Rotor**

4-place microplate rotor includes 4 microplate carriers. Use with up to 8 plates; 3,500 rpm @ 2,340 x g.

TX23556E \$2,837.00

Discount 2,539.00

**Aerosol
Containment
8 x 50 mL
Angle Rotor**

45° fixed angle rotor includes sealed shields and 45° angle rack. 5,100 rpm @ 4,190 x g.

TX23556D \$3,083.00

Discount 2,858.00

**4-place
Swinging
Bucket
Rotor**

Includes sealed buckets domes. 3,400 rpm @ 2,550 x g.

TX23556C \$2,284.00

Discount 2,130.00

**Sealed
4-place
Rotor**

GSA

Includes sealed buckets and sealed transparent windshield. 4,600 rpm @ 4,550 x g.

TX23556A \$3,791.00

Discount 3,535.00

**Centrac Adapters/
Microplate Carriers**

GSA

Color-coded Centrac adapters let your Centra GP8 Centrifuge accommodate many different tube sizes. Tube sites are numbered for positive sample identification. Ready to use—adapters do not require separate handles or cushions. Tubes are fully supported to reduce chance of breakage. Autoclavable.

Accepts Container	Max. in rotor A B	Catalog Number	Pack of 2	Discount Pack
0.25 & 0.4 mL; 0.5 and BD Microtainers; 4 to 5 mL (10 x 75/12 x 75 mm)	148 148	TX23557A	\$183.00	\$163.00
1.5 mL Microtubes; 7 to 10 mL (13 x 75/13 x 100mm)	108 108	TX23557B	195.00	174.00
4 mL Shell Files; 7 to 12 mL (16 x 75/16 x 100 mm); 10 to 18 mL (16 x 125/17 x 120 mL)	76 76 28 76	TX23557C	205.00	185.00
15 mL (Falcon/Corning Conical)	24 48	TX23557D	183.00	163.00
15/50 mL (Conical); Individually Sealed	12 12	TX23557E	174.00	162.00
50 mL (Round, Flat or Conical Bottom)	28 28	TX23557F	193.00	173.00
50 mL (Falcon/Corning Conical, #5805 Adapter)	20 20	TX23557G	261.00	173.00
50 mL (Round Bottom); 100 mL (Corning 8460-100)	— 8	TX23557H	298.00	265.00
140 mL (Bellco 3045-00140); 200 mL (NUNC 376813 conical); 175/225 mL (Falcon 2074/2075); 250 mL (Corning 25350 Conical); 250 mL (IEC 2050, 2051, 2261, 1625)	— 4 — 4 4 4	TX23557J	183.00	163.00
Microplate Carrier for (86 x 128 mm)	4 — — 8	TX23557K TX23557M	833.00 756.00	738.00 679.00