The Thermo Scientific Orion 2-Star pH benchtop meter is easy to use and value priced for your general purpose pH measurements in the lab.

Thermo Scientific Orion 2-Star

pH Benchtop Meter



Features and Benefits

- The Thermo Scientific Orion SMART STABILITY™ and SMART AVERAGING™ functions for automatically optimized accuracy, precision and response time
- Built-in pullout guide for quick access to calibration instructions
- Up to 3 point pH calibration
- Automatic buffer recognition of USA, NIST and DIN buffers
- E_H autocalibration for accurate ORP calibration and readings
- Data storage for up to 50 points with time and date stamp
- RS232 port for data downloading
- IP54-rated housing
- Includes universal power supply
- 3 year meter warranty

Thermo Scientific Orion 2-Star pH Benchtop Meter

pH	Range	0.000 to 14.999
	Resolution	0.1, 0.01, 0.001
	Relative Accuracy	± 0.002
	Auto-Buffer Recognition	4.01, 7.00, 10.01/ 1.68, 6.86, 9.18, 12.46
	Calibration Points	1, 2 or 3
	Selectable Calibration Sets	USA or EURO
	Custom Buffers Accepted	Yes
Millivolts	Range	± 1999.9
	Resolution	0.1
	Relative Accuracy	± 0.2 or 0.05% of reading, whichever is greater
Relative Millivolts	Range	± 1999.9
	Resolution	0.1
	Relative Accuracy	± 0.2 or 0.05% of reading, whichever is greater
	ORP Mode with E _H Autocalibration to Standard Hydrogen Electrode	Yes
Temperature	Range	- 5 to 105 °C
	Resolution	0.1 up to 99.9 °C, 1.0 over 99.9 °C
	Relative Accuracy	± 0.1 °C



Meter Onl	У		
Cat. No.	Package		
1111000	Includes meter, universal power adapter and user guide		
Meter Kits	S		
Cat. No. / Application	Sensor (s)	Accessories and Solutions	
1111001 / Fresh Water	9157BNMD Refillable epoxy pH/ATC Triode electrode	 Swing arm electrode stand pH 4.01, 7.00, and 10.01 buffers and pH electrode storage solution (pint bottles) 	
Meter wit	h Electrode and Stand		
Cat. No. / Application	Sensor (s)		
1111201 / Wastewater, Clog-Free	9165BNWP Sure-Flow refillable epoxy pH electrode		
1111202 / Fresh Water	• 9157BNMD Refillable epoxy pH/ATC Triode electrode		
1111203 / <i>TRIS</i>	9156APWP AquaPro low maintenance epoxy pH electrode		
1111204 / <i>Fresh Water</i>	9142BN Refillable glass pH electrode927007MD Stainless steel ATC probe		
1111205 / <i>Fresh Water</i>	9142BN Refillable glass pH electrode		
1111206 / <i>Fresh Water</i>	• 9145BN Low maintenance epoxy pH electrode		
1111207 / <i>Fresh Water</i>	• 9107BNMD Low maintenance epoxy pH/ATC Triode electrode		
1111208 / TRIS, Education	9107APMD AquaPro low maintenance epoxy pH/ATC Triode electrode		

©2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.



North America 166 Cummings Center Beverly, MA 01915 USA Toll Free: 1-800-225-1480 Tel: 1-978-232-6000 Dom. Fax: 1-978-232-6015 Int'l Fax: 978-232-6031

Europe

Denmark House, Angel Drove Ely, Cambridgeshire England, CB7 4ET Tel: 44-1353-666111 Fax: 44-1353-666001

Asia Pacific

Blk 55, Ayer Rajah Crescent #04-16/24, Singapore 139949 Tel: 65-6778-6876 Fax: 65-6773-0836



S-2STAR-E 0608 RevB

