



BIOMEDICAL LABORATORY  
INSTRUMENTS

Manufacturing & Servicing Laboratory  
Instruments Since 1975

## Rocker Platform

Not only the speed but also the pitch can be changed!



### Model RP-50

(Approx. Shipping Wt. 18.00 lbs.)

## Features

- The platform is interchangeable 12" x 12", 16" x 16", 20" x 20" Can work with the Orbital Mixer, and Battery Rocker (16" or smaller) Platforms.
- The noise free motor operates from 0 - 100 RPM
- Pitch from 2° to 40°
- The platform's working area is covered with a silicon rubber surface.
- Platform can be used with Screw In Flask Holders.

## Specifications

- 110V or 220V
- Max Weight Limit 20 lbs.
- Dimensions: 11" x 9" 4"

PO Box 1835 • Rockville, MD 20849-1835 • Ph. 301.208.0445 • Fax 301.208.0146  
Customer Service Available Mon.- Fri. 8am - 4:30pm (ET) • email: [elmeco@elmecoengineering.com](mailto:elmeco@elmecoengineering.com)